

May 22, 2024

Melatonin ELISA

DOI

dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1

daniel.dautan daniel¹, Per Svenningsson¹

¹Karolinska Institute Stockholm

ASAP Collaborative Rese...

Kaplift Protocols



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1**

Protocol Citation: daniel.dautan daniel, Per Svenningsson 2024. Melatonin ELISA. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbx35ylpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 15, 2024

Last Modified: May 22, 2024

Protocol Integer ID: 95309

Keywords: ASAPCRN, melatonin, sleep

Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: 020608



Abstract








Measurement of mouse plasma melatonin using ELISA Kit (Enzo Life Sciences, ENZ-KIT150-0001, NY, US) according to manufacturer instructions.

Materials

Enzo Life Sciences, ENZ-KIT150-0001












Melatonin Extraction





- 1 Mix  200 μL of plasma with an equal volume of cold Ethyl Acetate. Vortex gently.
- 2 Allow layers to separate on ice for  00:03:00 . Vortex again and incubated  On ice for  00:02:00 . 5m
- 3 Spin samples at 1000g for  00:10:00 at  4 °C . Transfer the organic layer to a new tube. 10m
- 4 Dry samples and resuspend in  220 μL of 1X stabilizer.

ELISA

2h

- 5 Add  100 μL of standards' working solutions and samples to provided 96-well plate in duplicates.
- 6 Immediately after add  50 μL of melatonin tracer to each well (except blanks) followed by  50 μL of 1X melatonin antibody (except blanks).
- 7 Cover plates with the provided plate sealer. Incubate for  01:00:00 at  37 °C with 500 rpm shaking. 1h
- 8 Decant the solution from each well. Add  400 μL of wash solution to each well.
- 9 Decant the solution from each well and pat dry against clean absorbent paper.
- 10 Repeat wash step 3 times.
- 11 Add  200 μL of melatonin conjugate solution to each well (except blanks). Cover plate with the sealer. Incubate for  00:30:00 at  Room temperature . 30m



- 12 Perform the wash step as described above. Add  200 μ L of substrate reagent to each well. Incubate for  00:30:00 at  37 °C protected from light. 30m
- 13 Add  50 μ L of stop solution to each well. Measure the optical density using a micro-plate reader with absorbance set to 450 nm.