



Aug 10, 2020

Conjugation of peptide fragment 579-601 of HIV-gp41 or fragments 308-331 or 421-438 of the HIV-gp120 with keyhole limpet hemocyanin (KLH).

Angel A Justiz-Vaillant¹¹University of the West Indies St. Augustine**1** Works for me dx.doi.org/10.17504/protocols.io.bjhukj6w[University of the West Indies](#) angel.vaillant@sta.uwi.eduAngel Justiz-Vaillant
University of the West Indies St. Augustine

ABSTRACT

Chemical synthesis facilitates the generation of peptides which are very difficult to express in bacteria, peptide/protein backbone modification, the incorporation of unnatural amino acids, and the production or synthesis of D-proteins.

The C-terminal cysteine can be added to the amino acidic sequences of HIV peptides (fragment 579-601 of the HIV-gp41 and fragments 308-331 or 421-438 of the HIV-gp120). These peptide fragments were dimerized by cysteine oxidation with dimethyl-sulfoxide [1] to facilitate their conjugation to keyhole limpet hemocyanin that acts as a carrier protein.

Reference:

1. Tam JP, Wu CR, Liu w, Zhang JW (1991) Disulfide bond formation in peptides by dimethyl sulfoxide. Scope and applications. J Am Chem Soc 113: 66576662.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Angel Justiz-Vaillant. Conjugation of peptide fragment 579-601 of HIV-gp41 or fragments 308-331 or 421-438 of the HIV-gp120 with keyhole limpet hemocyanin (KLH).. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjhukj6w>

DOI

dx.doi.org/10.17504/protocols.io.bjhukj6w

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Conjugation of peptide fragment 579-601 of HIV-gp41 or fragments 308-331 or 421-438 of the HIV-gp120 with keyhole limpet hemocyanin (KLH).. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjhukj6w>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Angel Justiz-Vaillant. Conjugation of peptide fragment 579-601 of HIV-gp41 or fragments 308-331 or 421-438 of the HIV-gp120 with keyhole limpet hemocyanin (KLH).. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjhukj6w>

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 10, 2020

LAST MODIFIED

Aug 10, 2020

PROTOCOL INTEGER ID

40212

GUIDELINES

The Protocol has a high level of reproducibility and has worked for many other HIV peptides.

MATERIALS

NAME	CATALOG #	VENDOR
Glutaraldehyde EM Grade 25%	G5882-50ML	Sigma Aldrich
10mg KLH (Keyhole Limpet Hemocyanin) (Immunological Grade)	786-088	G-Biosciences
Peptide 579-601 of HIV-gp41		
Peptide 308-331 HIV-gp120		
Peptide 421-438 HIV-gp120		

SAFETY WARNINGS

Side effects include skin irritation, nausea, headache, and shortness of breath.

BEFORE STARTING

Ensure that the HIV peptide requires a terminal cysteine, so a disulfide bond formation in the peptides by dimethyl sulfoxide is needed before the conjugation to KLH. For example, fragment 254-274 of the second conserved domain of HIV gp120 has a terminal cysteine, and this peptide can be conjugated to KLH without previous treatment with dimethyl sulfoxide.

- 1 These peptide fragment (579-601 from HIV-gp41) is dimerized by cysteine oxidation with dimethyl-sulfoxide. The HIV peptide is dissolved in 5% acetic acid to a final concentration of 5.1 mg/ml.
- 2 The pH of the medium is adjusted to 6 with 1 M (NH₄)₂CO₃.
- 3 Dimethyl-sulfoxide is added to 20% of the final volume, and after four hours at room temperature (RT), the solute is extracted.
- 4 Then, the peptide is dissolved in 3 ml 5% trifluoroacetic acid and precipitated with 35 ml cold ether.
- 5 The precipitate is dialyzed against 1.2 liters of deionized water, pH 7 at 4°C overnight.
- 6 Then, 1 mg of keyhole limpet hemocyanin (KLH) is diluted in 2.1 ml 0.1 M borate buffer (1.24 g boric acid, 1.90 g sodium tetraborate, pH 10, in 500 mL deionized water).
- 7 In a 20 ml glass tube, with a gentle stirring, 1.1 µmol of the HIV synthetic peptide (with C-terminal cysteine added) and 0.22 mL 0.3% glutaraldehyde solution (ACS reagent grade, pH 5.5, Sigma-Aldrich) are slowly mixed at RT and left to

stand for 1.50 hrs.

- 8 When a yellow coloration is observed this indicates that the conjugation process is successful.
- 9 To blocking the excess of glutaraldehyde, 0.26 ml of 1 M glycine (Sigma-Aldrich) is added.
- 10 The mix is left for 32 min at RT.
- 11 The HIV-hemocynin conjugate is then dialyzed against 1.3 liters 0.1 M of borate buffer, pH 8.4 through the night at 4°C.
- 12 Then the previous buffer is used to dialyze the preparations for 8 hrs at 4°C.
- 13 The dialysates is stored at 4°C until further use.