

May 22, 2024 Version 3

HTTPM : gDNA extraction V.3

DOI

dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v3

Antoine Champie¹, Amélie De Grandmaison¹, Sebastien Rodrigue¹

¹Université de Sherbrooke



Antoine Champie

Université de Sherbrooke

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v3

External link: <https://doi.org/10.1371/journal.pone.0283990>

Protocol Citation: Antoine Champie, Amélie De Grandmaison, Sebastien Rodrigue 2024. HTTPM : gDNA extraction. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7yzz3gwz/v3> Version created by **Antoine Champie**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 24, 2022

Last Modified: May 22, 2024

Protocol Integer ID: 100192

Abstract

Part of the HTTPM protocol dedicated to the extraction of gDNA from transposon mutated cell pellets.



Attachments



HDTM Protocol-2.pd...

298KB

Image Attribution

Make with BioRender.com



Materials

▪ **Homemade DNA lysis Buffer :**

A	B
Component	Amount for 1000ml of solution
CTAB 2%	20g
1,5M Guanidine HCl	143,2g
10mM Tris HCl	1,57g

Mix well and adjust pH to 8.0.

▪ **Homemade wash solution :**

A	B
Component	Amount for 1000ml of solution
Ethanol 100%	800ml
Tris HCl 1M pH 8,0	10ml
NaCl 4M	25ml
EDTA 0,5M	2ml

Mix well and adjust pH to 8.0.

▪ **Elution Buffer (Low TE Buffer):** 10 mM Tris-HCl (pH 8.0) + 0.1 mM EDTA

Solutions for plate regeneration, from this protocol : (1)<https://doi.org/10.1016/j.ab.2008.10.021>.

▪ **NaOH 1N + Triton X100 0,15% (v/v)**

A	B
Component	Amount for 1000ml of solution
Water	960ml



A	B
NaOH	40g
Triton X-100	1,5ml

Mix well and store in a base resistant container.

▪ **HCl 1.5N + Triton X100 0,15% (v/v)**







A	B
Component	Amount for 1000ml of solution
Water	873,5ml
HCl Stock (37%)	125ml
Triton X-100	1,5ml

Mix well and store in an acid resistant container.






DNA extraction

2h 5m

- 1 Prepare the lysis solution by adding  165 μL of proteinase K to  66 mL of homemade lysis buffer and mix well.
- 2 Add  600 μL of lysis solution to each well of the deep-well plate and resuspend the pellet.
- 3 Cover with an adhesive aluminum foil and incubate at  55 $^{\circ}\text{C}$ for  01:00:00 .
- 4 While still warm, add  260 μL of ethanol 100%, without overmixing.

Note

Overmixing will result in DNA agglomeration and difficulty with the extraction.

- 5 Transfer immediately to a deep-well plate fitted with an array of silica columns.
- 6 Centrifuge twice at  3270 x g, 00:10:00 .
- 7 Discard flowthrough and add  500 μL of wash solution.
- 8 Centrifuge at  3000 x g, 00:10:00
- 8.1 Repeat steps 7 and 8.
- 9 Discard flowthrough.

1h


10m



10m



- 10 Centrifuge at 3000 x g, 00:05:00 to eliminate traces of wash solution. 5m
- 11 Discard flowthrough.
- 12 Add a collector plate between the silica column array and the deep-well plate.
- 13 Add 50 μL of low TE to the silica matrix in each well.
- 14 Cover with an adhesive aluminum foil and incubate at 55 °C for 00:15:00 . 15m
- 15 Centrifuge at 3270 x g, 00:05:00 . 5m
- Silica array regeneration (Optional)** 1h 5m
- 16 Put the contaminated silica array on an empty deep-well plate.
Add 150 μL of 1N NaOH + 0.15%(v/v) Triton X-100 to each well.
- 17 Incubate at Room temperature for 00:05:00 5m
- 18 Centrifuge 3000 x g, 00:02:00 2m
- 19 Add 200 μL of 1,5N HCl+ 0,15% (v/v) Triton X-100 to each well.
- 20 Incubate at Room temperature for 00:30:00 30m
- 21 Centrifuge 3000 x g, 00:02:00 2m



22 Add  150 μL of 1N NaOH + 0,15%(v/v) Triton X-100 to each well.

23 Incubate at  Room temperature for  00:05:00

5m

24 Centrifuge  3000 x g, 00:02:00

2m

24.1 Collect the flowthrough in a beaker. Neutralize pH if needed and dispose of the flow through.

25 Add  600 μL of ddH₂O to each well.

26 Centrifuge  3000 x g, 00:05:00

5m

27 Silica columns array are ready to be reused.