



Sep 30, 2020

# Mouse Liver Single-Cell Dissociation with Cold Protease

Sanjay Subramanian<sup>1</sup>, Stacey S Huppert<sup>2,3</sup><sup>1</sup>[Division of Gastroenterology, Hepatology & Nutrition, Cincinnati Children's Hospital Medical Center, Cincinnati, OH, United States]

;

<sup>2</sup>[Division of Gastroenterology, Hepatology & Nutrition, Cincinnati Children's Hospital Medical Center, Cincinnati, OH, United States;<sup>3</sup>Department of Pediatrics, University of Cincinnati College of Medicine, Cincinnati, OH, United States]

1

Works for me

[dx.doi.org/10.17504/protocols.io.bhpcj5iw](https://dx.doi.org/10.17504/protocols.io.bhpcj5iw)

Sanjay Subramanian

## ABSTRACT

This protocol is used to dissociate murine liver cells at 4 °C in an effort to preserve transcriptomic expression profiles in cells to be used for single-cell RNA sequencing and single-cell ATAC sequencing. An expedited protocol length limits the time of exposure to dissociative solutions in order to retain natural cellular distribution and viability while limiting transcriptomic changes arising from injury response, especially in fragile cell populations. This protocol has been tested and works in mice aged 2-4 weeks, dissociating cells from the left lobe for single-cell RNA sequencing using the 10x Chromium system.

This project was supported in part by NIH P30 DK078392 (Gene Analysis Core) of the Digestive Diseases Research Core Center in Cincinnati, R01DK107553, and the Cincinnati Pediatric Cell Atlas Center funded by Cincinnati Children's Research Foundation to Stacey S. Huppert

DOI

[dx.doi.org/10.17504/protocols.io.bhpcj5iw](https://dx.doi.org/10.17504/protocols.io.bhpcj5iw)

## PROTOCOL CITATION

Sanjay Subramanian, Stacey S Huppert 2020. Mouse Liver Single-Cell Dissociation with Cold Protease.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bhpcj5iw>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 18, 2020

## LAST MODIFIED

Sep 30, 2020

## PROTOCOL INTEGER ID

38340

## GUIDELINES

### Bacillus licheniformis enzyme mix

*The Bacillus licheniformis mix is made up in DPBS at a concentration of 100mg/mL and stored at -80 °C in 100µL aliquots*

To make 10mL of 1mg/mL Bacillus licheniformis mix:

9900µL DPBS

100µL Bacillus licheniformis 100mg/mL mix

### 10x PBS (1L)

80g NaCl  
2g KCl  
11.5g Na<sub>2</sub>HPO<sub>4</sub> 7H<sub>2</sub>O  
2g KH<sub>2</sub>PO<sub>4</sub>  
1L diWater

### 1x PBS for Flush

Make up 5 mL of 1x PBS by diluting 500µL 10x PBS in 4.5mL MilliQ water

### Required Equipment

Refrigerated centrifuge  
Tube rotator (or similar device)

### MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">Protease from Bacillus Licheniformis</a>	P5380	<a href="#">Sigma</a>
<a href="#">3 mL syringes</a>	BD #309657	<a href="#">BD Biosciences</a>
<a href="#">Falcon® Cell Strainers, Mesh size: 70um; white</a>	087712	<a href="#">Thermo Fisher</a>
<a href="#">William's E Medium, no glutamine</a>	12551032	<a href="#">Thermo Fisher</a>
<a href="#">DPBS, calcium, magnesium</a>	14040133	<a href="#">Thermo Fisher</a>
<a href="#">100um Cell Strainer</a>	22363549	<a href="#">Fisher Scientific</a>
<a href="#">EqualFetal Bovine Serum</a>	EF-0500-A	<a href="#">Atlas Biologicals</a>
<a href="#">PrecisionGlide Needle 23G</a>	BD #305120	<a href="#">BD Biosciences</a>
<a href="#">BD Intramedic™ Polyethylene Tubing PE-50</a>	BD #427516	<a href="#">Fisher Scientific</a>

### ABSTRACT

This protocol is used to dissociate murine liver cells at 4 °C in an effort to preserve transcriptomic expression profiles in cells to be used for single-cell RNA sequencing and single-cell ATAC sequencing. An expedited protocol length limits the time of exposure to dissociative solutions in order to retain natural cellular distribution and viability while limiting transcriptomic changes arising from injury response, especially in fragile cell populations. This protocol has been tested and works in mice aged 2-4 weeks, dissociating cells from the left lobe for single-cell RNA sequencing using the 10x Chromium system.

This project was supported in part by NIH P30 DK078392 (Gene Analysis Core) of the Digestive Diseases Research Core Center in Cincinnati, R01DK107553, and the Cincinnati Pediatric Cell Atlas Center funded by Cincinnati Children's Research Foundation to Stacey S. Huppert

### BEFORE STARTING

Prepare essential solutions and leave on ice

Cool centrifuges to 4 °C

To make cannulation apparatus:

- Cut ~12 inch length of PE-50 tubing and make an angled cut at one end with a razor blade to serve as the cannulation tip.
- Carefully insert 23-gauge needle tip into opposite end of tubing, making sure not to poke a hole in tubing
- If a smaller cannula is needed to fit into the portal vein, the PE-50 can be warmed by rolling between fingers and then stretched by pulling. Additionally, PE-10 can be used for very small vessels.

Tissue Isolation 20m

1 Anesthetize mouse

- 1.1 Fill syringe with 1x PBS and twist on to needle with attached cannulation apparatus. Flick syringe and push plunger to ensure 1x PBS has filled the tubing and there are no air bubbles in the syringe or tubing
- 1.2 Open abdominal cavity and locate main trunk of portal vein. Make a small incision to allow for insertion of angled tip of cannulation apparatus without cutting all the way through portal vein
- 1.3 Make incision in inferior vena cava to allow blood and perfusing solutions to flow out

- 2 Cannulate the portal vein and flush 3-5mL ice cold 1x PBS through portal vein until liver tissue is sufficiently cleared of blood

 **5 mL ice-cold 1x PBS**

- 3 Switch PBS syringe out for one with *Bacillus licheniformis* mixture while keeping the portal vein cannulated. Ensure no bubbles are in the system and then perfuse 3-5mL ice cold 1mg/mL *Bacillus licheniformis* mixture through the portal vein


 **5 mL ice-cold *Bacillus licheniformis* enzyme mix**

- 4 Dissect out desired liver tissue and mince thoroughly with razor blade in petri dish with an additional 1mL *Bacillus licheniformis* mix; start timer for 1 hour

 **1 mL ice-cold *Bacillus licheniformis* enzyme mix**

#### Dissociation and Filtration


1h 5m





- 5 Transfer tissue-enzyme slurry to 5mL capped tube and incubate, rotate, or rock at 4 °C until timer reaches 1 hour  
 **01:00:00 at 4 °C**
- 6 Mince slurry again with razor blade, making sure to separate or scrape larger chunks to release the maximum number of cells
- 7 Pass slurry through a 100µM filter into a 50mL conical tube
- 8 Bring volume in tube up to 10mL by slowly passing ice cold Williams' E media with 20% FBS through filter to dislodge any sticking cells and stop the enzyme activity

 **10 mL Williams' E 20% FBS with slurry**

#### Centrifugation and Resuspension

15m

- 9 Centrifuge cells at 250g for 5min at 4 °C  
 **00:05:00 centrifuge 250g at 4 °C**



- 10 Discard supernatant and resuspend pellet in 5mL Williams E with 10% FBS  
 **5 mL resuspend Williams' E 10% FBS**
- 11 Filter suspension through a 70µM filter into a new 50mL conical tube; rinse filter with an additional 5mL of media  
 **5 mL Williams' E 10% FBS**
- 12 Repeat steps 9-10  
 **00:05:00 centrifuge 250g at 4 °C**  
 **5 mL resuspend Williams' E 10% FBS**



Final resuspension volume may be greater or less than stated and should be judged accordingly based on pellet size

Quality Control

10m

- 13 Analyze cell morphology using hemocytometer and Trypan blue  
 Viable hepatocytes tend to take up trace amounts of Trypan blue
- 14 Adjust cell concentration as needed (700-1200 cells/µL for 10x Chromium system)  
 **1000 cells/µL**