



FEB 29, 2024

T cell differentiation from mice spleen tissue

Ningbo Zheng^{1,2}, Weiyi Peng^{1,2}

¹Department of Biology and Biochemistry, University of Houston, Houston, TX, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network

ASAP Collaborative Research Network

LRRK2



sarfraz.ahmed

OPEN ACCESS



ABSTRACT

This protocol is for T cell differentiation from mice spleen tissue to investigate T cell function (including differentiation and proliferation) in vitro.

DOI:

dx.doi.org/10.17504/protocols.io.261gedy8ov47/v1

Protocol Citation: Ningbo Zheng, Weiyi Peng 2024. T cell differentiation from mice spleen tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261gedy8ov47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 29, 2024

Last Modified: Feb 29, 2024

MATERIALS

PROTOCOL integer ID: 95978

Funders Acknowledgement:

Aligning Science Across

Parkinson's (ASAP)

Grant ID: Grant ID: ASAP-000312

1. Pre coating of 96 well plate with antibodies using following reagents

A	B	C
Reagent	Final concentration	Amount
PBS	n/a	100mL/well
Anti-CD3e (0.5mg/mL)	2.5mg/mL	0.5mL/well
Anti-CD28 (1mg/mL)	1mg/mL	0.1mL/well

2. Prepare FACS buffer and T cell culture media (See Appendix 1 below)

A	B	C
FACS buffer		
Reagent	Final concentration	Amount
PBS	n/a	490mL
FBS	2%	10mL
Total	n/a	500mL

Note: Can be prepared in advance and stored at 4°C for at least 4 months. Filtration through a 0.2 µm vacuum filter is recommended.

A	B	C
T cell culture media		
Reagent	Final concentration	Final concentration
MEMα	n/a	429mL
FBS	10%	50mL
GlutaMAX 100×	2×	10mL
HEPES, 1M	0.02M	10mL

A	B	C
2-mercaptoethanol (55 mM)	55μM	0.5mL
Normocin, 1g (50mg/ml)	50μg /ml	0.5mL
Total	n/a	500mL

Note: Can be prepared in advance and stored at 4°C for up to 1 month.

2. Prepare T-cell differentiation cocktails [1-4] (Appendixes 2 below)

A	B	C
Th0 differentiation cocktail (2.5×)		
Reagent	Final concentration	Amount
T cell culture media	n/a	1mL
IL-2 (5×10 ⁵ U/mL)	100U/mL	0.5mL

A	B	C
Th1 differentiation cocktail (2.5×)		
Reagent	Final concentration	Amount
T cell culture media	n/a	1mL
IL-2 (5×10 ⁵ U/mL)	100U/mL	0.5mL
IL-12 (10μg/mL)	4ng/mL	1mL
Anti-IL-4 (10mg/mL)	10μg/mL	2.5mL

A	B	C
Th2 differentiation cocktail (2.5×)		
Reagent	Final concentration	Amount
T cell culture media	n/a	1mL

A	B	C
IL-2 (5×10 ⁵ U/mL)	100U/mL	0.5mL
IL-4 (10µg/mL)	10ng/mL	2.5mL
Anti-IFNγ (10mg/mL)	10µg/mL	2.5mL

A	B	C
Th9 differentiation cocktail (2.5×)		
Reagent	Final concentration	Amount
T cell culture media	n/a	1mL
IL-2 (5×10 ⁵ U/mL)	100U/mL	0.5mL
IL-4 (10µg/mL)	10ng/mL	2.5mL
Anti-IFNγ (10mg/mL)	10µg/mL	2.5mL
TGF-β (5µg/mL)	1ng/mL	0.5mL
IL-1β (10µg/mL)	10ng/mL	2.5mL

A	B	C
Treg differentiation cocktail (2.5×)		
Reagent	Final concentration	Amount
T cell culture media	n/a	1mL
IL-2 (5×10 ⁵ U/mL)	300U/mL	1.5mL
Anti-IFNγ (10mg/mL)	10µg/mL	2.5mL
TGF-β (5µg/mL)	5ng/mL	2.5mL
Anti-IL-4 (10mg/mL)	10µg/mL	2.5mL

Note: Cocktails can be prepared in advance and stored at 4°C for several days, seed 100mL in each well.

3. Reagents or Resources

A	B	C
MEMα	Fisher Scientific	Cat#32-571-101
FBS	R&D Systems	Cat#S11150
GlutaMAX 100×	Fisher Scientific	Cat#35050061
HEPES	Fisher Scientific	Cat#MT25060CI
2-mercaptoethanol	Fisher Scientific	Cat#21985023
Anti-CD3e	Cytek Biosciences	Cat#145-2C11
Anti-CD28	BD Biosciences	Cat#553294
IL-2	Prometheus Laboratories	Cat#NDC65483-116-07
IL-12	PeproTech	Cat#210-12
TGF-β1	BioLegend	Cat#763104
IL-1β	R&D Systems	Cat#401-ML/CF
IL-4	PeproTech	Cat#214-14
Anti-IL-4	BioXCell	Cat#BE0045
Anti-IFNγ	BioXCell	Cat#BE0055
Mouse CD4+ T-cell Isolation Kit	STEMCELL Technologies	Cat#19852
Mouse IFN-γ DuoSet ELISA	R&D Systems	Cat#DY485
Mouse IL-4 DuoSet ELISA	R&D Systems	Cat#DY404

A	B	C
Mouse IL-9 DuoSet ELISA	R&D Systems	Cat#DY409
Foxp3/transcription factor staining buffer set	Fisher Scientific	Cat#00-5523-00
Anti-CD4-eFluor 450	eBioscience	Cat#48-0042-82
Anti-CD25-PerCP	TONBO Biosciences	Cat#65-0251
Anti-Foxp3-PE	Fisher Scientific	Cat#12-5773-82
PMA	Sigma-Aldrich	Cat#P1585

Day 0

- 1 Prepare precoated 96-well plate using antibodies described in materials section.

Day 1

- 2 Prepare single-cell suspensions of mice spleen tissues as per steps below
- 3 Euthanize mice by CO₂ inhalation or other means of euthanasia.
- 4 Collect the spleen from the mice.

- 5 Put a 40 μ m vacuum filter on a 50mL tube. Mesh the spleen on the strainer using the back of a syringe.
- 6 Add 5mL of RPMI medium to the mashed spleen to collect splenocytes into a 50mL tube.
- 7 Spin down cells at 1500rpm for 3min. Discard supernatant.
- 8 Add 0.5mL ACK lysis buffer to the red cell pellet, resuspend, and wait 30 seconds to lyse RBC.
- 9 Directly add 5mL T cell medium.
- 10 Spin down cells at 1500rpm for 3min. Discard supernatant.
- 11 Prepare a new 50mL tube with a new 40 μ m vacuum filter.
- 12 Resuspend cells in 10mL T cell medium and filter through a new 40 μ m vacuum filter.

13 Count the cells.

Isolation of naïve CD4⁺ T cells from single-cell suspensions from spleen ti...

14 Naïve CD4⁺ T cells were isolated from single-cell suspensions from spleen tissue by negative selection using the EasySep Mouse CD4⁺ T-cell Isolation Kit (#19852, STEMCELL Technologies, Vancouver, Canada).

T cell differentiation on Day 1

15 Seed the cells into plate and culture with different conditions

16 Regulate the cell concentration at 5×10^5 /mL with T cell medium

A	B	C
Naïve CD4 ⁺ T cells	5×10^4 /well	100µL/well
T cell medium		50µL/well
Total	Total	150µL/well

17 Take out the pre-coated plate from 4°C.

18 Remove the supernatant and add 200µL PBS to the pre-coated plate and then remove PBS to wash the coated plate. Repeat 2 times.

19 Seed wells of a 96-well plate with 150µL cells cocktail (100µL cells and 50µL T cell media.

20 Seed 100µL differentiation cocktail.

21 Incubate for 5 days at 37°C.

Detect T cell differentiation by ELISA on Day 5 and Day 6 (Note: The chose...

22 Collect cells and count.

23 Regulate cell concentration at $1 \times 10^6/\text{mL}$.

24 Seed 100µL cells into a new 96-well plate (U bottle)

25 Add 100µL PMA cocktail and culture overnight. For non-stimulate cells, add 100µL T cell medium.

A	B	C
	Final concentration	Amount
T cells	$1 \times 10^6/\text{mL}$	100µL
T cell Medium	n.a	100µL

A	B	C
PMA (50µg/mL)	50ng/mL	0.2µL

- 26 Collect supernatant for ELISA. (Note: The supernatant can be frozen at -20°C.)

Detect Treg cell differentiation by Flow cytometry on Day 5

- 27 Spin down cells at 1500rpm for 3min at RT.

- 28 Remove the supernatant over the sink in one motion.

- 29 Add 150µL FACS buffer and spin down cells at 1500rpm for 3min at RT. Remove the supernatant.

- 30 Surface staining, add 50µL FACS Buffer mAbs cocktail:

A	B
FACS buffer	50µL
Anti-CD4-eFluor 450	0.3µL
Anti-CD25-PerCP	0.3µL

- 31 Incubate at 4°C for 30min in the dark.

- 32 Directly add 150uL FACS buffer and spin down cells at 1500rpm for 3min at RT.
- 33 Remove the staining buffer and thoroughly resuspend cells.
- 34 Add 200μL of Foxp3 Fixation/Permeabilization working solution to each well for 1 hour at RT in the dark. (4×stock solution must be diluted prior to use with the Fixation/Permeabilization Diluent dilute 1 part concentrate with 3 parts diluents, make fresh).
- 35 Centrifuge samples at 400-600g for 5min at RT.
- 36 Add 200μL 1×Permeabilization Buffer (make fresh) to each well and centrifuge samples at 400-600g for 5min at RT. Discard the supernatant.

- 37 Add 50uL 1×Permeabilization Buffer mAb cocktail (Foxp3):

A	B
1×Permeabilization Buffer	50μL
Anti-Foxp3-PE	0.3μL

- 38 Incubate at 4°C for 30min in the dark.

- 39** Directly add 200 μ L 1 \times Permeabilization Buffer to each well and centrifuge samples at 400-600g for 5min at RT. Discard the supernatant.
- 40** Add 200 μ L 1 \times Permeabilization Buffer for Flow Cytometry.