

Nov 24, 2024

Bulk DNA Extraction (Ding Lab)

DOI

dx.doi.org/10.17504/protocols.io.bsnhndb6

Reyka Jayasinghe¹, Li Ding¹, Feng Chen¹, Satok¹

¹Washington University, Saint Louis



Reyka Jayasinghe

Washington University (St. Louis)

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.bsnhndb6**

Protocol Citation: Reyka Jayasinghe, Li Ding, Feng Chen, Satok 2024. Bulk DNA Extraction (Ding Lab). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bsnhndb6>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: February 22, 2021

Last Modified: November 24, 2024

Protocol Integer ID: 47529

Abstract

Bulk DNA Extraction



- 1 Cut tissue (~5 mg) and place in a 1.5 mL microcentrifuge tube. Add 180 μ L Buffer ATL and 20 μ L proteinase K, mix by vortexing and incubate at 56°C until completely lysed (1–3 h).
- 2 Add 2 μ L RNase (100 mg/ml) and incubate for 30 min. at 37°C
- 3 Add 200 μ L Buffer AL. Vortex for 15 sec.
- 4 Incubate at 70°C for 10 min. Briefly centrifuge the tube to remove drops from the lid.
- 5 Add 200 μ L ethanol (96–100%). Vortex for 15 sec. Briefly centrifuge the tube to remove drops from the lid.
- 6 Pipet the mixture onto the QIAamp Mini spin column in a 2 mL collection tube. Centrifuge at 8,000 rpm for 1 min. at room temperature. Discard the flow-through and collection tube.
- 7 Place the QIAamp Mini spin column in a new 2 mL collection tube and add 500 μ L Buffer AW1. Centrifuge at 8,000 rpm for 1 min. at room temperature. Discard the flow-through and collection tube.
- 8 Place the QIAamp Mini spin column in a new 2 mL collection tube and add 500 μ L Buffer AW2. Centrifuge at 14,000 rpm for 3 min. at room temperature. Discard the flow-through and collection tube.
- 9 Place the QIAamp Mini spin column in a new 2 mL collection tube and centrifuge at full speed for 1 min. at room temperature.
- 10 Place the QIAamp Mini spin column in a new 1.5 mL microcentrifuge tube, add 50 μ L DNase free water and incubate at room temperature for 1 min. Centrifuge 8,000 rpm for 1 min. at room temperature.
- 11 Repeat step 10 with a further 50 μ L DNase free water.