



Shipping Paraffin Blocks to the Bodenmiller Lab for IMC Analysis

Marda Jorgensen¹

¹University of Florida

1 Works for me dx.doi.org/10.17504/protocols.io.bf5njq5e

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu

ABSTRACT

This SOP describes the methods used to ship paraffin-embedded blocks to the University of Zurich.

GUIDELINES

Shipments should be tracked via FedEx documentation

MATERIALS TEXT

Shipping box
Bubble wrap or alternative packing supplies
Paraffin blocks
Packing Tape
Fed Ex documentation
Customs Form

SAFETY WARNINGS

Human tissue derived paraffin blocks should be handled wearing gloves

BEFORE STARTING

Before submitting the box for shipping, please contact the Bodenmiller lab by email to confirm that they are prepared to receive samples.

- 1 Obtain a small box, large enough to accommodate a FedEx shipping label.
Wrap each paraffin block separately with protective wrapping.
Include paperwork identifying the block numbers for the blocks being shipped
Stuff the empty space in the box to secure the blocks in place. Seal the box with shipping tape.
- 2 Create a shipping label through the FedEx website. Select Priority Overnight Shipping. Enter the address in the document box below. Print out label.



Stefanie Engler
Winterthurestrasse 190
Bodenmiller Lab, Y11-J80
Zurich, CH8057
41446356630

- 3 In addition to the shipping label, a customs invoice must be created to ship internationally (template attached).

☐ [FedEx Customs Invoice.doc](#)

- 4 Place the customs invoice underneath the shipping label and have the shipping label barcode facing out. Deposit sealed box into a FedEx station.