



FEB 20, 2023

🌐 General procedure for long-term cryopreservation of Yeast

Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

A simple and generalized procedure for freezing yeast cells.

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023. General procedure for long-term cryopreservation of Yeast.

protocols.io

<https://protocols.io/view/general-procedure-for-long-term-cryopreservation-o-cpm6vk9e>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Feb 17, 2023

Last Modified: Feb 20, 2023

PROTOCOL integer ID:
77214

50% Glycerol solution

15m

1 Fill a 200 mL Duran bottle with 100 mL PBS



Materials:

⊗ PBS Gibco, ThermoFisher Catalog #20012019

- 2 Transfer  100 mL 100% Glycerol to Duran bottle

Materials:

 Glycerol **MP Biomedicals Catalog #194680**

- 3 Sterilize by autoclave at  121 °C for  00:15:00



15m



- 4 Incubate a  20 mL batch of Yeast overnight

Note

It is important that the yeast culture is saturated, but not oversaturated. An OD₆₀₀ between 0.4-0.6 is ideal.



- 5 Aliquot saturated culture into an appropriate container. For example:



 60 µL saturated culture in a  200 µL reaction tube

 600 µL saturated culture in a  1.5 mL reaction tube

 900 µL saturated culture in a  2.0 mL cryotube

- 6 Mix in 50% glycerol solution. For example:

 40 µL saturated culture in a  200 µL reaction tube

 400 µL saturated culture in a  1.5 mL reaction tube

 600 µL saturated culture in a  2.0 mL cryotube

- 7 Transfer to an ultra-low temperature or vapor-phase nitrogen freezer for long-term storage

Note

If possible, snap-freeze tubes in liquid nitrogen.