

Sep 10, 2020

# Splitting p0 (6wp) to p1 (T75)

In 1 collection

Andrea Argouarch<sup>1</sup><sup>1</sup>University of California, San Francisco

1 Works for me dx.doi.org/10.17504/protocols.io.8gjhtun

Neurodegeneration Method Development Community

Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)

Andrea Argouarch

University of California, San Francisco

## ABSTRACT

Protocol includes splitting outgrowth of dural cell from a 6 well plate into one T75 flask for expansion.

## DOI

[dx.doi.org/10.17504/protocols.io.8gjhtun](https://dx.doi.org/10.17504/protocols.io.8gjhtun)

## PROTOCOL CITATION

Andrea Argouarch 2020. Splitting p0 (6wp) to p1 (T75). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8gjhtun>

## COLLECTIONS ⓘ

**Dural Cell Isolation and Culturing - Collection**

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 19, 2019

## LAST MODIFIED

Sep 10, 2020

## PROTOCOL INTEGER ID

28907

## PARENT PROTOCOLS

Part of collection

[Dural Cell Isolation and Culturing - Collection](#)

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">Fetal Bovine Serum</a>	97068-091	<a href="#">Vwr</a>
<a href="#">Trypsin-EDTA (0.05%), phenol red</a>	25300062	<a href="#">Thermo Fisher</a>
<a href="#">Dumont #5 Forceps</a>	11251-30	<a href="#">Fine Science Tools</a>
<a href="#">GeneMate Cell Scrapers &amp; Lifter</a>	490000-254	<a href="#">Vwr</a>
<a href="#">DPBS, no calcium, no magnesium</a>	14190250	<a href="#">Thermo Fisher</a>
<a href="#">DMEM, high glucose, pyruvate</a>	11995073	<a href="#">Thermo Fisher</a>

## NAME

Penicillin-Streptomycin

## CATALOG #

15140122

## VENDOR

Gibco - Thermo Fisher

## Observations

- 1 After 3 weeks, split 6 wells of the 6 well plate (6wx6wp) into 1xT75 flask
  - a. Keep 6wp as backup by re-feeding the entire plate after trypsinization.
- 2 Cells will grow on the plastic and maybe on the glass coverslip itself

## Preparation

- 3 Turn off UV lights and clean hood with 70% ethanol
- 4 Clean items with 70% ethanol and bring into hood
  - a. Autoclaved # 5 Forceps



## Dumont #5 Forceps

by Fine Science Tools

Catalog #: 11251-30

- b. 10 cm dish for forceps, cell scraper, discarded coverslips, and discarded tissue
- c. Sterile cell scraper



## GeneMate Cell Scrapers &amp; Lifter

by Vwr

Catalog #: 490000-254

- d. DPBS <sup>-/-</sup>



## DPBS, no calcium, no magnesium

by Thermo Fisher

Catalog #: 14190250

- e. Sterile Filtered Media with PenStrep



## DMEM, high glucose, pyruvate

by Thermo Fisher

Catalog #: 11995073



#### Penicillin-Streptomycin

by Gibco - Thermo Fisher


Catalog #: 15140122



#### Fetal Bovine Serum

by Vwr

Catalog #: 97068-091

- f. Label 1xT75 flask with ID, date, and p1, **add** 15 mls  **15 mL** of media per flask
- g. 0.05% Trypsin








#### Trypsin-EDTA (0.05%), phenol red






by Thermo Fisher


Catalog #: 25300062

- i. Aliquot and warm before use

### Culturing

- 5 Aspirate old media
- 6 Rinse by adding 3 mls  **3 mL** of DPBS per well
- 7 Lift coverslips with forceps to allow adequate washing under the coverslip and gently swirl
- 8 Aspirate DPBS
- 9 Add 2 ml  **2 mL** of trypsin per well and lift coverslips briefly with forceps
- 10 Place in incubator  **37 °C** for ~5 mins  **00:05:00** until cells start to detach
  - a. Can also gently tap the side of the plate to detach cells from the coverslip
- 11 Add 2 mls  **2 mL** of media per well to inactivate trypsin

- 12 With forceps, flip over coverslip within each well
  - a. Can remove tissue and place in 10 cm dish to discard
- 13 Use cell scraper to gently scrap off cells in each well attached to the coverslip
- 14 Discard each coverslip in the 10 cm dish
- 15 Use cell scraper to scrap off remaining cells (cells that are attached to bottom of plate)
- 16 Check under microscope and mark areas with cells still attached with a pen
- 17 Re-scrape cell gently if needed
- 18 Collect cell suspension in a 50 ml conical
- 19 Rinse wells with additional media ~6mls  6 mL to collect any remaining cells
- 20 Spin at 1000 rpm  1000 rpm for 5 mins  00:05:00
- 21 Re-feed old 6wp plate for backup, 2 mls  2 mL per well
- 22 Aspirate supernatant, being careful to not aspirate the pellet
- 23 Tap the pellet to resuspend
- 24 Add media to a final volume of 5 mls  5 mL and pipette up and down to mix

- 25 Add cell suspension to flask
- 26 Place flask in incubator (5% CO<sub>2</sub>, 37°C)  37 °C
- 27 Observe and feed the next day to remove debris, then feed every 2-3 days until 90-100% confluent

#### Clean Up

- 28 Throw away biohazard materials properly
  - a. 10 cm dishes with glass coverslips should be thrown away in biohazard red sharp container
- 29 Clean surgical tools, wear waterproof lab coat and eye protection/PPE
  - a. Brush and clean with 409 soap water. Rinse with water, dry on kimwipe, rinse with 100% ethanol, and then dry completely with kimwipe to prevent rust or water marks.
  - b. Prep for next autoclaving cycle
- 30 Clean and sterilize hood with 70% ethanol and turn on UV
- 31 Update cell culture notes in lab notebook