

NOV 12, 2023

## (EM) and quantification

Leonardo A Parra-

Rivas<sup>1</sup>

<sup>1</sup>University of California, San Diego

Subhojit Roy Lab



Leonardo A Parra-Rivas University of California, San Diego

**ABSTRACT** 

Electron microscopy (EM) and quantification





DOI:

dx.doi.org/10.17504/protocol s.io.yxmvm3zm6l3p/v1

**Protocol Citation:** Leonardo A Parra-Rivas 2023. Electron microscopy (EM) and quantification. protocols.io https://dx.doi.org/10.17504/p rotocols.io.yxmvm3zm6l3p/v

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status: Working** We use this protocol and it's working

Created: Nov 12, 2023

Oct 12 2023

Last Modified: Nov 12,

2023

### **PROTOCOL** integer ID:

90812

## **Electron microscopy (EM)**

1

Hippocampal neurons were plated at a density of  $60,000 \text{ cells/cm}^2$ . DIV-3 neurons were infected with lentiviruses carrying  $\alpha$ -syn tagged at the C-terminus to mScarlet, or mScarlet alone (multiplicity of infection or MOI=2.5). The transduced neurons were cultured to maturity (DIV17-DIV21) before imaging. The  $\sim 100$  % infection efficiency of the mScarlet-tagged constructs was confirmed before using the coverslip for EM experiments.

- 2 EM embedding protocol for monolayer cultured cells was performed on ice. The coverslips were fixed in 2 % glutaraldehyde and in 0.1 M sodium cacodylate (SC buffer) buffer for at least 60min at 4 °C, washed 3 min with 0.1 M SC buffer (5 times), and postfixed in 1% osmium tetroxide (0sO4) diluted in 0.1 M SC buffer for 30 min to 1hr.
- The coverslips were then washed for 3 min (5 times) with 0.1 M SC buffer, followed by a water rinse, sequential ethanol dehydration (20%, 50%, 70%, 90%, 2x100%) 1min each, and subsequent incubation with dry acetone for 1 min at room temperature (2 times).
- The samples were then incubated with a Durcupan:Acetone solution (50:50%) for 1hr in rotation, followed by incubation with 100% Durcupan for 1 hr (2 times). The cells were then embedded using foil plates along with wood sticks in the embedding oven (overnight incubation).
- 60 nm ultrathin sections were cut with a diamond knife and mounted on 300 mesh grids. Grids were post-stained with 2% of uranyl acetate for 5 min and lead for 1 min. JEOL 1400 Electron Microscopy and Gatan digital Camera were used to obtain images.

# quantification

6

For quantification of inter-vesicular (IV) distances between vesicles, we considered all the vesicles in the synapse encompassing  $1\mu mx1\mu m$ . To mark each vesicle, a line was drawn across the center of the vesicle, recording its diameter and its position.

7 These parameters were then fed into our algorithm on MATLAB (RRID:SCR\_001622)( http://www.mathworks.com/products/matlab/), which gave the output of IV distances for a particular image. The IV distance for all images per condition was pooled together. This algorithm is available upon request.