



FEB 10, 2023

## S. O. C. medium

Andreas Sagen<sup>1</sup>

<sup>1</sup>University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

### ABSTRACT

S.O.C. Medium is used in the final step of bacterial cell transformation to obtain maximal transformation efficiency of *E. coli*. S.O.C.

### MATERIALS

LAF bench

Autoclave

Filter unit

Scale

### OPEN ACCESS

**Protocol Citation:** Andreas Sagen 2023. S. O. C. medium. [protocols.io](https://protocols.io/view/s-o-c-medium-cn5svg6e) <https://protocols.io/view/s-o-c-medium-cn5svg6e>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 08, 2023

**Last Modified:** Feb 10, 2023

**PROTOCOL integer ID:**  
76690

**Keywords:** SOB, SOC, Transformation, medium, recovery medium


## 500 mL 1 M Glucose solution

1 Fill the bottle with  400 mL distilled water

2 Measure  90.08 g Glucose

Materials:


 Glucose **Sigma-aldrich Catalog #G7021**

3 Add distilled water to a total of  500 mL

4 Filter sterilize with a 0.2 µm filter


## 250 mL 250 mM Potassium chloride

5 Fill the bottle with  200 mL distilled water

6 Measure  4.66 g Potassium chloride

Materials:

 Potassium chloride **Sigma-aldrich Catalog #P3911**


7 Add distilled water to a total of  250 mL

8 Autoclave solution at  121 °C for  00:15:00

15m


## 125 mL 2 M Magnesium chloride

9 Fill the bottle with  100 mL distilled water

10 Measure  50.825 g Potassium chloride

Materials:

 Magnesium chloride **Sigma-aldrich Catalog #M0250**





11 Add distilled water to a total of  125 mL

12 Autoclave solution at  121 °C for  00:15:00

15m

## 1 000 mL S. O. B. medium

13 Fill the bottle with  950 mL distilled water

14 Measure  20 g Tryptone,  5 g Yeast extract and  0.5 g Sodium chloride. Add  15 g Agar to create agar plates





Materials:


 Tryptone **Millipore Catalog #T9410**

 Yeast Extract **Sigma-aldrich Catalog #Y0875**

 Sodium chloride **Sigma-aldrich Catalog #S9625**


 Agar **Sigma-aldrich Catalog #A1296**

15 Add  1 mL Potassium chloride solution (  250 millimolar (mM) ) and  0.5 mL Magnesium chloride solution (  2 Molarity (M) )

16 Adjust pH to  7.0 with concentrated sodium hydroxide



Materials:

 Sodium hydroxide solution **Sigma-aldrich Catalog #S8263**

17 Add distilled water to a total of  1000 mL

18 Autoclave solution at  121 °C for  00:00:00

### 100 mL S. O. C. medium

19 Sterilely transfer  98 mL S. O. B. medium to a sterile  100 mL container

20 Sterilely add  2 mL  1 Molarity (M) Glucose solution