



Feb 12, 2021

# *Creolimax fragrantissima*\_culture method

Meritxell Antó<sup>1</sup><sup>1</sup>Institut de Biologia Evolutiva (CSIC-Universitat Pompeu Fabra)**1** Works for me This protocol is published without a DOI.

Meritxell Antó

SUBMIT TO PLOS ONE

## PROTOCOL CITATION

Meritxell Antó 2021. *Creolimax fragrantissima*\_culture method. **protocols.io**  
<https://protocols.io/view/creolimax-fragrantissima-culture-method-bsdjna4n>

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Feb 12, 2021

## LAST MODIFIED

Feb 12, 2021

## PROTOCOL INTEGER ID

47243

## Growing medium

- 1 BD Difco Marine Broth 2216 • Marine Agar 2216  
[https://legacy.bd.com/europe/regulatory/Assets/IFU/Difco\\_BBL/212185.pdf](https://legacy.bd.com/europe/regulatory/Assets/IFU/Difco_BBL/212185.pdf)

## Growth conditions

- 2 Temperature 12°C

## Cryopreservation and recovery (from liquid cultures)

- 3 Cryopreservation (10%DMSO in Growing medium)

1. Pre-grow the culture in appropriate growing medium until almost saturated but not too old.
2. Mix 1.35mL of the culture with 150 ul of DMSO in the screw-top cryovials.
3. Transfer the cryovials to the Freezing container, Nalgene® Mr. Frosty unit. Place it in the -80°C freezer and leave it overnight.
4. Remove cryovials from the Freezing container and store in appropriate boxes for -80°C.

## Recovery

1. Fill a flask with 10mL fresh Growing medium.
2. Take out the frozen stock from the freezer, carry it to the culture hood in a cryo box to prevent it from thawing.
3. With a pipette tip (or a loop), take a chunk of the frozen stock and resuspend it in the culture flask (~10uL is probably sufficient). Take the frozen stock back to the -80°C as rapidly as possible.
4. Let it recover for at least 3-4 days and proceed with another transfer in a fresh medium even if the cells are not completely grown it will help to dilute the remaining DMSO from the cryopreservation, the time will vary depending on the amount of the initial inoculum and the culture type.

\*NOTE: because of DMSO some precipitates may appear in the culture.

All protists must be stored in liquid nitrogen for long-term storage. Although storage at -80°C is sometimes sufficient for a temporary alternative.