



# Preparation of staphylococcal protein-A conjugated to horseradish peroxidase by the periodate method.

Angel A Justiz-Vaillant<sup>1</sup>

<sup>1</sup>University of the West Indies St. Augustine

1 Works for me [dx.doi.org/10.17504/protocols.io.bjk5kky6](https://dx.doi.org/10.17504/protocols.io.bjk5kky6)

University of the West Indies [angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)

Angel Justiz-Vaillant  
University of the West Indies St. Augustine

## ABSTRACT

This reagent can be used in ELISA, Western blotting and Dot blot to detect antigens and antibodies. It is important in the immunodiagnosis of infectious diseases and other problems. I find this useful in the detection of anti-HIV antibodies by ELISA.

## DOI

[dx.doi.org/10.17504/protocols.io.bjk5kky6](https://dx.doi.org/10.17504/protocols.io.bjk5kky6)

## PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Preparation of staphylococcal protein-A conjugated to horseradish peroxidase by the periodate method.. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjk5kky6>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 11, 2020

## LAST MODIFIED

Aug 11, 2020

## PROTOCOL INTEGER ID

40317

## GUIDELINES

All reagents but specially the enzyme and sodium periodate solution has to be prepared freshly before mixing it with the enzyme.

## MATERIALS

NAME	CATALOG #	VENDOR
Ammonium Sulfate		P212121
Sodium periodate	SB0875.SIZE.100g	Bio Basic Inc.
sodium borohydride	452882	Sigma Aldrich
Horseradish Peroxidase (HRP) type IV	P8375-25KU	Sigma Aldrich
Staphylococcal Protein-A		Sigma Aldrich

## MATERIALS TEXT

Pipettes

20ml to 1000 ml glass  
Scale  
Incubator  
Refrigerator  
Freezer  
Centrifuges

- 1 Horseradish peroxidase (500 µg in 50 µl NaCO<sub>3</sub>, pH 9.6) is mixed with freshly made sodium periodate solution (1.71 mg/ml) followed by incubation in the dark for 2 h.
- 2 Mix 500 µg of staphylococcal protein-A (SpA) with an equal amount (500 micrograms) of a mix of horseradish peroxidase-sodium periodate.
- 3 The mixture is incubated for 3 hours at 4°C with gentle agitation.
- 4 Forty µl of freshly prepared NaBH<sub>4</sub> solution (5 mg NaBH<sub>4</sub> /ml 0.1 mM NaOH) is then added to the preparation.
- 5 The preparation is incubated for 90 min at 4°C in the dark with gentle agitation.
- 6 Cold 50% saturated ammonium sulphate solution (pH 7.4) is added drop by drop in the ratio 1:1 (v/v).
- 7 The mixture is then centrifuged for 25 min at 4°C and recover the pellet at the bottom of the tube.
- 8 The pellets is re-suspended in 200 µl of PBS pH=7.4 and dialysed against 1L of PBS for 24 h with 3 buffer changes.
- 9 An equal volume of glycerol is added to the dialysate followed by 100 µl of bovine serum albumin, BSA (20 mg/ ml).
- 10 The conjugate is then stored at -20°C until further used.

11