



Jan 02, 2021

# Growth curves

Elizabeth Fozo<sup>1</sup><sup>1</sup>In-house protocol

1

Works for me

This protocol is published without a DOI.



Eadewunm

## ABSTRACT

Rules of thumb and directions for making growth curves

## PROTOCOL CITATION

Elizabeth Fozo 2021. Growth curves. **protocols.io**  
<https://protocols.io/view/growth-curves-bq49myz6>



## KEYWORDS

growth curves

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jan 02, 2021

## LAST MODIFIED

Jan 02, 2021

## PROTOCOL INTEGER ID

45953

## GUIDELINES

- Rules of thumbs
- Making an overnight culture
- In the morning
- Growth curve

## DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## ABSTRACT

Rules of thumb and directions for making growth curves

Rules of thumb

- 1 Never remove > 1/3 of a culture for a growth curve
- 2 Read OD<sub>600</sub> (optical density at 600 nm) of 200 microliters of culture added to 800 microliters of fresh media → 1,000 microliters of fluid in cuvette; = 1:5 dilution
  - Remember to set the factor equal to 5 when reading from the spec
- 3 Easier to measure 800 microliters of fresh media into cuvettes ahead of time; simply store in cuvette box
- 4 Easier to measure media into a 50 mL conical for filling fresh cuvettes; less likely to get a contaminated bottle
- 5 Easier to use same blank multiple times; suggest every time you refill conical with media for cuvettes, make a new blank with that and use with cuvettes from that mixture.
- 6 Place cuvettes into 10% bleach solution when finished \*This is a must for BSL 2 organisms such as *E. faecalis*.

#### Making an overnight culture

- 7 Measure 5mL appropriate media with appropriate antibiotics into fresh conical (make sure to label conical with date, strain!)
- 8 Take 1 colony from a plate with appropriate strain (use pipette to poke that colony)
- 9 Inoculate media
- 10 Grow overnight (usually start culture ~3-5 p.m. night before culture is needed)

#### In the morning:

- 11 Measure OD<sub>600</sub>
- 12 Use  $V_1C_1 = V_2C_2$  to make a 0.01 OD<sub>600</sub> dilution in media for growth curve
- 13 Place culture in flasks, mix, place in appropriate incubator, start count up timer.

Growth curve:

- 14 The first time point is at 30 min on count up timer.
- 15 Take time point every 30 min, write down data and graph as you go
- 16 Use the count-up timer to measure time. When it reaches 30 min, 1 hour, 1 ½ hours, etc. remove cultures from incubator, pipette into cuvettes (remember 1/5 dilution, e.g. 200 microliters culture in 800 microliters media)
- 17 Remember to mix culture gently before adding to cuvette (the culture can be clumpy at the bottom, and you want a homogenous solution)