



Sep 09, 2022

Immunocytochemistry (ICC)

mineechoi^{1,2}¹Queen Square UCL Institute of Neurology; ²The Francis Crick Institute

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.q26g74w79gwz/v1

Minee-Liane Choi

UCL Institute of Neurology, The Francis Crick Institute

ABSTRACT

This protocol describes how to do immunocytochemistry for primary and hiPSC-derived cells.

DOI

dx.doi.org/10.17504/protocols.io.q26g74w79gwz/v1

PROTOCOL CITATION

mineechoi 2022. Immunocytochemistry (ICC). **protocols.io**
<https://protocols.io/view/immunocytochemistry-icc-b9yfr7tn>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 25, 2022

LAST MODIFIED

Sep 09, 2022

PROTOCOL INTEGER ID

63207

Cell fixation

- 1 Cells are fixed in **4 % volume** paraformaldehyde (PFA) and stored in phosphate-buffered saline (PBS) until use.

Permeabilizing and blocking cells

- 2 Wash cells with PBS twice.
- 3 Incubate the cells in **0.2 % volume** Triton X-100, **5 % volume** bovine serum albumin^{1h} (BSA) for **01:00:00** at **Room temperature**.

5 % volume BSA (made in PBS) is used to block non-specific binding.

For ATTO 425 labelled Aptamer staining, cells are permeabilized with **0.25 % volume** Triton X-100 and blocked with **10 % volume** normal goat serum (NGS) for **00:20:00** followed by another **03:00:00** with **0.1 % volume** Triton X-100 and **10 % volume** NGS.

Incubate cells in primary antibodies

- 4 2h

Do not wash after permeabilising and blocking steps

Dilute primary antibody in **5 % volume** BSA and incubate cells at **4 °C**, **Overnight** or **01:00:00** at **Room temperature**.

The final volume should be sufficient to cover each coverslip around **170 µL** for 8-ibid

chambers, #80806). For 8-ibid, it recommends incubating the cells at room temperature for 🕒 **01:00:00** at 🌡 **Room temperature** if possible.

- 5 Wash cells with [M] **5 % volume** BSA for 🕒 **00:05:00** three times. 5m

Incubate cells in secondary antibodies 1h 5m

- 6 Dilute primary antibody in [M] **5 % volume** BSA and incubate cells at 🌡 **4 °C** 🕒 **Overnight** or 1h 5m
🕒 **01:00:00** at 🌡 **Room temperature** .

- 7 Wash cells with [M] **5 % volume** BSA for 🕒 **00:05:00** three times away from light. 5m

Add Hoechst ([M] **10 micromolar (μM)**) in the second wash and leave for 🕒 **00:15:00** .

- 8 Take away PBS and load anti-fading medium to cover cells.

For the short term, imaging in PBS is also fine.