


Aug 05, 2024

Endocytosis and internalization assay in primary neuronal culture



DOI

dx.doi.org/10.17504/protocols.io.4r3l22j5ql1y/v1

Arpine Sokratian¹, andrew.west west²

¹Duke Univeristy; ²Duke University

ASAP Collaborative Rese...

West lab protocols



Arpine Sokratian

Duke Univeristy

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4r3l22j5ql1y/v1

Protocol Citation: Arpine Sokratian, andrew.west west 2024. Endocytosis and internalization assay in primary neuronal culture. protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l22j5ql1y/v1>

Manuscript citation:

DOI: [10.1126/sciadv.adi8716](https://doi.org/10.1126/sciadv.adi8716)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2024

Last Modified: August 05, 2024

Protocol Integer ID: 93919

Keywords: ASAPCRN

Funders Acknowledgement:

ASAP

Grant ID: ASAP-020527

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

This protocol outlines a method for quantifying endocytosis in primary neuronal cultures, specifically focusing on different a-synuclein fibril structures. By utilizing pH-sensitive dye conjugated to experimental proteins/fibrils, we can identify endocytosis pathways and analyze intensity changes in the pHrodo-specific channel from real-time acquired images at various time points.

Materials

 EIPA **caymanchem Catalog #14406**

 Dyngo **Abcam Catalog #ab120689**


 Methyl-β-cyclodextrin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4555**


 Pitstop® 2 **Abcam Catalog #ab120687**

 Wortmannin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W1628-1MG**

Protocol materials

 Hoechst 33342 **Catalog #H3570** Step 5

 LysoTracker™ Green DND-26 - Special Packaging **Thermo Fisher Catalog #L7526** Step 5

 pHrodo™ Red Dextran, 10,000 MW, for Endocytosis **Thermo Fisher Catalog #P10361** Step 5

 EIPA **caymanchem Catalog #14406** Materials

 Dyngo **Abcam Catalog #ab120689** Materials

 Methyl-β-cyclodextrin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4555** Materials

 Pitstop® 2 **Abcam Catalog #ab120687** Materials

 Wortmannin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W1628-1MG** Materials

Preparation of endocytosis inhibitors

- 1 For the assay we use DIV7 primary hippocampal neuron culture plated in 48-well plates
- 2 Prepare stock solution of the endocytosis inhibitors in DMSO:

A	B
2.3 mM	Wortmannin
50 mM	Dyngo
10 mM	Pitstop 2
384 mM	methyl- β -cyclodextrin (M β CD)
50 mM	ethyl-isopropylamiloride

Recommended concentrations for stock solutions (stocks can be freezed at -80C for couple of days)

- 3 At the day of the experiment, thaw down the stock solutions in a water bath and dilute the drugs to reach concentrations:

	A	B	C	D	E	F	G
	Ethyl-isopropyl amiloride (EIPA)	50 uM	500uM: 10 ul of stock + 990ul PBS				
	Dyngo	10 uM	10 mM: 20ul of the drug + 80 ul (50% DMSO)	for 100 uM of Dyngo: 5 ul of the 10mM dyngo + 495 ul of PBS			
	Wortmannin	0.2 uM	10 uM: 4ul of the stock + 916 ul PBS	1uM: 200 ul of 10 uM + 800 ul of PBS			
	Pitstop 2	15 uM	150 uM: 15ul of P + 984 ul PBS				
	Methyl-β-cyclodextrin (MβCD)	2 mM	20 mM: 25 ul of stock + 475 ul PBS				

Drugs are prepared to be diluted to final concentration in the cell culture




4 Before adding protein/drug of interest, add endocytosis inhibitors 30 minutes before the treatment. As calculated for 48-well plates with 3 mL of the media, 30 ul of diluted drugs are sufficient. ** Note: Be careful adding the diluted drugs (should be RT). Place the plate back to the cell culture incubator.

5 After 30 minutes, add the protein/drug of interest supplement with


 Hoechst 33342 **Contributed by users Catalog #H3570**

As for pHrodo-conjugated fibrils, please, see the protocol describing the conjugation process (step 5):

As a control for endocytosis add lysotracker

 LysoTracker[®]; Green DND-26 - Special Packaging **Thermo Fisher Catalog #L7526**

and phrodo-10kDA dextran

 pHrodo[®]; Red Dextran, 10,000 MW, for Endocytosis **Thermo Fisher Catalog #P10361**

Protocol

NAME



Preparation of mouse and human α -synuclein fibrils

CREATED BY

Arpine Sokratian

PREVIEW

6 After two hours of incubation, start taking images of each well at the high-speed (high sensitive mode)

Continue acquiring images every two to four hours up to 48 hours of incubation.

7 Calculate the images using cell count (Hoechst) for normalization and calculate the signal for each time-point as a proportion to the maxima of the signal

Calculate the signal from LysoTracker after 30 minutes of incubation and compare the control (no endocytosis inhibitors) to the experimental reactions.

