



Oct 16, 2020

Bacterial Glycerol Stocks for Long-term Storage

Jiaxin Li¹¹South China University of Technology**1** Works for me dx.doi.org/10.17504/protocols.io.bm8xk9xn

Jiaxin Li

DOI

dx.doi.org/10.17504/protocols.io.bm8xk9xn

PROTOCOL CITATION

Jiaxin Li 2020. Bacterial Glycerol Stocks for Long-term Storage. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bm8xk9xn>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 10, 2020

LAST MODIFIED

Oct 16, 2020

PROTOCOL INTEGER ID

42999

- 1 Add 750 μ L of the overnight culture to 250 μ L of 60% glycerol in a 2 mL cryotube and gently mix.
- 2 Freeze the glycerol stock tube at -80°C.