



Feb 18, 2021

Prestwick screen protocol

Forked from [Syngenta divergent strain screen protocol](#)Ida Barlow¹¹Imperial College London**1** Works for me This protocol is published without a DOI.

Behavioural Genomics

Ida Barlow

SUBMIT TO PLOS ONE

ABSTRACT

Protocol for screening the Prestwick *C. elegans* drug library of 240 drugs at 3 concentrations and imaging under baseline and bluelight conditions using the Hydra (Loopbio) imaging rigs. All drugs at all concentrations were imaged with N2 in a single day of tracking. This protocol was repeated 3 times so that there were 12 replicates per drug per concentration.

PROTOCOL CITATION

Ida Barlow 2021. Prestwick screen protocol. **protocols.io**
<https://protocols.io/view/prestwick-screen-protocol-bsibncan>



FORK NOTE

FORK FROM

Forked from [Syngenta divergent strain screen protocol](#), Ida Barlow

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 18, 2021

LAST MODIFIED

Feb 18, 2021

PROTOCOL INTEGER ID

47395

GUIDELINES

Careful planning of how drugs to be arranged in plates and the number of strains is required before undertaking screening experiments. Using a google calendar to pre-plan timings and days is advised in order to efficiently manage the workload

MATERIALS TEXT

VIAFLO

96 channel pipette

Integra

VIAFLO 96



VIAFILL

reagent dispenser

Integra

VIAFILL



BEFORE STARTING

96 well plates containing the drug library at all doses need to be randomised by column using the OpenTrons robot to create 3 stock plates for each library plate.

Pick L4 worms for bleaching (-9 days from first day of tracking) 1h

- 1 Pick 10 x N2 L4s onto 8 x 90mm plates (pre-seeded with OP50) for each strain

Pour 96WPs (up to 3 days before tracking) 4h

- 2 Prepare 1L low peptone NGM and autoclave



Making low peptone NGM for imaging plates

by Ida Barlow

PREVIEW

RUN



- 2.1
 - Book the autoclave (notebook on top of the machine).
 - Take clean flasks from the glass kitchen.
 - Measure all the pre-autoclave reagents and add to the flask (Use a new weighing boat and spatula for each reagent. Also, the cholesterol is kept in the fridge.)
 - Once water is added mix thoroughly and label with autoclave tape ('Low peptone NGM Rm 5020'). Make sure the bottle is not screwed completely when placing it inside the autoclave machine.

- 2.2
 - Turn ON the autoclave
 - Make sure that the autoclave's probe bottle is the same size as the largest bottle you use and fill it with water.
 - Place the temperature probe in it.
 - Fill up the autoclave with water until it reaches the grill.
 - Place the bottles in the autoclave and make sure that the cap is not screwed completely.
 - Check the waste flask is not too full
 - Use 'media' program.
 - Press START.
 - It will take about 2 hours for a 500ml bottle to autoclave and about 2.5 hours for 1L or larger bottles.

- When autoclave is complete, remove the probe flask

- 2.3
 - Make sure to wear gloves as the flask will be hot
 - Let the agar to cool to around 55°C, ie the bottle is cool enough to hold for a second with a gloved hand.
 - Add the post autoclave reagents.
 - Mix it well and start pouring onto imaging plates (See Protocol for plate pouring)
 - Try not to shake the bottle too much while mixing to avoid air bubbles.
 - The agar needs to be warm to be poured without blocking the tubings, so try to pour as quickly as possible and if not poured immediately put the bottle on a waterbath set to 60°C until being used.
- 3 Once agar has cooled to around **65 °C**, add the salts and dispense agar into 40 x square well 96 well plates using VIAFILL dispenser. Dispense **200 µl** per well. Once cooled, store agar side up in an airtight container at **4 °C**

Dry plates, bleach worms (-5 days)

4h

- 4 Dry 15 x 150mm plates in cabinet dryer (setting 2) or 3 hours
- 5 Seed 15 x 150mm plates with OP50 and leave to dry overnight at room temperature
- 6 Bleach worms prepared for day 1 of tracking



Bleach synchronisation of *C. elegans*
by Ida Barlow

PREVIEW

RUN



- 6.1 Wash hermaphrodites off plate with several ml of M9 solution and transfer to 15ml falcon tube (Fisher Scientific-Falcon 352096)
- 6.2 Fill falcon tube up to 15ml with M9 solution
- 6.3 Centrifuge for 2 minutes at 1500 rpm (RCF:210, ascending 9; descending 7) – program 1

Program 1 retains the worms as pellets and the bacteria is suspended as the supernatant

The descending is slow as the worm pellet is loose at this stage which we don't want to break
- 6.4 Remove supernatant using a plastic Pasteur pipette taking care not to disturb pellet
Leave at least 0.5ml M9 to avoid disturbing the pellet
- 6.5 Fill the tube with M9 up to 15ml
- 6.6 Spin program 1

6.7 Repeat steps 4-6

6.8 On final wash remove as much supernatant as possible and add M9 upto 4ml

6.9 Add 4ml 2X Bleach solution (From here onwards try to work as quickly as possible to avoid over-exposure of the worms to the bleach)

USE FRESHLY PREPARED BLEACH EVERYTIME

2X Bleach solution:

5% Sodium hypochlorite solution - 4ml
Sterile water - 3.5 ml
1M NaOH solution - 2.5 ml
TOTAL - 10 ml

6.10 Vortex on maximum setting for 4 min (no more as this will damage the eggs)

Makesure the vortex forms

After vortexing, top up the tube with M9 till 15ml

6.11 Centrifuge for 2 mins at 2500rpm (RCF:590, ascending 9; descending 7) – program 2

(Always check the program on the centrifuge before using it)

6.12 Remove supernatant by pouring into waste bottle – pellet should be compact and yellow in colour at bottom of falcon, but be careful not to lose

6.13 Add 15ml M9

6.14 Centrifuge at program 2

6.15 Repeat steps 12-14 four more times

The number of washes is crucial here as we need to get rid of all the bleach

6.16 After final wash add 15ml M9 and store eggs/larvae in the falcon on the rotator that is constantly spinning at 20°C, until feeding

L1 arrested larvae can be starved for up to 5 days before refeeding

6.17 Centrifuge larvae on program 2 to pellet

6.18 Remove supernatant with plastic Pasteur pipette

The pellet is loose here so make sure not to disturb it

6.19 Add 15ml M9, spin to wash

6.20 On final wash leave 0.5ml M9 in falcon

6.21 Resuspend the pellet by gently tapping the tube/flicking it

6.22 Place droplet containing larvae onto seeded plate and allow to grow to desired developmental state (ie. 2 days for L4s, 2.5 days for young adults)

Use glass pipette to place the droplet onto seeded plate, avoid using plastic pipette as larvae will stick to it

Development times at 20°C:

- 2 days for L4s
- 2.5 days for young adults

Note:

- If you feed larvae within 12hrs of bleaching then they develop faster than the longer arrested ones
- It is a good practice to bleach in two tubes in parallel
- If you drop the tube at any point of the process, make sure to transfer the contents into a new tube as the dropped tube may get cracked resulting in loss of worms during centrifugation/vortexing
- Any unused larvae can be topped up with M9 and stored spinning in the rotator to be re-used
- Use clean autoclaved rubber bulbs for the refeeding everytime to avoid contamination
- Put the used bulb in the box labelled 'Used Teets'

Stages	Grown at 20 C from L1	Grown at 25 C from L1
L1 division	11.7hrs	9hrs
Mid L1	16.9hrs	13hrs
First L2 division	22.1hrs	17hrs
Between L2 divisions	23.4hrs	18hrs
Second L2 divisions	24.3hrs	19hrs
Mid L2	29.9hrs	23hrs
L3 division	32.5hrs	25hrs
Mid L3	37.7hrs	29hrs
L4 division	42.9hrs	33hrs
Mid L4	49.4hrs	38hrs
Early adult	55.9hrs	43hrs
Adult	62.4hrs	48hrs

Table of Development times for different temperatures

Refeed L1s (-3 days)

4h

- 7 At 11:00, spin L1s for day 1 of tracking at 2500rpm. Remove supernatant and using glass pipette, drop 4 small droplets around the edges of the plate (off food) onto 14 x 150mm plates

Allow to grow at \uparrow 20 °C

Dry plates and dispense drugs onto imaging plates using VIAFLO (-1 days from imaging)


3h

- 8 In the morning, take 40 poured 96 well plates from the cold room, and weigh three random plates without their lids
- 9 Place in cabinet dryer (setting 1.5-2) and allow to dry for 2-3 hours with lids off
- 10 Weigh 3 random plates and verify that at least 3-5% reduction in weight
- 11 Pre-label dried imaging plates (square well) with the imaging run and drug plate information, so that every plate on a single day of imaging has a unique plate id, for example P01_sh01_01 where:


P01-09 - library plate number (out of 4 library plates)
sh01-04 - shuffle number (out of the 4 shuffled stock plates)
01-07 - imaging run number for that day
- 12 Remove the shuffled library plates (stock) plates from the \uparrow -20 °C freezer, allow to thaw at room temperature and spin to collect contents at bottom of wells
- 13 Pre-label the skirted 96 well plates with library plate IDs (P01-L09_sh01-sh04) to make up diluted drug plates
- 14 Dispense \square 7 μ l water into each well of the prelabeled dilution plates using multichannel pipette and reagent reservoir
- 15 Using VIAFLO (hedgehog) dispenser on BG_STOCK custom program, premix drug in drug library plates in slot A, and then transfer \square 0.5 μ l drug in DMSO to the dilution plates prefilled with water in slot B. Repeat for all drug library plates

Double check the dispense volumes before making up the diluted plates

Z-heights have been configured in this program to prevent pipette and plate crashes

- 16 Using VIAFILL (octopus) dispenser, dispense  5 µl water onto 5 x predried imaging plates

5 plates at a time to prevent the agar absorbing all the liquid before the drug is dispensed into the water droplet

- 17 Using VIAFLO in custom program BG_AGAR, with correct drug library plate in slot B, transfer  3 µl of diluted drug and water mixture onto the correctly labelled imaging plate in slot A.


Repeat until all imaging plates have had drug dispensed onto them.

Z-heights have been configured in this program to prevent pipette tips from piercing the agar

- 18 Prepare 1:10 dilution of OP50 in M9 in a small bottle:

 5 mL OP50

 45 mL M9

- 19 Using VIAFILL dispenser, seed all the imaging plates with  5 µl per well

- 20 Place lids on each plate and leave drugged and seeded imaging plates overnight at room temperature in the dark (with box on top)

Imaging

- 21 Wash worms off half (7) of the 150mm plates with M9 buffer using pasteur pipette into 15ml falcons

- 22 Spin at 1500rpm for 2 minutes to pellet the worms

- 23 Remove supernatant and fill with M9

- 24 Repeat steps 28-29

- 25 After final wash, fill falcon with M9 and transfer contents from 15ml falcon to 3 x 50ml and fill up 30ml with M9

- 26 Use COPAS wormsorter to dispense 3 worms per well (pure, no double) into 5 imaging plates at a time. Use pre-made YYYYMMDD_wormsorter.csv to plan and determine which plates to dispense onto
- 27 Allow liquid to dry off by placing imaging plates in **20 °C** incubator with lid off for 30 minutes, then replace lids and keep in **20 °C** incubator
- 28 Expose worms to drug for 4 hours in total as calculated from the middle wormsorter time and allow worms to acclimate for 30 minutes in the cave prior to imaging

Example:
wormsorter start time 10 :00
wormsorter end time 10:40
middle wormsorter time 10:20

cave time 11:50
Imaging start time 12:20

- 29 Imaging on hydra using protocol script (5 mins prestim; 6 mins bluelight with 60 sec OFF, [10sec ON, 90sec OFF] x 3 times; 5 mins postsim)
- 30 After 15 plates have been filled with this batch of N2s, repeat steps 21-25 with the other batch of N2s on 150mm plates and record which batch was used in each plate in the wormsorter.csv file.