



May 15, 2022

# DNA extraction from dermatophytes using the Qiagen DNEasy™ UltraClean Microbial kit (REF: 12224-50)

Khalid El Moussaoui<sup>1</sup><sup>1</sup>Université de Liège[dx.doi.org/10.17504/protocols.io.q26g74d79gwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g74d79gwz/v1)Khalid El Moussaoui  
Université de Liège

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

This protocol describes the steps necessary to extract and purify genomic DNA from dermatophytes (and more specifically from dermatophytes of the genus *Trichophyton*).

DOI

[dx.doi.org/10.17504/protocols.io.q26g74d79gwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g74d79gwz/v1)

Khalid El Moussaoui 2022. DNA extraction from dermatophytes using the Qiagen DNEasy™ UltraClean Microbial kit (REF: 12224-50). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.q26g74d79gwz/v1>



---

 protocol ,

May 15, 2022

May 15, 2022

62631

Perform the DNA dosage directly after extraction and not after a freeze/thaw cycle. Be careful when preparing the medium : work under sterile conditions as much as possible to avoid contaminating the liquid medium.

Qiagen DNEasy UltraClean Microbial kit (REF: 12224-50)

Sabouraud Dextrose Broth (REF : Merck S3306)

:



#### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

### Medium preparation 35m








- 1 Dissolve **30 g** of **Sabouraud dextrose broth Merck Millipore Sigma Catalog #S3306** in **1 L** of **MilliQ water Contributed by users** and let mix on the heated magnetic stirrer for **00:05:00** (temperature and mixing speed knob at mid-step). 5m
- 2 Cover the flask with glass wool and aluminium foil. Autoclave it at **121 °C** for **00:30:00**. 30m

### Cultivation of the strains

- 3 After allowing to cool, transfer  **25 mL** of this medium into a tube. Label the tube with the strain number.
- 4 Using a sterile swab (or a sterile inoculation loop), gently collect the primary culture and dip the swab (or the sterile inoculation loop) into the tube containing the culture medium (prepared in the previous step). Close the tube halfway to allow gas exchange.
- 5 Allow to grow in the incubator at  **30 °C** until a sufficient flocculate is formed (requires at least 96 hours). Incubation time varies from strain to strain but flocculate should be visible after 5 days. If this is not the case, repeat the cultivation step.









#### Preliminary steps




3m

- 6 Using a Pasteur pipette, carefully remove the flocculate from the tube containing the previously cultured dermatophyte strain. Transfer this flocculate to a PowerBead tube containing glass beads, let's call it primary tube. Add  **300 µL** of PowerBead solution and  **50 µL** of SL solution to this tube.
- 7 Cool this tube to  **-196 °C** in liquid nitrogen for  **00:01:00**. Then, heat this tube in a water bath at  **56 °C** for  **00:01:00**. Finally, run this tube through the cell disruptor at maximum speed for  **00:01:00**. This constitutes 1 cycle of 3 steps. You must repeat this cycle 3 times. The recovered mixture is referred to as primary lysate in the following steps.

#### DNA extraction














6m 50s

- 8 Centrifuge the tube at  **10000 x g** for  **00:00:30** at room temperature. Gently transfer the supernatant to a clean collection tube (provided in the kit) and discard the PowerBead tube. 30s
- 9 Add  **100 µL** of IRS solution to the supernatant and vortex for  **00:00:10**. Incubate at  **4 °C** for  **00:05:00**. 5m 10s
- 10 After that, centrifuge the tube at  **10000 x g** for  **00:01:00** at room temperature. Gently transfer the supernatant to a new collection tube and discard the tube containing the pellet. 1m


- 11 Add  **900 µL** of SB solution to the tube containing the supernatant from the previous step<sup>10s</sup> and vortex for  **00:00:10**. Load approximately  **700 µL** of this suspension into a silica membrane chromatography column (provided in the kit).

#### DNA purification

4m

- 12 Centrifuge the column at  **10000 x g** for  **00:00:30** at room temperature. Keep the column and discard the flows-through. Repeat until the entire volume from step 11 is loaded into the column. <sup>30s</sup>
- 13 Add  **300 µL** of CB solution into the column and centrifuge it at  **10000 x g** for  **00:00:30** at room temperature. Keep the column and discard the flows-through. <sup>30s</sup>
- 14 Centrifuge the column alone (empty) to remove the last residues of CB solution. The conditions are identical to the previous step :  **10000 x g** for  **00:00:30** at room temperature. <sup>30s</sup>
- 15 Place the column in a new collection tube. Add  **50 µL** of EB solution to the center of the silica membrane. Let stand for  **00:01:00** at room temperature. Let stand for  **00:01:00** at room temperature and then centrifuge the column at  **10000 x g** for  **00:00:30** to elute the DNA. <sup>2m 30s</sup>
- 16 Discard the column and keep the flows-through which is the purified DNA. Store DNA at  **-80 °C** to ensure stability.

#### Spectrophotometric dosage

- 17 To determine the purity and concentration of the DNA, a NanoDrop dosage was performed. For this purpose, a negative control was prepared beforehand. This control will have undergone all the extraction steps but will not contain any material from dermatophytes.
- 18 Launch the computer program and select the "nucleic acid" mode. Make sure the sample deposit spot is clean and dry. If necessary, clean it with the wipes provided for this purpose. Then drop  **2 µL** of the negative control and click on the "blank" box.
- 19 Proceed in the same way to measure the sample containing the DNA, but click on "measure" instead of "blank". There is no need to redo a blank between measurements.

