



Version 1

Nov 23, 2020

DNA extraction for human microbe samples. V.1

In 3 collections

Lilan Hao¹¹BGI-Shenzhen, Shenzhen 518083, China.

1

Works for me

dx.doi.org/10.17504/protocols.io.bcmriu56

Lilan Hao

ABSTRACT

This protocol is used to clarify the process of total DNA extraction for human microbe samples.

DOI

dx.doi.org/10.17504/protocols.io.bcmriu56

PROTOCOL CITATION

Lilan Hao 2020. DNA extraction for human microbe samples.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bcmriu56>

COLLECTIONS ⓘ

The Female Urinary Microbiota Protocols Collection **Protocols for the female urinary microbiota in relation to the reproductive tract microbiota** **Protocols for "The female urinary microbiota in relation to the reproductive tract microbiota."**

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 17, 2020

LAST MODIFIED

Nov 23, 2020

PROTOCOL INTEGER ID

33169

PARENT PROTOCOLS

Part of collection

[The Female Urinary Microbiota Protocols Collection](#)[Protocols for the female urinary microbiota in relation to the reproductive tract microbiota](#)[Protocols for "The female urinary microbiota in relation to the reproductive tract microbiota."](#)

ABSTRACT

This protocol is used to clarify the process of total DNA extraction for human microbe samples.

- 1 200µl samples were mixed with 500 µl TE buffer and 20 µl of 10 mg/ml lysozyme solution.

- 2 Incubate at 37 °C for 90min. Invert the tube several times to mix every 15 minutes. Centrifuge at 13,000 rcf/min for 10 min at 4 °C and discard the supernatant.
- 3 Add 300 µl of 20 mg/ml proteinase K and 5µl of 10% SDS. Vortex to resuspend cells.
- 4 Incubate at 55 °C in a shaking water bath for 120 min to effect complete lysis. Invert the tube several times to mix every 15 minutes.
- 5 Centrifuge at 13,000 rcf/min for 10 min at 4 °C, then transfer the supernatant to a nuclease-free 2 ml microfuge tube.
- 6 Add 700 µl phenol-chloroform-isoamylalcohol (25:24:1) and vortex to mix well, centrifuge at 13,000 rcf/min for 10 min at 4 °C, then transfer the supernatant to a nuclease-free 2 ml microfuge tube.
- 7 Add 700 µl chloroform-isoamylalcohol (24:1) and vortex to mix well. Centrifuge at 13,000 rcf/min for 10 min at 4 °C, then transfer the supernatant to a nuclease-free 1.5 ml microfuge tube.
- 8 Add 1/10 volume of 3M sodium acetate, 2 µl of 5 mg/ml glycogen and 800 µl cold isopropanol. Incubate at -20 °C for the night.
- 9 Centrifuge at 13,000 rcf/min for 10 min at 4 °C, then discard the supernatant.
- 10 Washed with 800 µl of 70% ethanol. Centrifuge at 13,000 rcf/min for 10 min at 4 °C, then discard the supernatant.
- 11 Wash the precipitate with a second 800 µl of 70% ethanol and centrifuge as above. discard the supernatant.
- 12 Let sit for 15 min to dry the precipitate. Add 1×TE buffer to resuspend the precipitate.
- 13 Add 1µl RNase A to samples and invert tube several times to mix. Incubate at 37 °C for 30 minutes.