



JUL 13, 2023

OPEN  ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.8epv59e16g1b/v1

Protocol Citation: Mark Ellisman, Mason Mackey, Eric Bushong, Marianne Martone, Thomas Deerinck 2023. Ringer's Solution.

protocols.io
<https://dx.doi.org/10.17504/protocols.io.8epv59e16g1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Apr 20, 2022

Last Modified: Jul 13, 2023

PROTOCOL integer ID:
61105

Keywords: ringers, perfusion, mice, Lidocaine, Heparin, salts, solution

Ringer's Solution

Mason

Mark Ellisman¹, Mackey¹, Eric Bushong¹,
Marianne Martone¹, Thomas Deerinck¹

¹National Center for Microscopy and Imaging Research, University of California, San Diego, La Jolla CA

NCMIR@UCSD



NCMIR

University of California, San Diego

ABSTRACT

Standard protocol for preparing 100mL of Ringer's Solution for mice and rat perfusions.

GUIDELINES

Make sure you bubble for at least 5 minutes before adding calcium chloride to Ringer's solution to prevent precipitation. If precipitation occurs, discard and remake entire solution.

MATERIALS

- 100mL Erlenmeyer flask
- 1000 μ L Pipette
- Serological Pipettes

SAFETY WARNINGS



Always wear appropriate PPE.

- 1 1.0 mL (Na_2HPO_4 , 18 g/L, at room temperature)
- 2 9.9 mL (NaCl, 79.8 g/L, 4 °C fridge)
- 3 1.0 mL (KCl, 37.5 g/L, 4 °C fridge)
- 4 1.0 mL ($\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$, 20 g/L, 4 °C fridge)
- 5 2.5 mL (NaHCO_3 , 50 g/L, 4 °C fridge)
- 6 0.200g of Dextrose
- 7 Bring solution up to 97ml with double distilled water (ddH_2O).
- 8 1.0 ml Lidocaine hydrochloride
- 9 1.0 ml heparin (prevents blood clotting)
–bubble 95% air: 5% CO_2 at 35 °C

10 1.0 mL ($\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$, 30 g/L, 4 °C fridge)