



Jan 29, 2022

Library construction for human placenta bulk RNAseq

Scott Lindsay-Hewett¹¹University of California, San Diego

1

dx.doi.org/10.17504/protocols.io.b4fyqtpw

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu

Scott Lindsay-Hewett

This protocol describes the generation of stranded RNA-seq libraries from placenta total RNA. Since RNA quality can be an issue with total RNA isolated from placental tissue, it is advisable to enrich mRNA using a ribodepletion method rather than polyA selection.

Prior to ribodepletion, it is especially important to ensure that the input RNA is free of contaminating DNA. This protocol, therefore, begins with DNase treatment of samples, using Ambion's DNA-free DNase Treatment and Removal Reagents. Next, ribodepletion and library construction are performed using the KAPA RNA HyperPrep Kit with RiboErase (HMR). Libraries are indexed with KAPA Unique Dual-Indexed Adapters to enable multiplexed sequencing on an Illumina instrument.

DOI

dx.doi.org/10.17504/protocols.io.b4fyqtpwScott Lindsay-Hewett 2022. Library construction for human placenta bulk RNAseq. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.b4fyqtpw>

protocol ,

Jan 28, 2022

Jan 29, 2022

57560

Veriti 96-Well Thermal Cycler

Applied Biosystems 4375786

Qubit
Fluorometer

Invitrogen Q33228 [↗](#)



Centrifuge
Benchtop Centrifuge

Eppendorf 5405000441 [↗](#)
Any benchtop centrifuge will suffice



Set of micropipettes with rack: 100-1000
μl, 20-200 μl, 2-20 μl, and 0.5-10 μl
Pipettor set

Pipetman QP-1001-07 [↗](#)
Can use equivalent Pipettors



Vortex mixer

Any xx

Magnetic Stand
Magnetic Stand

Thermo Scientific MR02 [↗](#)
Any magnetic rack that fits your tubes will suffice.



Bioanalyzer
Bioanalyzer

Agilent G2991AA [↗](#)
Any bioanalyzer will suffice.



[↗](#) Nuclease-free Water Contributed by users

[↗](#) 1M Tris pH 8.0 Sigma

Aldrich Catalog #T2694

[↗](#) 200 Proof Ethanol pure Sigma

Aldrich Catalog #E7023

[↗](#) Qubit RNA BR Assay Kit Thermo Fisher

Scientific Catalog #Q10211

[↗](#) Qubit dsDNA HS Assay kit Thermo Fisher

Scientific Catalog #Q32854

[↗](#) Agilent High Sensitivity DNA Kit Agilent

Technologies Catalog #5067-4626

[↗](#) DNA-free DNA Removal

kit Invitrogen Catalog #AM1906

[↗](#) KAPA RNA HyperPrep Kit with RiboErase (HMR) Kapa

Biosystems Catalog #KK8560 (24 libraries) or KK8561

[↗](#) KAPA Unique Dual-Indexed Adapter Kit Kapa

Biosystems Catalog #KK8727

DNA removal

- 1 To remove any possible contaminating DNA prior to ribodepletion, follow Ambion's protocol (Publication # 1906M, Revision E) for their DNA-*free* DNase Treatment and Removal Reagents.

Follow instructions for routine DNase treatment.

Tip: when transferring supernatant during the final step, be sure not to draw up any of the DNase Inactivation Reagent. Take off a volume that you are comfortable with. There is practically zero loss of RNA using this method, only what you leave behind. If you need to maximize your RNA yield you can always increase the starting volume by adding nuclease-free water, but be sure that the RNA will still be at least **100 ng/ul** for library construction in the next step.

Tip: The DNase Inactivation Reagent settles very quickly, so it is important to follow the instruction to resuspend it immediately before use, and several times during incubation with the reaction mix.

Quantitate DNase-treated total RNA using Qubit RNA Broad Range Assay.

Ribodepletion and library construction 10m 30s

- 2 Follow KAPA's protocol (KR1351 - v4.21) for their KAPA RNA HyperPrep Kit with RiboErase (HMR). In our hands, this protocol yields great results with partially degraded total RNA samples from placenta, some with RIN scores as low as 2. ^{10m 30s}

 **KAPA RNA HyperPrep Kit with RiboErase (HMR) KR1351 - v4.21.pdf**

For HuBMAP samples, the following parameters were used:

Input total RNA: **1000 ng**

Fragmentation: tailored to RIN score (RIN score 6.0 and above, **00:04:00** **85 °C** ; RIN score 5.0-5.9, **00:03:30** **85 °C** ; RIN score 4.9 and below, **00:03:00** **85 °C**)

Library amplification: 8 cycles

For indexing, and to enable efficient multiplexed sequencing, use the KAPA Unique Dual-Indexed Adapter Kit , and dilute adapters to **7 micromolar (μM)** . Using UDIs reduces "index hopping" which can be an issue when using combinatorial dual indexes, especially on patterned flow cells like the NovaSeq.


 **KAPA Unique Dual-Indexed Adapter Kit KR1736 - v3.20.pdf**

Be extremely careful when pipetting out of the adapter plate - you do not want to cross-contaminate wells. Make sure to centrifuge the adapter plate before carefully peeling off the foil cover. Promptly replace with a fresh foil cover when done. Use provided KAPA Adapter Dilution Buffer to dilute adapters from **15 micromolar (μM)** to **7 micromolar (μM)** .

If using a thermocycler to incubate the adapter ligation reaction at **20 °C** , do not pre-heat the lid as this may cause the tube/plate temperature to be too warm, allowing DNA ends to breathe and reducing adapter ligation efficiency.

Quality control


- 3 Quantitate libraries using the Qubit DNA High Sensitivity Assay, and check library distribution by running the DNA High Sensitivity Assay on an Agilent Bioanalyzer. A typical trace is shown below.



- 4 If multiplexing samples, first perform a balancing run to ensure equal representation of all samples in the pool. For the balancing run, prepare an "equal volume" pool by combining **2 μ L** each library together. Run the pool on a MiSeq instrument using a MiSeq Reagent Kit Nano. Based on the proportion of reads assigned to each index during the Nano run, prepare a balanced pool that will yield an equal read depth for all samples in the pool.
- 5 Prior to submitting for sequencing, quantitate the pool using the Qubit High Sensitivity DNA assay. Determine the average fragment size for each library from the Bioanalyzer traces.

Tip: set a region from 200bp to 1000bp in the region table tab (fragments over ~1kb don't cluster efficiently on the flow cell) and the Agilent software will calculate the average fragment size.

Determine the average fragment size in the balanced pool, and use the following formula to determine the nM concentration:


$$\frac{(\text{concentration in ng/}\mu\text{l})}{(660 \text{ g/mol} \times \text{average library size in bp})} \times 10^6 = \text{concentration in nM}$$

Submit the pool to your sequencing facility, noting the nM concentration.

For HuBMAP bulk RNA-seq samples, the multiplexed pool was sequenced on a NovaSeq 6000 S4 lane using a 100bp paired-end run configuration.