



Mikro SARS-CoV-2 Multiplex Protocol V.1

Toby Overmaat¹

¹Mikro

Version 1

Sep 10, 2020

In Development

dx.doi.org/10.17504/protocols.io.bk8dkzs6

Toby Overmaat

DOI

dx.doi.org/10.17504/protocols.io.bk8dkzs6

PROTOCOL CITATION

Toby Overmaat 2020. Mikro SARS-CoV-2 Multiplex Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk8dkzs6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 10, 2020

LAST MODIFIED

Sep 10, 2020

PROTOCOL INTEGER ID

41957

MATERIALS

NAME	CATALOG #	VENDOR
Ethanol, Absolute, Molecular Biology Grade	BP2818500	Thermo Fisher Scientific
Buffer AVL	19073	Qiagen
Buffer AVE	1020953	Qiagen
UltraPlex 1-Step ToughMix (4X)	95166-01K	Quantabio
Nuclease-Free Water	129117	Qiagen

STEPS MATERIALS

NAME	CATALOG #	VENDOR
Ethanol, Absolute, Molecular Biology Grade	BP2818500	Thermo Fisher Scientific
Buffer AVE	1020953	Qiagen
Buffer AVL	19073	Qiagen
UltraPlex 1-Step ToughMix (4X)	95166-01K	Quantabio

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors,

BEFORE STARTING

The protocol is designed for a single sample, utilizing a single well.

It is designed for a custom RT-qPCR device, but may be adapted for other equipment.

Sample Collection

- 1 Remove the lid of sample collector. Direct the sample provider to expel saliva into the collector until 0.5 mL has been collected. Replace lid.



Teeth brushing, mouthwash, drinks, food, and nasal sprays should be avoided for half an hour before sample collection.

Heat Pre-Treatment

- 2 Load sample for RNA extraction.

2.1 Apply heat-treatment to sample

Step	Action	Specification	Time
1	Heat	95°C	2 min
2	Cool	Room temperature	2 min
3	Centrifuge	500 x g	2 min

2.2 Collect cellular debris.

Add  100 µl of



Ethanol, Absolute, Molecular Biology Grade

by Thermo Fisher Scientific

Catalog #: BP2818500

Mix for  00:00:30

RNA Extraction

- 3 Immediately proceed with RNA extraction. Do not store sample.

3.1 Prepare 500 µl of



Buffer AVL

by Qiagen

Catalog #: 19073

3.2 Add 100 µl of heat-treated sample to buffer.

Mix for  00:00:15

3.3 Incubate at room temperature for 00:07:00

3.4 Centrifuge at 500 x g for 00:02:00

Then remove buffer.

4 Proceed with RNA elution.

4.1 Add 60 µl of



Buffer AVE

by Qiagen

Catalog #: 1020953


4.2 Centrifuge at 1000 x g for 00:02:00


RT-qPCR

5 Prepare 8 µl multiplex primer-probe-water mix in advance, containing the following:

Component	Spec	Concentration	Volume (1 reaction)
2019-nCoV_N1-F	US CDC	200nM	0.6 µL
2019-nCoV_N1-R	US CDC	400nM	1.2 µL
2019-nCoV_N1-P	US CDC	200nM	0.6 µL
2019-nCoV_N2-F	US CDC	200nM	0.6 µL
2019-nCoV_N2-R	US CDC	400nM	1.2 µL
2019-nCoV_N2-P*	Custom probe (JUN-QSY)	200nM	0.6 µL

RP-F	US CDC	200nM	0.6 µL
RP-R	US CDC	400nM	1.2 µL
RP-P*	Custom probe (VIC-MGB)	200nM	0.6 µL
Nuclease-free water			0.8 µL



- 6 Prepare  **140 µl** of master mix in advance, using



UltraPlex 1-Step ToughMix (4X)
by Quantabio
Catalog #: 95166-01K

with the following quantities:

Component	Volume (1 reaction)
UltraPlex 1-Step ToughMix	15 µL
Primer-probe-water mix	8 µL
Nuclease-free water	117 µL

- 7 Add  **60 µl** of previously purified RNA sample to  **140 µl** of master mix.
Load onto qPCR device.

- 8 Proceed with thermocycling, under the following conditions:

Step	Temperature	Time
1	50°C	10 min
2	95°C	1 min
3	95°C	3 sec
4	60°C	30 sec
Repeat 40 Cycles for Step 3-4		