



Version 1 ▼

Aug 04, 2022

SNP Genotyping and ApoE Genotyping V.1

Huw Morris^{1,2}, Manuela MX Tan^{1,3}, Donald G Grosset⁴, Nigel M Williams⁵

¹Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815.;

³Department of Neurology, Oslo University Hospital, Oslo, Norway;

⁴Institute of Neurological Sciences, Queen Elizabeth University Hospital, Glasgow, UK;

⁵Institute of Psychological Medicine and Clinical Neurosciences, MRC Centre for Neuropsychiatric Genetics and Genomics, Cardiff University, Cardiff, UK

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.by9ypz7w

h.morris

ABSTRACT

This protocol details the steps for DNA extraction from a human blood sample, quality control, and SNP and APOE genotyping.

The protocol has been adapted from the **PRoBaND SNP Genotyping and ApoE Genotyping Protocol**.

The overall protocol for PRoBaND /Tracking Parkinson's is published:

Malek, N., Swallow, D. M. A., Grosset, K. A., Lawton, M. A., Marrinan, S. L., Lehn, A. C., Bresner, C., Bajaj, N., Barker, R. A., Ben-Shlomo, Y., Burn, D. J., Foltynie, T., Hardy, J., Morris, H. R., Williams, N. M., Wood, N., & Grosset, D. G. (2015). Tracking Parkinson's: Study Design and Baseline Patient Data. *Journal of Parkinson's Disease*, 5(4), 947–959. <https://doi.org/10.3233/JPD-150662>

DOI

dx.doi.org/10.17504/protocols.io.by9ypz7w

EXTERNAL LINK

<https://doi.org/10.1136/jnnp-2021-328365>

PROTOCOL CITATION

Huw Morris, Manuela MX Tan, Donald G Grosset, Nigel M Williams 2022. SNP Genotyping and ApoE Genotyping . **protocols.io**
<https://protocols.io/view/snp-genotyping-and-apoe-genotyping-by9ypz7w>



FUNDERS ACKNOWLEDGEMENT

The Michael J. Fox Foundation for Parkinson's Research (MJFF) and the Aligning Science Across Parkinson's (ASAP) initiative.

Grant ID: ASAP-000478

Parkinson's UK

Grant ID: H-1703

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

<https://content.iospress.com/articles/journal-of-parkinsons-disease/jpd150662>

KEYWORDS

Parkinson's disease, Blood, DNA, Genotyping, Apolipoprotein E (APOE), Single Nucleotide Polymorphisms (SNPs), ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 20, 2021

LAST MODIFIED

Aug 04, 2022

OWNERSHIP HISTORY

Oct 20, 2021  Oke Avwenagha

Oct 26, 2021  h.morris

PROTOCOL INTEGER ID

54296

GUIDELINES

The 1000 genomes reference panel:

1000 Genomes Project Consortium, Auton, A., Brooks, L. D., Durbin, R. M., Garrison, E. P., Kang, H. M., Korbel, J. O., Marchini, J. L., McCarthy, S., McVean, G. A., & Abecasis, G. R. (2015). A global reference for human genetic variation. *Nature*, 526(7571), 68–74. <https://www.nature.com/articles/nature15393>

MATERIALS TEXT

BD Vacutainer EDTA Blood Collection Tubes:

<https://www.camlab.co.uk/bd-vacutainer-edta-blood-collection-tubes>

Illumina HumanCore Exome Array

<https://www.illumina.com/products/by-type/microarray-kits/infinium-exome.html>

SAFETY WARNINGS


The process of venepuncture has the potential to expose research and clinical personnel and other health care professionals and participants to blood from other people, thereby putting them at risk from blood borne pathogens.

In order to minimise risk to the participants, research and clinical personnel and other health care professionals it is imperative to ensure that personnel performing a venepuncture have received appropriate training and education to ensure they are competent. The procedure, the use of PPE and disposal of sharps etc., must be in accordance institution policies, approved risk assessments and procedures.

BEFORE STARTING

Check the planned procedure is covered by existing Research Ethics Committee.
Check that you have the correct PPE, supplies and consumables.

Sample Collection and Genotyping

- 1 At study entry (baseline), collect a  **4 mL** blood sample in an ethylenediaminetetraacetic acid (EDTA) vacutainer for DNA extraction.



The process of venepuncture has the potential to expose research and clinical personnel and other health care professionals and participants to blood from other people, thereby putting them at risk from blood borne pathogens.

In order to minimise risk personnel performing a venepuncture must have received appropriate training and deemed competent. This procedure, the use of PPE and the

use and disposal of sharps etc., must be in accordance institution policies, approved risk assessments and procedures.

2 Perform DNA extraction.

DNA extraction was carried out at [K biosciences](#), UK

3 Genotype DNA samples using the Illumina HumanCore Exome array with custom content following manufacturers instructions : <https://www.illumina.com/products/by-type/microarray-kits/infinium-exome.html>

This covered approximately 250,000 common variants, 250,000 rare variants, and over 27,000 custom variants that have been implicated in neurological and psychiatric disorders.

Purcell, S., Neale, B., Todd-Brown, K., Thomas, L., Ferreira, M. A. R., Bender, D., Maller, J., Sklar, P., de Bakker, P. I. W., Daly, M. J., & Sham, P. C. (2007). PLINK: a tool set for whole-genome association and population-based linkage analyses. *American Journal of Human Genetics*, 81(3), 559–575. <https://doi.org/10.1086/519795>

Genotyping was performed Cardiff University: <https://www.cardiff.ac.uk/mrc-centre-neuropsychiatric-genetics-genomics>

Quality Control and Principle Component Analysis

4 Perform standard quality control procedures in PLINK v1.9: <https://www.cog-genomics.org/plink/>

5 Remove individuals:

- 5.1 With low overall genotyping rates ($< 98\%$), related individuals (Identity-By-Descent PIHAT > 0.1) and heterozygosity outliers (> 2 SDs away from the mean)
- 5.2 Whose clinically reported biological sex did not match the genetically determined sex.
- 6 Conduct principle component analysis on a linkage disequilibrium (LD) pruned set of variants (removing SNPs with an $r^2 > 0.05$ in a 50kb sliding window shifting 5 SNPs at a time) after merging with European samples from the HapMap reference panel.
- 7 Remove:
 - 7.1 Individuals who are > 6 SDs away from the mean of any of the first 10 principal components.
 - 7.2 Variants if they had a low genotyping rate ($< 99\%$), Hardy-Weinberg Equilibrium p -value $< 1 \times 10^{-5}$, and minor allele frequency $< 1\%$.

1,000 Genomes Project and APOE Genotyping

- 8 Following quality control, input genotypes separately to the 1,000 Genomes Project reference panel (phase 3 release 5) using the Michigan Imputation Server (<https://imputationserver.sph.umich.edu>).
- 9 Retain only variants with imputation quality > 0.8 , to keep only high quality calls to merge across the cohorts.
- 10 Tracking Parkinson's data was lifted over to genome build hg38 using liftOver (<https://genome.ucsc.edu/cgi-bin/hgLiftOver>).

Twenty genetic principal components are generated from a linkage-pruned SNP set (removing SNPs with an $r^2 > 0.02$ in a 1000kb sliding window shifting 10 SNPs at a time).

Extreme outliers are removed according to the first 5 principal components (> 6 SDs away from the mean).

- 11 Recalculate the genetic principal components after removing outliers, as extreme outliers can substantially affect the calculation of genetic principal components.
 - 11.1 These first 5 new principal components are included as covariates in future GWAS to adjust for population substructure
 - 11.2 Additional outliers who are > 6 SDs away from the mean of any of the first 5 principal components are excluded.
- 12 Restrict genotyping of APOE to samples and SNPs that survive these QC procedures based on two common SNPs of the APOE gene: 388 T > C (rs429358) and 526C > T (rs7412).
- 13 The three haplotypes ($\epsilon_2(388\text{ T}-526\text{ T})$, $\epsilon_3(388\text{ T}-526\text{ C})$, $\epsilon_4(388\text{ C}-526\text{ C})$) and six genotypes (ϵ_2/ϵ_2 , ϵ_2/ϵ_3 , ϵ_2/ϵ_4 , ϵ_3/ϵ_3 , ϵ_3/ϵ_4 , ϵ_4/ϵ_4) formed by these SNPs are then used to determine ϵ_4 status in each sample.