



Aug 04, 2021

Integration of a cargo brick

In 1 collection

Carolyn N Bayer¹, Maja Rennig¹, Anja Ehrmann¹, Morten Norholm¹¹Technical University of Denmark

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bvk9n4z6

SEGA

Carolyn Bayer
Technical University of Denmark

ABSTRACT

SEGA, the Standardized Genome Engineering Architecture, is a comprehensive strain collection that enables genome engineering by combining only two reagents: a DNA fragment that can be ordered from a commercial vendor and a stock solution of bacterial cells followed by incubation on agar plates. Recombinant genomes are identified by visual inspection using green-white colony screening akin to classical blue-white screening for recombinant plasmids. The modular nature of SEGA allows precise multi-level control of transcriptional, translational, and post-translational regulation. The SEGA architecture simultaneously supports increased standardization of genetic designs and a broad application range by utilizing well-characterized parts optimized for robust performance in the context of the bacterial genome

This protocol describes the process of integrating a SEGA cargo brick. A cargo brick is integrated using *tetA* counterselection. This protocol also applies to integration of other bricks using *tetA* counterselection, e.g. splitting *tetA*.

DOI

dx.doi.org/10.17504/protocols.io.bvk9n4z6

PROTOCOL CITATION

Carolyn N Bayer, Maja Rennig, Anja Ehrmann, Morten Norholm 2021. Integration of a cargo brick.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bvk9n4z6>

COLLECTIONS ⓘ

 **SEGA protocol collection**

KEYWORDS

SEGA, genome engineering, recombineering, E. coli, synthetic biology

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 08, 2021

LAST MODIFIED

Aug 04, 2021

PROTOCOL INTEGER ID

50561

PARENT PROTOCOLS

Part of collection

[SEGA protocol collection](#)

BEFORE STARTING

Transform a SEGA strain with pSIM19 (Spectinomycin resistance). From now on cultures have to be kept at **30 °C** to retain the plasmid (temperature-sensitive *ori*)

preculture and DNA fragment- Day 1

1d

1d

- 1 Prepare a PCR product of the cargo brick and purify it from an agarose gel.
- 2 Setup a preculture of the strain with pSIM19 in LB medium supplemented with Spectinomycin (**0.05 mg/ml**) and incubate at **250 rpm, 30°C** overnight

Recombineering- Day 2

1d

- 3 Prepare:

3.1 Cold sterile water

3.2 Cold Glycerol (**15 % volume**)

3.3 Pre-chilled centrifuge and tabletop centrifuge at **4 °C**

3.4 M9 agar plates supplemented with (**50 Micromolar (µM)**) NiCl₂

Additionally, specific inducers to enable green-white screening can be added

- 4 Inoculate **50 mL** LB-Medium supplemented with Spectinomycin (**0.05 mg/ml**) with **500 µl** of the preculture from step 3



4.1 Incubate at **250 rpm, 30°C** until cultures reached an OD₆₀₀ of 0.5

- 5 Induce expression by transferring the culture to a shaking water bath at 🔄 **150 rpm, 42°C, 00:20:00**
 - 6 Transfer culture to prechilled 📄 **50 mL** falcon tubes and put on ice for ⌚ **00:15:00** 15m
 - 7 Spin the culture down at 🌀 **4000 x g, 4°C, 00:05:00** and discard the supernatant 5m
 - 8 Add 📄 **1 mL** of ice cold water, resuspend and transfer to a 1.5 ml tube
 - 9 Spin at 🌀 **11000 x g, 4°C, 00:00:30** in a tabletop centrifuge 30s
 - 10 Wash pellet twice with 📄 **1 mL** ice cold water ➡ **go to step #9**
 - 11 Resuspend the pellet in 📄 **600 µl** cold glycerol (**[M]15 % volume**)
 - 11.1 Unused cells can be stored at 🧊 **-80 °C**
- This is not possible for *E. coli* Nissle
- 12 Electroporate 📄 **50 µl** of cells with 📄 **200 ng** of purified PCR product from step 2 or 📄 **2 µl** of a **[M]100 Micromolar (µM)** single stranded oligonucleotide
 - 13 Recover cells 🔄 **800 rpm, 30°C, 01:00:00** in a tabletop shaker using SOC medium.
 - 14 Transfer the cells into 📄 **5 mL** LB medium supplemented with Spectinomycin
 - 15 Incubate at 🔄 **250 rpm, 30°C** overnight

Plating- Day 2


1d

16 Wash  **1 mL** of the recovered cells twice with sterile water. Centrifuge at  **11 000 rpm, 20°C, 00:00:30** 30s

17 Make a dilution series and plate  **100 µl** of the 1:10 - 1:1000 dilution on M9 agar supplemented with  **50 Micromolar (µM)** NiCl_2

the appropriate inducer can be added to the agar plates to enable green-white screening

18 incubate the plates at  **30 °C** for  **48:00:00** to  **72:00:00** 5d

Cells can be incubated at  **37 °C** if it is the last step of recombination and pSIM19 is no longer needed in the cells

Screening- Day 4-5

1d

19 Screen for positive colonies by "green-white screening" on a blue-light table and perform colony PCR on the colorless colonies to identify the correct recombinants.

Streaking out of the colonies at least twice is necessary to ensure a clean culture