



OCT 31, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.kqdg3x3m7g25/v1

Protocol Citation: Ester Kalef-Ezra, Ben Harvey, Katherine Roper, Christos Proukakis 2023. Section 1: Enzymatic DNA Fragmentation (Manually). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg3x3m7g25/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Oct 03, 2023

Section 1: Enzymatic DNA Fragmentation (Manually)

In 1 collection

Ester Kalef-Ezra^{1,2}, Ben Harvey³, Katherine Roper³, Christos Proukakis^{1,2}

¹Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

³Diagnostics and Genomics Group, Agilent Technologies LDA UK Ltd

ASAP Collaborative Research Network

University College London



Ester Kalef-Ezra
 University College London

ABSTRACT

This protocol details manual enzymatic DNA Fragmentation prior to Section 2: NGS library preparation for sequencing.

ATTACHMENTS

[861-2221.pdf](#)

Enzymatic DNA Fragmentation (Manually)

5m

- 1 Group the samples scWGA-products of interest for library preparation based on their size length to require the same fragmentation time. The recommended fragmentation times per amplicon are suggested on Table 1.

Table 1. Recommended fragmentation time based on amplicon size.

A	B
scWGA products peak size (bp)	Fragmentation Time (min)
400-450	3
450-500	4
500-550	5
550-650	6
650-750	7
750-900	8
1000-1100	10
>10.000	15

- 2 Calculate the amount of input material needed from each scWGA-product to have input of 10-200 ng DNA in a total volume of 7 µL .

Note



Note: If possible, use the maximum amount of input DNA available within the recommended range.

- 3 Thaw the vial of 5X SureSelect Fragmentation Buffer On ice , vortex, then keep On ice .

4 Thaw a PCR plate cooler @  Room temperature for  00:05:00 . 5m

5 Pre-label PCR tubes/strips to be used for the enzymatic DNA fragmentation reaction.

6 Place the prelabelled PCR tubes/strips/plate in the PCR plate cooler.

7 Dilute  10-200 ng of each scWGA-product with dH₂O up to a final volume of  7 µL in the PCR tubes/strips.

Note

Note: From now on, keep the samples on ice if not stated differently.



8 Pre-program a thermal cycler as on Table 2. If required, use a reaction volume setting of  10 µL and lid temperature @  60-65 °C .



Table 2. Thermal cycler program for Enzymatic Fragmentation (with the heated lid ON).

A	B	C
Step	Temperature (°C)	Time (min)
1	37	The max need based on Table 1
2	65	5
3	4	hold


9 Prepare the appropriate volume of Fragmentation master mix by combining the reagents as on Table 3 in 0.2ml or 1.5 ml tube.



Table 3. Fragmentation master mix.

A	B
Reagent	1 reaction
5X SureSelect Fragmentation Buffer	2 μ l
SureSelect Fragmentation Enzyme	1 μ l
Total	3 μ l

- 10** Mix well by pipetting up and down 20 times or seal the tube and vortex at high speed for 5–10 sec.



- 11** Spin briefly to remove any bubbles and keep  On ice .

- 12** Add  3 μ L of the Fragmentation master mix to each sample well containing  7 μ L of input DNA.



- 13** Mix well by pipetting up and down 20 times or cap the wells and vortex at high speed for 5–10 seconds.





- 14** Spin the samples briefly and place on the PCR plate cooler.

- 15** Place the tubes/strip in the thermal cycler and start the 'Enzymatic Fragmentation' program.

Note

Note: Place the strips in the thermal cycler in an order to perform the recommended enzymatic fragmentation time based on Table 1. As an example, place first a strip with samples that require 10 min fragmentation and 2 min later a strip with samples that require 8 min fragmentation.

16 When the program reaches the  4 °C Hold step, remove the samples from the thermal cycler.

17 Add  40 µL of nuclease-free water to each sample and place the samples  On ice .



Note

Notes.

1. The 50 µl reactions are now ready for NGS sequencing library preparation (Section 2).
2. This is not a stopping point in the workflow, and analysis of the enzymatically fragmented samples is not required prior to library preparation. Moreover, electrophoretic analysis of the fragmented samples may produce misleading results due to the presence of agents that affect DNA fragment migration.