



Apr 20, 2020

Isolation of SARS-Cov2 RNA from Humans Without High Demand Reagents

In 1 collection

Joseph Patterson¹, Allyson Cole-Strauss¹, John Beck¹, Caryl Sortwell¹, Jack Lipton¹

¹Michigan State University

1

Works for me

dx.doi.org/10.17504/protocols.io.be8ujhww

Joseph Patterson
Michigan State University



ABSTRACT

Viral RNA isolation kits used in PCR based tests for the novel coronavirus (SARS-Cov2) are in short supply. Our group sought to identify a method to isolate viral RNA without the use of a kit or other supplies in high demand by clinical labs. Using a TRIzol based RNA extraction, followed by a glycogen precipitation, we are able to isolate a sufficient quality and quantity of RNA for quantitative PCR, and/or droplet digital PCR. This method has been confirmed to work as well as a viral RNA isolation kit, shown to work with SARS-Cov2 RNA-spiked into human samples, and in detecting SARS-Cov2 in known positive cases of COVID-19.

GUIDELINES

Samples should be processed for RNA extraction (at least up until they can be frozen at -80 °C) within 48 hours of collection.

MATERIALS TEXT

- TRIzol (Ambion 15596026)
- GlycoBlue (Thermo Fisher #AM9515)
- 100% Isopropanol
- 100% Chloroform
- 80% ethanol
- Nuclease Free Water
- DNase/RNase free 1.5 mL microcentrifuge tubes
- Phasemaker tubes (Invitrogen A33248)

SAFETY WARNINGS

Human samples should be handled with care, and sample preparation performed in at least a BSL-2 lab.

- 1 Transfer the swab into a DNase/RNase free 1.5 mL microcentrifuge tube.
•NOTE: The swabs we use can be pulled off their stick for easier processing.
- 2 Add 500 µL of TRIzol to the microcentrifuge tube.
- 3 Close tube and shake for 2 min. by hand.
- 4 Remove the swab and briefly centrifuge.
•NOTE: The swab should be removed from the TRIzol, and squeezed against the side of the tube to retain as much sample as possible.
- 4.1 Samples can be frozen at -80 °C at this point if needed, then thawed for use when ready.

- 5 Centrifuge Phasemaker tubes (Invitrogen A33248) for 30 s at 16,000xg
- 5.1 If Phasemaker tubes are not available, substitute the steps below for steps 6-11:
 - 1) Add 100 μ L of 100% chloroform.
 - 2) Vortex for 30 s.
 - 3) Incubate for 3 min. at room temperature.
 - 4) Centrifuge at 12,000 x g for 10 min. at 4 °C.
 - 5) Transfer supernatant to a new DNase/RNase free 1.5 mL microcentrifuge tube.
 - 6) Continue with protocol at step 7, adding 250 μ L of 100% isopropanol to each tube.
- 6 Transfer the entire sample to a Phasemaker tube and incubate at RT for 5 min
- 7 Add 100 μ L of 100% chloroform to the phasemaker tube.
- 8 Shake the tube for 15 seconds by hand (DO NOT VORTEX).
- 9 Incubate for 10 min. at room temperature.
- 10 Centrifuge for 5 min. at 16,000xg at 4°C.
- 11 Transfer the aqueous phase (clear) to a new DNase/RNase free 1.5 mL microcentrifuge tube, and dispose of the Phasemaker tube.
- 12 Add 250 μ L of 100% isopropanol
- 13 Add 2 μ L of GlycoBlue coprecipitant
- 14 Vortex for 30 s.
- 15 Centrifuge at 12,000 x g for 20 min. at 4 °C.
- 16 Discard the supernatant, keep the pellet.
- 17 Add 80% ethanol to wash the pellet.
- 18 Immediately remove all of the ethanol.
•NOTE: Once most of the ethanol is removed, the tubes can be briefly centrifuged and a pipet used to remove as much of the remaining ethanol as possible.
- 19 Allow pellet to dry for 3 min. at room temperature.
- 20 Resuspend pellet in nuclease free water.
- 20.1 For both droplet digital PCR and qPCR protocols, 25 μ L nuclease free water is an appropriate starting volume, but can be adjusted if needed.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited