

Aug 14, 2024

proteomics (mitochondria)

DOI

dx.doi.org/10.17504/protocols.io.36wgqny15gk5/v1

Livia Hecke Morais¹, Baiyi Quan¹, Tsui-fen Chou¹

¹California Institute of Technology



Livia Hecke Morais

California Institute of Technology

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.36wgqny15gk5/v1**

Protocol Citation: Livia Hecke Morais, Baiyi Quan, Tsui-fen Chou 2024. proteomics (mitochondria). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqny15gk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 14, 2024

Last Modified: August 14, 2024

Protocol Integer ID: 105242

Keywords: ASAPCRN

Funders Acknowledgement:
ASAP

Abstract

Here we describe proteomics experiment with isolated mitochondrial extracts at the Proteome exploration laboratory (PEL) at Caltech by Baiyi Quan, Jeff Jones and Tsui-Fen Chou in collaboratin with Livia Hecke Morais and Sarkis Mazmanian.



Sample Preparation

1

1. Proteins extracted from mitochondria are reduced using a final concentration of 5 mM of TCEP (tris(2-carboxyethyl)phosphine) at room temperature for 10 min.
2. The proteins are further alkylated using a final concentration of 20 mM of CAA (chloroacetamide) at room temperature for 15 min.
3. A 25 ul aliquot of the lysate is acidified using 2.5 ul of 12% phosphoric acid. Using a pH test paper to make sure solution pH< 2.
4. Combine 165 ul of the S-trap buffer (90% MeOH with 100 mM triethylammonium bicarbonate) with the acidified lysate.
5. Transfer the combined colloid into the S-trap micro (Protifi, NY).
6. Centrifuge the S-trap micro at 4000 g for 1 min at room temperature.
7. Using 150 ul of the S-trap buffer to wash the S-trap and centrifuge at 4000 g for 1 min. Repeat the step two more times.
8. Add 20 ul of 100 mM triethylammonium bicarbonate containing 20 ug of Trypsin into S-trap micro. Allow the incubation stay overnight (14-16 hr) at 37 C.
9. After overnight incubation, directly adding 40 ul of the 50 mM triethylammonium bicarbonate and centrifuge at 4000 g for 1 min.
10. Add 40 ul of 2% formic acid in water to the S-trap, and centrifuge at 4000 g for 1 min.
11. Add 40 ul of 50% acetonitrile in water to the S-trap, and centrifuge at 4000 g for 1 min.
12. Speedvac the sample.
13. The sample is reconstituted in 20 ul of 2% acetonitrile and 0.2% formic acid in water and 500 ug of the peptide is used for LC-MS/MS analysis.

LC-MS/MS

- 2 1. An aliquot of sample containing 500 ug of the peptides is subjected to the LC-MS/MS analysis. The sample is separated on an Aurora UHPLC Column (25 cm × 75 µm, 1.6 µm C18, AUR2-25075C18A, Ion Opticks) using an Easy-nLC 1200 liquid chromatography system. The gradient settings follows **Table 1**.

Time	Duration	Flow (nl/min)	%B
0:00	0:00	350	3
1:00	1:00	350	3
73:00	72:00	350	19
101:00	28:00	350	29

121:00	20:00	350	41
124:00	3	350	95
131:00	7	350	98

Table 1. LC gradient for the sample

Mobile Phase A: 0.2% formic acid, 2% acetonitrile, and 97.8% water.

Mobile Phase B: 0.2% formic acid, 80% acetonitrile, and 19.8% water.

1. The sample is analyzed on a Thermo Q-Exactive HF mass spectrometer using a data-dependent acquisition method. Detailed parameters of the scans are listed in **Table 2**.

Global settings	
Ion source type	NSI
Spray voltage	2000 V
Ion transfer tube temperature	300 C
Polarity	Positive
MS1 scan settings	
Resolution	60000
AGC target	3e6
Maximum IT	15 ms
Scan range	375-1500 m/z
MS2 scan settings	
Resolution	30000
AGC target	1e5
Maximum IT	45 ms

Loop count	12
Isolation window	1.2 m/z
NCE	28
Spectrum data type	Centroid
Fixed first mass	100 m/z

Table 2. MS settings for the run.

Proteomic Data analysis

3. 1. The raw data generated by mass spectrometer is analyzed using Proteome Discoverer 2.5. The data is searched using the mouse proteome achieved from UniprotKB on 10/26/2020 (swissprot + trembl). The processing and consensus workflow are set according to **Figure 1** and **2**. The parameters are listed below in **Table 3**. All the parameters that are not mentioned are left defaulted.




Figure1. Processing workflow of the PD analysis.




Figure
2. Consensus workflow of the PD analysis.

4

SequestHT settings	
Enzyme name	Trypsin (Full)
Max. missed cleavage	2
Min. peptide length	6
Max. peptide length	144

Precursor mass tolerance	10 ppm
Fragment mass tolerance	0.02 Da
Max. equal modification	3
Dynamic modification	Oxidation/ +15.995 Da (M)
Dynamic modification (protein terminus)	Acetyl/ + 42.011 Da (N-Terminal)
Dynamic modification (protein terminus)	Met-loss/ - 131.040 Da (M)
Dynamic modification (protein terminus)	Met-loss+Acetyl/ - 89.030 Da (M)
Static modification	Carbamidomethyl/ + 57.021 Da (C)
Percolator	
Target/Decoy selection	Concatenated
Validation based on	q-Value
Target FDR (Strict)	0.01
Target FDR (Relaxed)	0.05
Feature Mapper settings	
Maximum RT shift (min)	2
Mass tolerance	5 ppm

Table 3. Parameters for PD searching

1. The protein list exported from PD results are further analyzed using R scripts based on the TidyProteomics package (jeffsocal.github.io/tidyproteomics/articles/overview.html).