



Nov 20, 2020

Primers for bark beetles and associated fungi

Andrew J. Johnson¹, [Demian F Gomez](#)¹, Jiri Hulcr¹¹University of Florida1 [Works for me](#) dx.doi.org/10.17504/protocols.io.bnvcme2w

Bark Beetle Mycobiome Research Coordination Network

ABSTRACT

This protocol describes the primers used for bark beetles and associated fungi.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Seltmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. *Symbiosis* 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

DOI

dx.doi.org/10.17504/protocols.io.bnvcme2w

DOCUMENT CITATION

Andrew J. Johnson, Demian F Gomez, Jiri Hulcr 2020. Primers for bark beetles and associated fungi.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bnvcme2w>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 23, 2020

LAST MODIFIED

Nov 20, 2020

DOCUMENT INTEGER ID

43652

ABSTRACT

This protocol describes the primers used for bark beetles and associated fungi.

This protocol is part of the Bark Beetle Mycobiome (BBM) Research Coordination Network. For more information on the BBM international network: Hulcr J, Barnes I, De Beer ZW, Duong TA, Gazis R, Johnson AJ, Jusino MA, Kasson MT, Li Y, Lynch S, Mayers C, Musvuugwa T, Roets F, Seltmann KC, Six D, Vanderpool D, & Villari C. 2020. Bark beetle mycobiome: collaboratively defined research priorities on a widespread insect-fungus symbiosis. *Symbiosis* 81: 101–113 <https://doi.org/10.1007/s13199-020-00686-9>.

Diluting primers

Reconstituting lyophilized primers to create 100X stock:

Resuspend to 100 uM/L concentration: for 28.5 nM of lyophilized primer, add 285 ul water (10 times more, in uL).

Dilute to make working solution :

To 10uL of the working stock above, add 90 uL of water, to make 10 uM/L (or 10 pM/uL) solution.
Working concentration: Use 1uL of the working solution above in 25uL reaction (10/25=0.4pM per rxn).

Primers we use

GENE	PRIMER CODE	SEQ	REFERENCE
16S	1492R	GGTTACCTTGTTACGACTT	
16S	27F	AGAGTTTGATCCTGGCTCAG	
16S bact rev	338R	TGCTGCCTCCCGTAGGAGT	Amann et al. 1995. Microbiol. Rev. 59, 143–169.
16S bact+archaea for	515F	GTGCCAGCMGCCGCGGTAA	Bergmann et al 2011 (Fierer lab)
16S bact+archaea rev	806R	GGACTACVSGGGTATCTAAT	Bergmann et al 2011 (Fierer lab)
28S	A4053 (alt_name: A1)	TCKGKTTCAAGACGGGGT C	Whiting et al. 1997
28S	A4221 (alt_name: A160)	CGCCTCTTCTCGCAATGAGA	Sequeira et al, 2000
28S	A4285 (alt_name: A247)	CCTGACTTCGTCCTGACCAG GC	Sequeira et al, 2000
28S	A4394 (alt_name: A335)	TCGGARGGAACCAGCTACTA	Whiting et al. 1997
28S	D2F1	ACTGTTGGCGACGATGTTCT	Jordal et al. (2008)
28S	D3R2	TCTTCGCCCCTATACCC	Jordal et al. (2008)
28S	f1094	GGATCGTCGCAAGACGGAC AGAAG	Sequeira et al, 2000
28S	f116 (alt_name: e)	CTGGTTGATCCTGCCACGT	Hamby & Zimmer 1988
28S	f1403	CGGAAGGATTGACAGATTG AGAG	Sequeira et al, 2000
28S	f420	GGCGACGCATCTTCAAATG TCTG	Sequeira et al, 2000
28S	r1094 (alt_name: P)	CGCTTTCGTAAACGGTT	Hamby & Zimmer 1988
28S	r1138	CGCCTTCGAACCTCTAAC	Sequeira et al, 2000
28S	r1626	GGCATCACAGACCTGTTATT GCTCAATCTC	Sequeira et al, 2000
28S	r1856 (alt_name: Q)	CAGCGAGGATGGCTAACTTA	Hamby & Zimmer 1988
28S	r803	CCACCGGCAGGACGTCTC	Sequeira et al, 2000
28S	S3690F	GAGAGTTMAASAGTACGTG AAAC	Dowton & Austin 1998
28S	S4072 (alt_name: S1)	GACCCGTCTTGAAMCAGG A	Whiting et al. 1997
Ap1098rev2		ATATTRTTKGGCARYTGDCC KCCC	[note: TmD 61.1 °C] [primer_reference:Danforth et al. 2006]
Ap787for2		TGCTTYGARCCDAGYCTHGA TTAYTG	[note: TmD 60.0 °C] [primer_reference:Danforth et al. 2006]
Ap835rev1		GCATHACYTCHCCACRCTY TTC	[note: TmD 59.8 °C] [primer_reference:Danforth et al. 2006]
ApCADfor1		GGWTATCCCGTDATGGCBM GWGC	[note: TmD 62.1 °C] [primer_reference:Danforth et al. 2006]

ApCADfor2mod		GATGGGAYCTNRGNAARTT YC	[note: TmD 53.1 °C] [primer_reference:Danforth et al. 2006]
ApCADfor3		CTCHGTKGARTTYGATTGGT GYGC	[note: TmD 60.9 °C] [primer_reference:Danforth et al. 2006]
ApCADfor4		TGGAARGARGTBGARTACGA RGTGGTYCG	[note: TmD 63.1 °C] [primer_reference:Danforth et al. 2006]
ApCADfor4		TGGAARGARGTBGARTACGA RGTGGYCG	This is not good, do not use.
ApCADrev1mod		GCCATYRCYTCBCCYACRCT YTTCAT	[note: TmD 62.2 °C] [primer_reference:Danforth et al. 2006]
ApCADrev1mod		GCCATYRCTBCCTACRCTY TTCAT	This is not good, do not use.
ApCADrev4a		GGCCAYTGNGCNGCCACYG TGTCTATYTGYYTNACC	[note: TmD 68.6 °C] [primer_reference:Danforth et al. 2006]
ArgK	ArgK-forB2	GAYTCCGGWATYGGWATCT AYGCTCC	Dole et al. 2010
ArgK	ArgK-revB1	TCNGTRAGRCCCATWCGTC TC	Danforth et al. 2005
ArgK	ArgK-revB2	GTATGYTCMCCRCGRGTAC CACG	Dole et al. 2010
cellobiohydrolase	fungCBH1f		
cellobiohydrolase	fungCBH1r		
CO1	1495b	AACAAATCATAAAGATATTG GRAC	Smith and Cognato, 2014
CO1	A2237R	CCGAATGCTTCTTTTTTACC TCTTTCTTG	
CO1	A2411R	GCTAATCATCTAAAACTTT AATTCCWGTWG	
CO1	HCO2198	TAAACTTCAGGGTGACCAAA AAATCA	Hebert et al. 2003
CO1	LC01490	GGTCAACAAATCATAAAGAT ATTGG	Hebert et al. 2003
CO1	LepF1	ATTCAACCAATCATAAAGAT ATTGG	Hebert et al. 2004
CO1	LepR1	TAAACTTCTGGATGTCCAAA AAATCA	Hebert et al. 2004
CO1	r750	GAAATTATNCCAATTCCTGG	Smith and Cognato, 2014
CO1	S1718F	GGAGGATTTGGAATTGAT TAGTTCC	
EF1a	A1043	GTATATCCATTGGAATTG ACCNGGRTGRTT	Normark et al. 1999
EF1a	S149	ATCGAGAAGTTCGAGAAGG AGGCYCARGAAATGGG	Normark et al. 1999
efTu	52F		Ludwig et al
ITS (end of SSU)	ITS1	TCCGTAGGTGAACCTGCGG	
ITS (in 5.8S, for)	ITS3	GCATCGATGAAGAACGCAGC	
ITS (in 5.8S, rev)	ITS2	GCTGCGTTCTTCATCGATGC	
ITS (rev, start of LSU)	ITS4	TCCTCCGCTTATTGATATGC	

ITS fungal (like ITS1, basidio)	ITS1-F	CTTGGTCATTTAGAGGAAGTAA	
LSU (28S)	JH-LSU-91	AGTAACGGCGAGTGAAG	probably universal for fungi, not for insect, position 91 in <i>S. cerevisiae</i> 28S
LSU (28S)	LR0R	ACCCGCTGAACTTAAGC	identical in insect and fungi
LSU (28S) rev	JH-LSU-377	TTTTCAAAGTGCTTTTC	almost identical to LR2 but with one bp more specific to filamentous ascos, position 377 in <i>S. cerevisiae</i> ; doesn't work very well
LSU (28S) rev	JH-LSU-638	GGTCCGTGTTTCAAGAC	reverse complement to Vilgalis's LR3R (which is not the same as reverse complement of LR3)
LSU (28S) rev	JH-LSU-643	CTCCTTGGTCCGTGTTTC	tested with <i>Raffaelea</i> spp., inconsistent product
LSU (28S) rev	JH-LSU-830	GTTTCCTCTGGCTTCACCC	tested with <i>Raffaelea</i> spp., inconsistent product
LSU (28S) rev	LR2	TTTTCAAAGTTCTTTTC	White et al., from yeast, 1bp difference from filamentous ascos
LSU (28S) rev	LR3	CCGTGTTTCAAGACGGG	identical for beetles and fungi
SSU (18S) general	NS1	GTAGTCATATGCTTGTCTC	eukaryotes?
SSU (18S) general	NS4	CTTCCGTCAATTCTTTAAG	eukaryotes?
SSU fungal	AU2	TTTCGATGGTAGGATAGDGG	Vandenkoornhuyse et al. 2002, Science 295:2051
SSU fungal	AU4	RTCTCACTAAGCCATTC	Vandenkoornhuyse et al. 2002, Science 295:2051
SSU fungal	JH-SSU-304	CCCTATCAACTTTCGATG	upstream of AU2, 2bp difference from Curculionidae, can use with NS2 for 240bp conserved product or with JH-SSU-822 for ~500bp of variable product
SSU fungal (rev)	JH-SSU-1646	CTAAGCCATTCAATCGGTA	6bp upstream of AU4, 2bp difference from Curculionidae
SSU fungal (rev)	JH-SSU-822	CTAGAAACCAACAAAATAG	2bp upstream of SR4, 4bp difference from Curculionidae
SSU rev	NS2	GGCTGCTGGCACCAGACTTGC	550 in <i>S. cerevisiae</i> ; eukaryotes?

R= A G (purine) D= G A T

Buying primers

ITS1 example CType: DNA Oligos – modified & non-modified

Oligo Name: ...

Synthesis Scale: 0.025 µmol

Purification: Desalt

Sequence (5' to 3'): ...

Format: Lab-ready



Citation: Andrew J. Johnson, Demian F Gomez, Jiri Hulcr (11/20/2020). Primers for bark beetles and associated fungi.
<https://dx.doi.org/10.17504/protocols.io.bnvcme2w>

This is an open access protocol distributed under the terms of the **Creative Commons Attribution License** (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Number of Tubes per Oligo: 1
Concentration (μM): NA