



Oct 14, 2020

## Transformation

Hung Liang Pai<sup>1</sup>, Huan Jui Chang<sup>1</sup><sup>1</sup>Chung Shan Medical University

1

Works for me

This protocol is published without a DOI.



Hung Liang Pai

Chung Shan Medical University

### ABSTRACT

This protocol is designed to transform plasmids from distribution kits of iGEM 2017, 2018, 2019 as well as those synthesized by Mission Biotech with yT&A backbone into E.coli DH5- $\alpha$  competent cells. Chloramphenicol is used for selection in transformation when the plasmids are retrieved from iGEM distribution kits, while ampicillins is utilized when the plasmid are derived from Mission Biotech.

### ATTACHMENTS

[Transformation.pdf](#)

### PROTOCOL CITATION

Hung Liang Pai, Huan Jui Chang 2020. Transformation . **protocols.io**  
<https://protocols.io/view/transformation-bg4ajyse>

### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Jun 03, 2020

### LAST MODIFIED

Oct 14, 2020

### PROTOCOL INTEGER ID

37730

### ATTACHMENTS

[Transformation.pdf](#)

### GUIDELINES

The ratio between competent cells and plasmid is 10 : 1. We will apportion 20  $\mu$ L of competent cells for this protocol.

### SAFETY WARNINGS

Make sure that all the steps should be operated in laminar flow (except heat shock).

### ABSTRACT

This protocol is designed to transform plasmids from distribution kits of iGEM 2017, 2018, 2019 as well as those synthesized by Mission Biotech with yT&A backbone into E.coli DH5- $\alpha$  competent cells. Chloramphenicol is used for selection in transformation when the plasmids are retrieved from iGEM distribution kits, while ampicillins is utilized when the plasmid are derived from Mission Biotech.

### BEFORE STARTING

Follow the "preparation" section of protocol.

#### Preparation

Sterilize styrofoam box with alcohol.

- 1
- 2 Turn on the UV light of the laminar flow for 🕒00:10:00 to 🕒00:15:00 in order to sterilize it.
- 3 Add ice into the styrofoam box and sterilize it again.
- 4 Place the styrofoam into the laminar flow.
- 5 Set the dry bath incubator at 🌡️42 °C

#### Protocol

- 6 Apportion 20uL of E-coli DH5α into an eppendorf.



From step 6 to 8, perform them on ice and in the laminar flow.

- 7 Add 2uL of target plasmid into the eppendorf and vortex
- 8 Ice the eppendorf for 3-5 min.
- 9 Heat shock the eppendorf for 1 min at 42°C in a dry bath incubator.
- 10 Ice the eppendorf for 3-5 min.
- 11 Add 100uL of antibiotic-free LB media into the eppendorf.
- 12 Incubate the eppendorf in 🌡️37 °C , shaking for 🕒01:00:00
- 13 Spread the transformed bacteria on agar plate containing antibiotics.

14 Incubate the agar plate for 🕒 **Overnight (16-18 hours)**