


Oct 06, 2022


HMW DNA extraction for Long Read Sequencing using CTAB

Benoît Vacherie¹, Karine Labadie¹, Cyril Falentin²¹Genoscope, Institut de Biologie François Jacob, Commissariat à l'énergie atomique (CEA) Université de Paris-Saclay Evry 91057 France.;²IGEPP, INRAE, Institut Agro, Univ Rennes, Le Rheu, 35653, France.

2 Works for me

 Sharedx.doi.org/10.17504/protocols.io.bp2l694yzlqe/v1

HPM/Genoscope

 Benoît Vacherie
Genoscope/CEA

ABSTRACT

High Molecular Weight DNA extraction protocol for Long Read sequencing.

Extraction is performed from flash frozen leaves stored at -80°C.

The protocol is adapted for an extraction of 1g of leaves.

This protocol is based on the protocol provided by Oxford Nanopore Technologies, Oxford, UK (ONT), "High molecular weight gDNA extraction from plant leaves" provided by the ONT community in March 2019, with slightly modifications.

This protocol involves a conventional CTAB extraction followed by purification using commercial Qiagen Genomic tips (QIAGEN, MD, USA). DNA fragment size selection is performed using Short Read Eliminator (Circulomics, MD, USA).

This protocol is particularly adapted for plant leaves, but also works with many other organisms (microalgae, insects...).

DOI

dx.doi.org/10.17504/protocols.io.bp2l694yzlqe/v1

PROTOCOL CITATION

Benoît Vacherie, Karine Labadie, Cyril Falentin 2022. HMW DNA extraction for Long Read Sequencing using CTAB. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l694yzlqe/v1>



KEYWORDS

extraction, high molecular weight, plant, leaves, DNA, Long read sequencing

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 15, 2022

LAST MODIFIED

Oct 06, 2022

PROTOCOL INTEGER ID

70072

GUIDELINES

Use only wide bore tips.

Work in a chemical hood when using 2-mercaptoethanol.

Allow the DNA to resuspend for at least of 24 hours before proceeding with QC.

MATERIALS TEXT

Reagents :

[Tris HCl](#) **P212121**

[EDTA \(0.5 M\), pH 8.0](#) **Life**

Technologies Catalog #AM9260G

[CTAB \(Hexadecyltrimethylamm onium bromide\)](#) **BB**

Biotech Catalog #CB0108-100g

[Sodium chloride](#) **P212121**

[PEG-](#)

8000 Promega Catalog #V30111

[Liquid nitrogen](#) **Contributed by users** Step 3

[2-mercaptoethanol](#) **Sigma**

Aldrich Catalog #M6250 Step 2

[RNase](#)

A Qiagen Catalog #19101 Step 6

[Chloroform](#) **Contributed by users** Step 9

[Isopropanol](#) **Contributed by users** In 2 steps

 [Genomic tip](#)

100G Qiagen Catalog #10223 Step 19

 [Buffer](#)

G2 Qiagen Catalog #1014636 Step 17

 [Buffer](#)

QBT Qiagen Catalog #19054 Step 19

 [Buffer](#)

QC Qiagen Catalog #19055 Step 23

 [Buffer](#)

QF Qiagen Catalog #19056 In 2 steps

 [Ethanol 70%](#) **Contributed by users**

Consumables :

 [MBP™ Wide Bore Pipette Tips Thermo](#)

Fisher Catalog #02707600

 [Falcon 50mL Conical Centrifuge Tubes Fisher](#)

Scientific Catalog #14-959-49A

 [15 ml conical tubes](#) **Contributed by users**

 [1.5ml Eppendorf DNA LoBind tubes](#) **Contributed by users**

Equipment :

 [Porcelain Mortar, 145mL Thermo](#)

Fisher Catalog #CP1782004

 [Eppendorf ThermoMixer](#)

C pipette.com Catalog #2231000667 Step 7

 [Vortex Mixer](#) **Contributed by users**

 [Glass pasteur pipettes](#) **Contributed by users**

preparation of reagents

- 1 **Extraction Buffer** : Prepare 20 mL of extraction buffer (per 1g of leaves) in a 50 mL tube

Reagents	Final concentration
Tris-HCl pH8	100 mM
EDTA	20 mM
CTAB	2%
NaCl	1.4. M
PEG 8000	1 %
H2O	QSP 20 ml

DNA Extraction 3h

- 2 Add 50 µl of **2-Mercaptoethanol** to 20 ml of **extraction buffer** and preheat to 65 °C (approx. 15 min).

📄 20 mL Extraction buffer

🔗 2-mercaptoethanol Sigma

📄 50 µL Aldrich Catalog #M6250

🔥 65 °C

- 3 Put a **mortar in ice**.
Cool the mortar and pestle with **liquid nitrogen** until the bubbling stops.

🔗 Liquid nitrogen Contributed by users

- 4 **Grind 1g of frozen sample** to a fine powder (approx. 2 min), without adding liquid nitrogen.^{2m}

📄 1 g ⌚ 00:02:00



Before grinding



After grinding

5 **Transfer the powder** to the pre-warmed extraction buffer

6 **Add 40 µl of RNase A** (100mg/ml) and vortex the tube 5 sec.

5s

 [RNase](#)

 **40 µL** [A Qiagen Catalog #19101](#)

 **00:00:05 Vortex**

7 **Incubate 1h at 65°C** with intermittent agitation (300 rpm every 10 min)

1h

 **01:00:00 65°C**

 [Eppendorf ThermoMixer](#)

[C pipette.com Catalog #2231000667](#)

8 Let the tube cool for 5 min at RT

5m

 **00:05:00 RT**

9 **Add 1 volume (20ml) of chloroform** and vortex 2x5sec.

 **20 mL**  [Chloroform Contributed by users](#)

10 **Centrifuge** with low acceleration and deceleration

10m


 **5500 x g, 4°C, 00:10:00 , Acc/Dec : 6/3**



Before centrifugation



After centrifugation

- 11 Gently **transfer the aqueous phase** to a new 50 ml tube.
- 12 Add **0.7 volume of isopropanol** and mix gently by inversion 10 times
⊗ Isopropanol Contributed by users
- 13 Put the tube at **-80°C for 15 minutes.** 15m
🕒 00:15:00 -80°C
- 14 

*If a DNA medusa appears after this step, recover the medusa using a Pasteur pipette (without breaking the tip of the pipette).
Wash the medusa in 3 successive baths of 70% ethanol then resuspend the medusa in 100µl of 1X TE.
Incubate for 1 hour at 55°C then store the tube at 4°C before quality control.
Continue the protocol with the rest of the tube (without the medusa).
If no DNA medusa is visible continue the protocol.*
- 15 **Centrifuge** the tube with low acceleration and deceleration 30m
⚙️ 5500 x g, 4°C, 00:30:00 , Acc/Dec : 6/3
- 16 Carrefully **remove the supernatant** without resuspending the pellet.
Remove the remaining liquid by turning the tube upside down on a paper towel (make sure the pellet does not come off).
- 17 Gently **resuspend with the pipette the DNA pellet with 9.5 ml of G2 buffer** (QIAGEN Genomic-tip 100/G).
⊗ Buffer
G2 Qiagen Catalog #1014636
- 18 **Incubate 15 min at 50°C.** *The sample can be stored overnight at 4°C at this stage.* 15m
🕒 00:15:00 50°C

Génomique tip purification

1h

- 19 **Equilibrate a QIAGEN Genomic-tip 100/G** column with **4 ml of QBT buffer**.
[⌕ Genomic tip](#)
[100G Qiagen Catalog #10223](#)
[⌕ Buffer](#)
[QBT Qiagen Catalog #19054](#)
- 20 **Preheat 5 ml of QF buffer** to 55°C.
[⌕ Buffer](#)
[QF Qiagen Catalog #19056](#)
- 21 Allow all the QBT buffer to drain by gravity flow into a 50 ml tube.
- 22 **Add the sample** in G2 buffer (9.5ml) on the column and let it to enter the resin by gravity flow.
- 23 **Wash the column twice** with 7.5 ml of QC buffer.
[⌕ Buffer](#)
[QC Qiagen Catalog #19055](#)
- 24 Put the column on a 15 ml tube.
- 25 **Elute the DNA** with 5 ml of QF buffer preheated to 55°C.
[⌕ Buffer](#)
[QF Qiagen Catalog #19056](#)

DNA précipitation

3h

- 26 Add **0.7 volume of isopropanol** to the eluate.
Mix gently by inverting 20 times.
[⌕ Isopropanol Contributed by users](#)

- 27 Incubate 15 minutes at RT. 15m
🕒 00:15:00 RT
- 28 **Centrifuge** with low acceleration and deceleration 30m
🌀 5500 x g, 4°C, 00:30:00 , Acc/Dec : 6/3
- 29 Carefully discard the supernatant, gently resuspend the pellet with 1 ml of cold 70% ethanol.
📦 70% Ethanol Contributed by users
- 30 **Transfer the resuspended DNA** to a 1.5 ml DNA LoBind tube.
📦 DNA LoBind Tube 1.5ml
Eppendorf Catalog #022431021
- 31 **Centrifuge** with low acceleration and deceleration 10m
🌀 5000 x g, 4°C, 00:10:00 , Acc/Dec : 6/3
- 32 Carefully discard the supernatant. 10m
Air dry the pellet at RT for about 10 minutes.
🕒 00:10:00 RT
- 33 **Resuspend the pellet** with 50-100 µl of 1X TE buffer.
📦 Buffer TE 1x Contributed by users
- 34 **Allow the DNA to resuspend** for 2 hours at 55°C or ON at RT. 2h
🕒 02:00:00 55°C or 🕒 Overnight RT
- 35 **Store the DNA at 4°C.**

Sample quality control

- 36 Quantify your sample with a **Qubit DNA HS assay**.











Check the purity of the sample with a **Nanodrop** (measurements of 260/280 and 260/230 absorbance ratios).

Estimate the molecular weight of the sample with a **Tapestation** and/or a **Femto pulse** and/or a **Pippin Pulse**.

- 37 Depending on the DNA concentrations and DNA length profiles, deplete short DNA molecules using **SRE size selection kits** (SRE XS, SRE or SRE XL kits).

Results

- 38 **QC results** obtained on different species.

species	DNA quantity (µg)	Extraction yield (µg DNA / g organism)	DNA length (kb) estimation with				Nanopore N50 reads
			Nanodrop 260/280	Nanodrop 260/230	Tapestation	Femto pulse	
<i>Arabidopsis thaliana</i>	14,8	24,7	1,97	1,73	> 60 kb 	72 kb 	30 kb
<i>Lactuca sativa</i>	38,5	38,5	1,96	2,30	> 60 kb 	149 kb 	43 kb
<i>Brassica napus</i>	169	181	1,88	1,92	> 60 kb 		45 kb
<i>Fagus sylvatica</i>	33,8	254	1,84	2,07	> 60 kb 	100 kb 	30 kb
<i>Quercus robur</i>	98	183	1,98	1,64	> 60 kb 		48 kb
<i>Ectocarpus sp.</i>	3,2	70	1,88	1,81	> 60 kb 		22 kb
<i>Tenebrio molitor</i>	22,4	370	1,88	2,17	> 60 kb 		35 kb