

Purification of
Human K560-GFP
molecular motor

May 23, 2022

Purification of Human K560-GFP molecular motor

Andrea M. Dickey¹

¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093

1



dx.doi.org/10.17504/protocols.io.bp2l61xrdvqe/v1



Mariusz Matyszewski

K560-GFP purification protocol from the Reck-Peterson Lab based on protocol from Nicholas et al. 2014. Edited for protocols.io by Andrea Dickey and Mariusz Matyszewski.

DOI

dx.doi.org/10.17504/protocols.io.bp2l61xrdvqe/v1

Andrea M. Dickey 2022. Purification of Human K560-GFP molecular motor.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bp2l61xrdvqe/v1>



single molecule, kinesin, purification, microtubule

protocol ,

Mar 03, 2022

May 23, 2022

59048

All protein purification steps should be performed at **4 °C** unless noted otherwise.

Materials:

cOmplete EDTA-free protease inhibitor cocktail tablets

lysozyme

PMSF

Buffers:

Lysis buffer:

- [M]50 millimolar (mM) Tris pH 7.5
- [M]300 millimolar (mM) NaCl
- [M]5 millimolar (mM) MgCl₂
- [M]0.2 Molarity (M) sucrose
- [M]1 millimolar (mM) DTT
- [M]0.1 millimolar (mM) Mg-ATP
- [M]0.5 millimolar (mM) Pefabloc

Wash buffer:

- [M]50 millimolar (mM) Tris pH 7.5
- [M]300 millimolar (mM) NaCl
- [M]5 millimolar (mM) MgCl₂
- [M]0.2 Molarity (M) sucrose
- [M]10 millimolar (mM) imidazole
- [M]1 millimolar (mM) DTT
- [M]0.1 millimolar (mM) Mg-ATP
- [M]0.5 millimolar (mM) Pefabloc

Elution buffer:

- [M]50 millimolar (mM) Tris pH 8.0
- [M]300 millimolar (mM) NaCl
- [M]5 millimolar (mM) MgCl₂
- [M]0.2 Molarity (M) sucrose
- [M]250 millimolar (mM) imidazole
- [M]5 millimolar (mM) bME
- [M]0.1 millimolar (mM) Mg-ATP

Storage buffer:

- [M]80 millimolar (mM) PIPES pH 7.0
- [M]2 millimolar (mM) MgCl₂
- [M]1 millimolar (mM) EGTA
- [M]0.2 Molarity (M) sucrose
- [M]1 millimolar (mM) DTT
- [M]0.1 millimolar (mM) Mg-ATP




Glycerol cushion:

- [M]80 millimolar (mM) PIPES pH 7.0
- [M]2 millimolar (mM) MgCl₂
- [M]1 millimolar (mM) EGTA
- [M]60 % volume glycerol
- [M]20 micromolar (μM) Taxol
- [M]1 millimolar (mM) DTT

BRB80:

- [M]80 millimolar (mM) PIPES pH 7.0
- [M]2 millimolar (mM) MgCl₂
- [M]1 millimolar (mM) EGTA
- [M]300 millimolar (mM) KCl
- [M]7.5 millimolar (mM) Mg-ATP














Expression

- 1 pET17b-Kif5b(1-560)-GFP-His should be transformed into BL21(DE3)RIPL cells.
- 2 Make enough LB for at least  7.5 L of culture.
- 3 Grow an overnight starter culture.
- 4 Transfer starter culture into LB. Make sure to add proper antibiotics (Ampicillin for the plasmid, and Chloramphenicol for the cells we used).
Allow it to grow at  200 rpm, 37°C until OD₆₀₀ reaches 0.6-0.8.
- 5 Chill cells and incubator to  18 °C
- 6 Add [M]0.5 millimolar (mM) IPTG and allow it to grow at  200 rpm, 18°C, 18:00:00

7 Harvest and freeze cell pellets.

Purification

2h 35m

- 8 Resuspend  **7.5 L** worth of pellets in  **120 mL Lysis buffer** supplemented with 1 30m
cOmplete EDTA-free protease inhibitor cocktail tablet (Roche) per  **50 mL of Lysis buffer** .
Also add  **1 mg/mL lysozyme** .
Incubate  **On ice** for  **00:30:00** .
- 9 Lyse the resuspension by sonication.
- 10 Add  **0.5 millimolar (mM) PMSF** to the sonicate, and clarify by centrifugation 1h
 **40000 rcf, 4°C, 01:00:00** in Type 70 Ti rotor (Beckman).
- 11 Incubate the supernatant with Ni-NTA agarose beads incubated with the **Wash buffer**. 1h
Nutate for  **01:00:00**
- 12 Apply to gravity column and wash with  **100 mL Wash buffer** .
- 13 Resuspend the beads in elution buffer, incubate for  **00:05:00** and elute in  **0.5 mL** 5m
increments.
- 14 Combine peak fractions and buffer exchange on PD-10 desalting column equilibrated with **Storage buffer**.
- 15 Peak fractions of the motor solution were either flash-frozen at  **-80 °C** until further use or immediately subjected to microtubule bind and release purification (see next steps).

- 16 A total of **1 mL** of motor solution from previous steps is used. Incubate with **1 millimolar (mM) AMP-PNP** and **20 micromolar (μM) Taxol** **On ice** in the **dark** for **00:05:00** and subsequently warm to **Room temperature** . 5m
- 17 Polymerize bovine brain tubulin (about **00:30:00** at **37 °C**) and centrifuge (TLA120.2 rotor at **80000 rpm, 00:12:00** at **Room temperature**) through a **glycerol cushion**. Resuspend the pellet. 42m
- 18 Incubate the microtubule solution with the resuspended microtubules in the **dark** for **00:15:00** at **Room temperature** . 15m
- 19 Put the incubated solution on top of **glycerol cushion** and centrifuge in a TLA120.2 rotor at **80000 rpm, 00:12:00** at **Room temperature** . 12m
- 20 Wash the final pellet with **BRB80** and incubate in **100 μL of release buffer** for **00:05:00** at **Room temperature** . 5m
- 21 Centrifuge the kinesin solution **72000 rpm, 00:07:00** in TLA100 rotor at **Room temperature** . 7m
- 22 Collect the supernatant and supplement with **660 millimolar (mM) sucrose** and flash freeze.

A typical kinesin prep in the lab yielded **0.5 micromolar (μM)** to **1.5 micromolar (μM)** K560-GFP dimer.