



VERSION 3

FEB 23, 2024

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v3](https://dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v3)

**Protocol Citation:** Rebecca Andrews, Joanne Lachica, Steven F. Lee, Sonia Ghandi 2024. Single-molecule Immunofluorescence Tissue Staining Protocol for Oligomer Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v3> Version created by Joseph S Beckwith

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

# Single-molecule Immunofluorescence Tissue Staining Protocol for Oligomer Imaging V.3

Rebecca Andrews<sup>1</sup>, Joanne Lachica<sup>2</sup>, Steven F. Lee<sup>1</sup>, Sonia Ghandi<sup>2</sup>

<sup>1</sup>University of Cambridge; <sup>2</sup>UCL

ASAP Collaborative Research Network

Becs' Workspace



Rebecca Andrews

Univeristy of Cambridge

## ABSTRACT

This protocol details background fluorescence quenching and immunofluorescence staining of human brain tissue for oligomer imaging.

## ATTACHMENTS

[kb3ib25np.pdf](#)

## GUIDELINES

- Use only clean bottles, flasks, magnetic stirrers, tweezers, weighing spatulas, measuring cylinders – everything should be cleaned, dried and covered if left on the side before next use.
- Everything should be handled with clean tweezers – gloves should not touch the samples, solutions and ideally anything placed into the solutions where the slides are.

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 23, 2024  
**Last Modified:** Feb 23, 2024  
**PROTOCOL integer ID:** 95660

**Keywords:**  
immunofluorescences, oligomer imaging


MATERIALS


**Materials and Reagents**

- Microtome
- Glass slides
- Xylene solution
- 100% alcohol
- Methanol
- Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) solution
- Citrate buffer pH6
- Milli Q water
- Pressure cooker
- PBS
- Goat Serum 10%
- AlexaFluor antibody
- 0.1% Sudan black solution
- Vectashield
- Overslip



Immunofluorescences staining protocol for oligomer imaging

10m

1 Cut  8 µm tissue sections on a microtome and load onto glass slides.

2 Dry slides  Overnight at  37 °C – **cover over the top.**

10m

3 Before staining commences keep slides for a few hours but ideally  Overnight at  60 °C .

10m

4 De-wax sections through three pots of xylene solution. Use each fresh pots of xylene each time.

4.1 De-wax section through pot of xylene solution for 00:02:00 . (1/3)

2m

4.2 De-wax section through pot of xylene solution for 00:02:00 . (/23)

2m

4.3 De-wax section through pot of xylene solution for 00:02:00 . (3/3)

2m

5 Take sections through two pots of 100% alcohol.

Note

Use fresh pots each time – methylated spirits.

5.1 Take sections through pot of 100% alcohol for 00:02:00 . (1/2)

2m

5.2 Take sections through pot of 100% alcohol for 00:02:00 . (2/2)

2m

6 Put slides into methanol + hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) solution (100 ml: 1 ml) for 00:10:00 at

10m

Room temperature **fresh pot each time.**

#### Note


This process will block any staining of endogenous peroxidase in the tissue sections.



7 Perform necessary antigen retrieval pretreatments by pressure cooking in citrate buffer.

7.1 Pressure cook sections in citrate buffer  6 for  00:10:00 at pressure (wait for it to release high pressure air) in **cleaned pressure cooker**. 10m

7.2 Cool under running **Milli Q** water – *never* under tap water.





8 Block non-specific antigen/antibody binding by placing sections in **PBS and Goat Serum 10%** for  00:30:00 . **NB Goat Serum is chosen for Goat-Raised antibodies**. 30m

9 Apply primary antibody for  01:00:00 at  Room temperature . 1h


10 Wash in PBS with fresh buffer (**at least filtered if not cell culture grade**).



10.1 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (1/3) 5m

10.2 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (2/3)

5m

10.3 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (3/3)


5m

11 Apply secondary AlexaFluor antibody for  01:00:00 at  Room temperature in the dark.


1h

12 Wash in PBS.




12.1 Wash  00:05:00 in PBS in the dark. (1/3)



5m

12.2 Wash  00:05:00 in PBS in the dark. (2/3)

5m

12.3 Wash  00:05:00 in PBS in the dark. (3/3)

5m

13 Add **filtered (0.22 um)** 0.1% Sudan black solution (0.1% sudan black/70% ethanol) for  00:10:00 at  Room temperature in the dark.

10m



**14** Wash 2-3 times in 30% ethanol.



**15** Mount section with Vectashield and coverslip (**Plasma cleaned slides**).

**16** Take for imaging.

