



Dec 22, 2021

Total Starch Enzymatic Digestion

Lynn Doran¹, Amanda P. De Souza¹

¹Realizing Increased Photosynthetic Efficiency (RIPE)

1



dx.doi.org/10.17504/protocols.io.b26tqhen

Burgess Lab UIUC



Lynn Doran

Realizing Increased Photosynthetic Efficiency (RIPE)

Enzymatic digestion of total soluble starch to glucose in plant tissue extracts for preparation for quantification via the GOD-POD Method (NZYtech).

DOI

dx.doi.org/10.17504/protocols.io.b26tqhen

Lynn Doran, Amanda P. De Souza 2021. Total Starch Enzymatic Digestion .

protocols.io

<https://dx.doi.org/10.17504/protocols.io.b26tqhen>



Realizing Increased Photosynthetic Efficiency (RIPE) that is funded by the Bill & Melinda Gates Foundation, Foundation for Food and Agriculture Research, and the U.K. Foreign, Commonwealth & Development Office

Grant ID: OPP1172157

protocol

Pak S. Chow, Simon M. Landhäusser, A method for routine measurements of total sugar and starch content in woody plant tissues, Tree Physiology, Volume 24, Issue 10, October 2004, Pages 1129–1136, <https://doi.org/10.1093/treephys/24.10.1129>

Starch, Soluble Starch, Enzymatic Digestion, Plant Tissue, GOD-POD

protocol ,

Dec 21, 2021

Dec 22, 2021

56243

Reagents

- α -amylase enzyme, 1000 U/mL
- Amyloglucosidase enzyme, 3260 U/mL
- MOPS Buffer 10mM, pH 6.5

2.09 g of MOPS per 1 L of water (0.4598 g for 220 mL water), pH adjusted to 6.5 using NaOH.

Store refrigerated for up to 2 months.

- Acetate Buffer 100mM, pH 4.5

6 mL of acetic acid per 1 L of water (1.5 mL per 250 mL water), pH adjusted to 6.5 using NaOH.

Store refrigerated for up to 2 months.

- Ice

Materials

- Pipette tips

Equipment

- Graduated cylinder
- Water bath
- Floating tube holder
- Single channel pipette
- Ice bucket

Extract and dry total starch pellet from plant tissue per [Extraction of Non-Structural Carbohydrates \(Total Soluble Sugars + Starch\) in Plant Tissues](#).

- 1 Prepare fresh daily 120 U/mL α -amylase in MOPS buffer. 1 mL per sample will be needed. Initial concentration of α -amylase is 1000 U/mL. Use $C_1V_1 = C_2V_2$ to calculate the volume of α -amylase and MOPS buffer to use.

α -amylase comes at 1000 U/mL concentration in a 20 mL bottle. If using all of the bottle, 20 mL of α -amylase in 146.7 mL of MOPS will result in 120 U/mL α -amylase in MOPS buffer.





- 2 Prepare fresh daily 30 U/mL amyloglucosidase in acetate buffer. 1 mL per sample will be needed. Initial concentration of amyloglucosidase is 3260 U/mL. Use $C_1V_1 = C_2V_2$ to calculate

the volume of amyloglucosidase and acetate buffer to use.

To prepare 250 mL of 30 U/mL amyloglucosidase in acetate buffer, dilute 2.3 mL amyloglucosidase (3260 U/mL) into 247.7 mL acetate buffer.

- 3 Heat the water bath to **75 °C**.
- 4 Add **500 µL** of α-amylase in MOPS buffer (120 U/mL) to each sample tube.
- 5 Vortex to suspend all solids. Flicking the tube may help dislodge the dried starch pellet.
- 6 Incubate for **00:30:00** at **75 °C** in the water bath. 30m
- 7 Add another **500 µL** of α-amylase in MOPS buffer (120 U/mL) to each sample tube.
- 8 Vortex to suspend all solids.
- 9 Incubate for **00:30:00** at **75 °C** in the water bath. 30m
- 10 Cool down the water bath to **50 °C**, tubes can be stored at room temp on the counter while water bath cools.

Leaving the lid off, especially if it is a shaking water bath, will help cool. For faster cooling, remove some water from the water bath and refill with cool distilled water. Do not use tap water as mineral buildup or heated chlorine could damage the water bath.

-
- 11 Add  500 µL amyloglucosidase in acetate buffer (30 U/mL) to each sample tube.
 - 12 Vortex to suspend all solids.
 - 13 Incubate for  00:30:00 at  50 °C in the water bath. 30m
 - 14 Add another  500 µL amyloglucosidase in acetate buffer (30 U/mL) to each sample tube.
 - 15 Vortex to suspend all solids.
 - 16 Place the tubes in ice to stop the reaction, until cool to the touch.
 - 17 Proceed to total starch (as glucose) quantification by NZYtech GOD-POD method or store the samples at -20°C up to one month.