




JAN 22, 2024

S-Basal Medium

 In 1 collection

Adrien Assie¹, Buck Samuel¹

¹Baylor College of Medicine

Adrien Assie: We are not the author of this protocol. It was first described in Lewis and Fleming (1995). Protocol available on wormbook.

Samuel Lab

OPEN  ACCESS



Adrien Assie
Baylor College of Medicine

DISCLAIMER

We are not the author of this protocol. It was first described in Lewis and Fleming (1995). Protocol inspired from the one available on the Wormbook website.

ABSTRACT

Large quantities of *C. elegans* can be grown in a liquid medium. Liquid cultures of *C. elegans* are usually grown on S Medium using concentrated *E. coli* OP50 as a food source. This is the S-Basal Medium recipe

External link:

http://www.wormbook.org/chapter_s/www_strainmaintain/strainmaintain.html

Protocol Citation: Adrien Assie, Buck Samuel 2024. S-Basal Medium. **protocols.io**
<https://protocols.io/view/s-basal-medium-zy9f7z6>

MANUSCRIPT CITATION:

Lewis, J.A. and Fleming, J.T. (1995). In: Methods in cell biology, Vol. 48, H.F. Epstein and D.C. Shakes, eds. (San Diego: Academic Press), p. 3.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited









Protocol status: Working

Created: Apr 10, 2019

Last Modified: Jan 22, 2024

PROTOCOL integer ID: 22273

Keywords: C. elegans, S medium, liquid culture

- 1 Start with  700 mL water
- 2  5.9 g NaCl ( 100 millimolar (mM))
- 3  50 mL of  1 Molarity (M) Potassium Phosphate Buffer, pH 6.0
OR
 1 g K₂HPO₄ and  6 g KH₂PO₄
- 4 Adjust water to  1000 mL
- 5 Autoclave