



# RCA of Gotcha V.1

Chia-Hsien Shih<sup>1</sup>

<sup>1</sup>Chung Shan Medical University

1



1 ▼

Oct 17, 2021

[dx.doi.org/10.17504/protocols.io.by5hpy36](https://dx.doi.org/10.17504/protocols.io.by5hpy36)

CSMU\_Taiwan



Chia-Hsien Shih  
Chung Shan Medical University

This protocol is to test that GotCha is functional and can works as designed.

DOI

[dx.doi.org/10.17504/protocols.io.by5hpy36](https://dx.doi.org/10.17504/protocols.io.by5hpy36)

Chia-Hsien Shih 2021. RCA of Gotcha. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.by5hpy36>



protocol ,

Oct 17, 2021

Oct 17, 2021









54153

## Preparation

- 1 Add **5  $\mu$ L** of GotCha(functional beads) into eppendorf
- 2 Centrifuge for **15000 rpm, 00:05:00** and remove supernatant. Make sure that eppendorf should put on DynaMag when removing supernatant. <sup>5m</sup>


## Protocol

- 3 Add **3  $\mu$ L** of 10X phi29 polymerase reaction buffer into eppendorf with GotCha

- 4 Add  **20.4 µL** RNase-free water
- 5 Add  **3 µL** of 100nM miRNA
- 6 Add  **3 µL** of 2mM dNTPs
- 7 Add  **0.6 µL** of 10U/ul phi29 polymerase
- 8 Pipetting to mix well
- 9 Incubate for  **02:00:00** at  **Room temperature**
- 10 Add  **1.5 µL** of 20X evagreen dye
- 11 Add  **1.5 µL** of 0.5M EDTA buffer to suspend the reaction

2h

#### Measuring

- 12 Load  **20 µL** of reaction solution into 384-well plate
- 13 Measure the fluorescence excitation and emission intensity

