

Jun 11, 2024

## Filter trap assay for the detection of alpha-synuclein aggregation

DOI

**[dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1)**

Cole S Sitron<sup>1</sup>, Victoria A Trinkaus<sup>1</sup>, F Ulrich Hartl<sup>1</sup>

<sup>1</sup>Department of Cellular Biochemistry, Max Planck Institute of Biochemistry

ASAP Collaborative Rese...



**Cole Sitron**

Max Planck Institute of Biochemistry

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1)**

**Protocol Citation:** Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Filter trap assay for the detection of alpha-synuclein aggregation. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 17, 2024

**Last Modified:** June 11, 2024

**Protocol Integer ID:** 101394

**Funders Acknowledgement:**  
**Aligning Science Across**  
**Parkinson's**  
**Grant ID:** ASAP-000282

## Abstract

This protocol describes an assay that detects aggregated alpha-synuclein in cell lysates.

## Materials



### Cell culture

Cells with alpha-synuclein aggregates, induced as in [[dx.doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1](https://doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1)], growing in 12-well plates.

### Reagents

- TrypLE™ Express **Thermo Fisher Scientific Catalog #12605010**
- Pierce&trade; Rapid Gold BCA Protein Assay Kit **Thermo Fisher Catalog #A53225**

### Buffers


- 10% FBS (  Fetal Bovine Serum **Gibco - Thermo Fischer Catalog #10270106** ) in 1X PBS (diluted from  PBS (10X), pH 7.4 **Thermo Fisher Scientific Catalog #70011051** )
- 1X PBS
- Lysis Buffer:

A	B
Tris-HCl pH 7.5	50 mM
NaCl	150 mM
Triton-x100	1% (v/v)
cOmplete protease inhibitor tablet	1 mini tablet per 15 ml buffer (Sigma-Aldrich cat. no. 4693159001)
Phos-STOP phosphatase inhibitor tablet	1 tablet per 15 ml buffer (Sigma-Aldrich cat. no. 4906845001)
Benzonase ((Max Planck Institute of Biochemistry Core Facility)	7.5 U/ml

- Wash Buffer:

A	B
Tris-HCl pH 7.5	50 mM
NaCl	150 mM
Triton-x100	1% (v/v)

### Consumables

- 0.2 micrometer pore size cellulose acetate membrane (e.g. GE cat. no. 10404131)
-  PARAFILM® M Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7793
































**Equipment**

Equipment	
Bioruptor® Plus sonication device	NAME
Sonication device	TYPE
Bioruptor®	BRAND
B01020001	SKU
<a href="https://www.diagenode.com/en/p/bioruptor-plus-sonication-device">https://www.diagenode.com/en/p/bioruptor-plus-sonication-device</a>	LINK

- Slot blot vacuum manifold (e.g. Hoefer PR648)
- Fume hood



## Harvesting cells

- 1 Prepare lysis buffer and place  On ice .
- 2 After  24:00:00 of exposure to PFFs/PBS, wash wells 1X with  1 mL warm PBS.   

- 3 Withdraw PBS and treat wells with  300  $\mu$ L of TrypLE. 
- 4 Quench dissociation reagent with  300  $\mu$ L of 10% FBS and transfer the contents of the wells to 1.5 ml tubes. 
- 5 Centrifuge  1500 x g for  00:03:00 .   

- 6 Withdraw the supernatant and wash cells with  1 mL PBS. 
- 7 Centrifuge  1500 x g for  00:03:00 .   

- 8 Withdraw supernatant and resuspend in  100  $\mu$ L of lysis buffer and incubate  On ice for  00:05:00 .   

- 9 Sonicate the cells in the Bioruptor for 3 cycles of  00:00:30 on and  00:00:30 off at high power. 
- 10 Centrifuge the lysate at  500 x g for  00:05:00 and place the clarified lysate in a new 1.5 ml tube.   


**Note**


This step is crucial to ensure the removal of high-diameter material (e.g. unlysed cells) that may clog the filter paper and non-specifically trap proteins on it.

**Note**



If desired, this is a convenient stopping point, wherein you can snap-freeze your lysate in liquid nitrogen and proceed with the rest of the assay at a later time.

- 11 Determine the concentration of protein in each sample by Pierce Rapid Gold BCA Protein Assay.




**Filter Trap****10m**

- 12 Soak the cellulose acetate membrane in wash buffer for >  00:10:00 .

**10m**



- 13 For each well that you run, prepare  1020  $\mu\text{L}$  of the sample containing 10.2 - 51 ug of lysate, with the remaining volume consisting of lysis buffer. 

- 14 Assemble the slot blot vacuum manifold apparatus with the cellulose acetate membrane.


- 15 Wash each well with  300  $\mu\text{L}$  wash buffer, inspecting each well to ensure that the buffer can drain through it.  

**Note**

Avoid using any non-functioning wells.

- 16 Load  1000  $\mu\text{L}$  (10 - 50 ug) of sample to each well. 



- 17 Once the sample has drained through the wells, wash each well 3 times with  300  $\mu$ L of wash buffer.
- 18 When all of the washes have drained through the wells, trim the membrane and proceed with standard immunodetection methods.

