



Apr 27, 2021

# Cryopreservation of microalgal cultures

Victoria Jackson<sup>1</sup><sup>1</sup>Living Systems Institute, University of Exeter**1** Works for me This protocol is published without a DOI.Victoria Jackson  
Living Systems Institute, University of Exeter

## ABSTRACT

Protocol for long-term cryopreservation of microalgal cell cultures.

## PROTOCOL CITATION

Victoria Jackson 2021. Cryopreservation of microalgal cultures. **protocols.io**  
<https://protocols.io/view/cryopreservation-of-microalgal-cultures-buixnufn>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Apr 27, 2021

## LAST MODIFIED

Apr 27, 2021

## PROTOCOL INTEGER ID

49463

### Preparing the cryoprotectant 10m

- 1 Make a **10 % (v/v)** solution of **DMSO** in the appropriate culture medium for the strains being preserved.

Eg. **45 mL culture medium** medium and **5 mL DMSO** to make **50 mL** cryoprotectant.

Work in the fume hood when using DMSO.

- 2 Filter sterilise the cryoprotectant using a **0.2 µM filter**.

### Preparing the cultures 15m


- 3 Prepare the cryotubes with a 1:1 ratio of cryoprotectant:culture

Eg. **1 mL cryoprotectant** and **1 mL culture** for a 2 mL cryotube.

- 4 Once the cryoprotectant and culture have been mixed in the cryotubes, allow **00:10:00** for the cryoprotectant to penetrate. <sup>10m</sup>




## Progressive freezing

16h

- 5 Place the cryotubes into a Mr Frosty freezing container (filled to the line with isopropyl alcohol/isopropanol/propan-2-ol).
- 6 Place the freezing container in a  $-80\text{ }^{\circ}\text{C}$  freezer  Overnight .
- 7 Remove the tubes from the freezer and flash freeze them in liquid nitrogen.
- 8 Store the tubes either in liquid nitrogen or in a  $-150\text{ }^{\circ}\text{C}$  freezer for as long as required.

## Culture recovery

2w

- 9 To recover a cryopreserved culture, thaw the cryotubes in a water bath pre-heated to  $40\text{ }^{\circ}\text{C}$  to prevent ice crystals from re-forming as they thaw.
- 10 Add the defrosted culture to  20 mL fresh culture medium .
- 11 Place the culture at the appropriate temperature but **keep it in the dark** for  24:00:00 , for example by wrapping<sup>1d</sup> the culture flask in aluminium foil.
- 12 After  24:00:00 , remove the aluminium foil. Monitor the culture to determine if recovery has been successful.<sup>1d</sup>