



# Total lactate dehydrogenase calorimetric enzyme activity measurement V.1

Renuka Sriram<sup>1</sup>

<sup>1</sup>University of California, San Francisco

Version 1

Jun 09, 2021

1

Works for me




Share

[dx.doi.org/10.17504/protocols.io.bvhmn346](https://dx.doi.org/10.17504/protocols.io.bvhmn346)

Renuka Sriram

## ABSTRACT

Lactate Dehydrogenase is a key cytoplasmic enzyme that confers the observed Warburg effect in tumors. Its a tetrameric protein composed of combinations of one or two isomers. This protein catalyzes the reversible reaction



Source: Wikipedia

This reaction is monitored spectrophotometrically by observing the reduction in absorbance of NADH (which is maximal at 420nm) as pyruvate is converted to lactate.

This protocol, details the steps involved in such a calorimetric assay of cells and tissues.

## DOI

[dx.doi.org/10.17504/protocols.io.bvhmn346](https://dx.doi.org/10.17504/protocols.io.bvhmn346)

## PROTOCOL CITATION

Renuka Sriram 2021. Total lactate dehydrogenase calorimetric enzyme activity measurement.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bvhmn346>

## KEYWORDS

Lactate dehydrogenase, calorimetric, lactate, pyruvate, NADH

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 02, 2021

## LAST MODIFIED

Jun 09, 2021

## PROTOCOL INTEGER ID

50445

## MATERIALS TEXT

Material	Source	Catalogue No.	Storage
<b>Cell Lysis Buffer (10X)</b>	Cell Signaling Technology	9803S	-20°C
<b>NaCl</b>	ThermoFisher	BP358-1	Room temperature
<b>NADH</b> ( $\beta$ -Nicotinamide adenine dinucleotide, reduced disodium salt)	Sigma	N6785-10VL	-20°C
<b>Sodium pyruvate</b> (11mg/mL, 100x)	Sigma	S8636	2-8°C
<b>Tris Base</b> (Trizma Base)	Sigma	93362(BioUltra, for molecular biology, $\geq 99.8\%$ (T)) or 93352 ( $\geq 99.0\%$ (T))	Room temperature
<b>Tris-HCl</b> (Trizma hydrochloride)	Sigma	T6666-50G	Room temperature
DTT (1,4-Dithiothreitol)	Sigma	10197777001	2-8°C
EDTA (Ethylenediaminetetraacetic acid)	Sigma	1233508	
Bradford Assay Kit	ThermoFisher	23200	
Triton X-100	Thermofisher	28314	room temperature

## Apparatus



1. If working with cells - Cell counter

We use [Bio-Rad TC20™ Automated Cell Counter](#)



or hemocytometer

hemocytometer

## 2. Stator-rotor homogenizer



or bead lyser



We use Qiagen TissueLyser LT



3. Plate reader

We use Tecan M200 Infinite

4. Multichannel pipette (20-200ul)
5. Micropipettes (10 - 1000 ul)
6. Stir plate
7. Water bath (@30C)
8. Centrifuge
9. Vortex mixer

#### Other consumables



1. 12 well reservoir
2. 96 well flat bottom clear plate
3. Pipette tips
4. 15 ml conical tubes
5. Eppendorf tubes (1.2 - 2 ml)

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Pre-preparation

- 1 Reaction buffer – Tris buffer @ **pH7.2**, **[M]200 Milimolar (mM)** **☒Sodium chloride P212121**,  
**[M]200 Micromolar (μM)**  
**☒β-Nicotinamide adenine dinucleotide reduced disodium salt hydrate (NADH) Sigma Aldrich Catalog #N8129**
  - 1.1 In a **☒200 mL** beaker, measure and add **☒140 mL** **☒dH2O Contributed by users**
  - 1.2 Place magnetic stir bar in beaker, place onto stir plate and set to **☒200 rpm**
  - 1.3 Measure and add the following:
    - i. **☒1.57 g** HCl Tris
    - ii. **☒0.15 g** Tris Base
    - iii. **☒1.64 g** of NaCl
    - iv. Use solution to dissolve NADH and add **☒20 mg** (10 vials of **☒2 mg** each)
  - 1.4 Turn off stir plate and remove stir bar
  - 1.5 Measure and record buffer pH (should be **pH7.2**, adjust with appropriate base/acid)
  - 1.6 Store in **☒-20 °C** freezer as **☒12 mL** aliquots (for 12 dilutions of pyruvate). This solution will be stable for upto 6 months.
- 2 Lysis buffer – **[M]50 Molarity (m)** Tris buffer **pH8.2**, **[M]2 Milimolar (mM)** DTT, **[M]2 Milimolar (mM)** EDTA, 1% Triton x-100.

- 2.1 In a **50 mL** conical tube, measure and add **20 mL** of dH<sub>2</sub>O
- 2.2 Measure and add the following:
- 0.71 g** HCl Tris
  - 0.67 g** Tris Base
  - 0.15 g** EDTA
  - 0.06 g** (or 59.326  $\mu$ L) of DTT
  - 2 mL** Triton
- 2.3 Store as **200  $\mu$ L** aliquots in **-20 °C** freezer (is stable for up to a year). This is 10 times concentrated, so for the assay add **1800  $\mu$ L** dH<sub>2</sub>O.
- 3 Store **500  $\mu$ L** aliquots of sodium pyruvate ( **11 mg/mL** ) in the **-20 °C** freezer. This solution is stable for up to a year when frozen.

#### Set-up

- 4 Take from the freezer, one tube each of the pyruvate and reaction buffer and place in **30 °C** water bath. Thaw the lysis buffer and add **1800  $\mu$ L** of dH<sub>2</sub>O and place on ice
- 5 Prepare 12 pyruvate concentrations using pyruvate and reaction buffer by carefully pipeting into pre-labeled conical tubes

- 5.1 Example of dilutions (can be subject to change based on sample type and amount of buffer needed):

Pyruvate Concentration (mM)	Vol. of Pyruvate stock sol'n to be added ( $\mu$ L)	Vol of Reaction Buffer to be added ( $\mu$ L)
0.05	1.25	2498.75
0.1	2.5	2497.5
0.2	5	2495
0.3	7.5	2492.5
0.4	10	2490
0.5	12.5	2487.5
0.6	15	2485
0.7	17.5	2482.5
0.8	20	2480
0.9	22.5	2477.5
1	25	2475
2	50	2450
TOTAL	188.75	29811.25

Pyruvate Dilution table

- 5.2 Calculate the amount of pyruvate + reaction buffer dilution needed for the entire assay -  $147 \mu\text{l} \times 3$  (triplicate)  $\times$  (# samples)
- 5.3 Pipette calculated amount +  $150 \mu\text{l}$  (to ensure enough is available) into each corresponding 12 well multichannel reservoir- from lowest to highest pyruvate concentration.
- 5.4 Keep the multi-well reservoir in the  $30^\circ\text{C}$  water bath

10m

- 6 TECAN spectrophotometer on and set temperature at  $30^\circ\text{C}$  .
  - a. Ensure the sample plate map of the appropriate method is set for the assay.
  - b. Set-up the protocol for a dynamic  $00:10:00$  readout at 420 nm and for the appropriate 96 well plate
- 7 Set-up the bead homogenizer in the cold room
- 8 Set -up buckets of ice and dry-ice as needed

#### Cell preparation

- 9 Trypsinize and collect cells:
  - 9.1 Aspirate the medium from cell culture vessel, and rinse with PBS. Aspirate the PBS, and add 0.25% trypsin in PBS.
  - 9.2 Place the culture vessel in the incubator (4-10 minutes based on cell type), and after the cells detach, add media (containing serum to inactivate the trypsin).
  - 9.3 Transfer the cells to a  $15 \text{ mL}$  tube and keep some in an eppendorf tube for cell counting.
  - 9.4 Centrifuge  $15 \text{ mL}$  tubes at  $300 \times g$  for  $00:05:00$  at  $4^\circ\text{C}$  .

5m


#### 10 Cell Counting

If using a cell counter, use the average of three measurements to determine the number of cells.



10.1

10.2 If using a hemocytometer (count the cells in 1,2,3 and 4 squares in both chambers and take the average). cells/ml = average count per square  $\times 10^4$



Hemocytometer

10.3 Total cells = cells/ml  $\times$  dilution (if used)  $\times$  total original volume of cell suspension from which sample was taken)

10.4 Re-suspend the cells in lysis buffer to obtain a concentration of  $10^7$  cells/ml

#### Tissue homogenization




- 11 Weigh and record tissue sample ( **6-10 mg** is sufficient) carefully using sterile forceps (keep frozen).
- a. Record any additional information about the sample (sample name/description, locus of tissue taken, date, etc.).
  - b. Place tissue into a screw cap vial (compatible with the bead homogenizer) and place on ice.

- 12 Add **500  $\mu$ l** of lysis buffer into the vial with tissue and place on ice

- 13 Two options for homogenizing tissue
- Option 1 - Using bead homogenizer -
- a. Place a bead (3-5 mm stainless steel) into the bead homogenizer compatible vial (from step 11b)
  - b. Operate the machine for **00:02:00** at 20-30 Hz

2m 50s

Option 2 - Using stator-rotor -

- a. Label plastic tubes EtOH and 0.1% DMPC (dimethyl pyrocarbonate) , one set for each sample) for cleaning.
- b. Place homogenizer into EtOH for  00:00:10 ; wipe excess with Kimwipe.
- c. Place homogenizer into DMPC H2O for  00:00:10 and wipe excess with Kimwipe.
- d. Use homogenizer to grind tissue at medium to high settings for approximately  00:00:30 or until tissue is fully homogenized.
- e. repeat steps b and c before moving to the next sample to clean the homogenizer

14 Centrifuge homogenate at  4 °C for  00:05:00 at  7000 rpm

5m

15 Transfer supernatant and place into a separate eppendorf tube labeled with sample name.

#### Estimation of tissue homogenate dilution

16 Prepare varying Tissue/Cell lysate : Lysis Buffer (TL:LB) dilutions (subject to change with varying sample type):

a. 2X (  100 µl TL :  100 µl LB)

b. 5X (  40 µl TL :  160 µl LB)

c. 10X (  10 µl TL :  190 µl LB)

Keep tissue lysate on ice until ready to pipette into 96 well plate.



17 Determine appropriate tissue/cell dilution to be used for the assay.

a. Use  0.6 Milimolar (mM) pyruvate concentration and plate triplicates of different sample dilutions.

b. In a 96-well plate, add  3 µl of each TL:LB dilution in triplicate.

c. Add  147 µl of the  0.6 Milimolar (mM) pyruvate concentration to each well.


18 Place the 96-well plate in the spectrophotometer and measure immediately.

19  147 µl Once the measurement has completed in  00:10:00 , Calculate a linear regression for each well to obtain the kinetic output slope (OD/min) and determine the correct dilution by identifying the dilution with the most linear slope.

10m

#### LDH activity measurement

20 

Pipette  3 µl of appropriate TL:LB dilution into three rows (technical replicates) and 12 columns (corresponding to the 12 pyruvate dilutions) on the 96 micro well plate.

Note: Ensure that there are no bubbles

21 

Using a  200 µl multichannel pipettor, carefully pipette  147 µl of each 12 pyruvate dilutions in triplicate; assay

immediately.  
Note: Ensure that there are no bubbles

22 Place the 96-well plate in the spectrophotometer and measure immediately.

23 Once the read has completed in  00:10:00 , click Export to excel and save the file. Eject the plate and clean up <sup>10m</sup>

#### Data Analysis

24 Do linear regression for each well (change in absorbance over the 10 minutes) to obtain the slope OD/min. (Note: the TECAN Infinity200M accounts automatically for the  $\pm 1$  cm path length. So the measurement is actually OD/min/cm)

25 LDH enzyme ctivity is calculated as  $\mu\text{M NADH/min} = \text{OD/min (slope)} * \text{Solution Volume in each well ( } \boxed{150 \mu\text{l}} \text{ )} / ( \text{volume of sample lysate ( } \boxed{3 \mu\text{l}} \text{ )} * e )$   
where  $e = \text{extinction coeff} = 0.00622 \text{ uM}^{-1}\text{cm}^{-1}$   
Note: take the absolute value of OD/min (omit the negative sign)

26 Take the average (of the three technical replicates) activity per pyruvate concentration subject to a threshold of  $R_{sq} > 0.9$

27 Plot Tissue Activity **Michaelis-Menten** curve  
a. X = Pyruvate Concentration (mM)  
b. Y = Average Activity

28 For simpler linear curve

28.1 plot Lineweaver-Burke:  
a. X =  $1/\text{Pyruvate Concentration}$   
b. Y =  $1/\text{Average Activity}$

28.2 Obtain trend-line linear fit in the Excel program and calculate  
a.  $K_M = \text{slope} / \text{intercept} ; [\text{mM}]$   
b.  $V_{\text{max}} = 1 / \text{intercept} ; \mu\text{M NADH/min}$

Use Bradford assay kit (<https://www.thermofisher.com/order/catalog/product/23200#/23200>) to determine protein concentration to express  $V_{\text{max}}$  (LDH activity) in  $\mu\text{mols NADH/min/mg}$  of protein