



Jan 18, 2022

Protocols from: Evolutionary analyses of visual opsin genes in frogs and toads: diversity, duplication, and positive selection

Ryan K Schott¹, Leah Perez², Matthew A Kwiatkowski², Vance Imhoff³, Jennifer Gumm³

¹York University; ²Stephen F. Austin State University; ³US Fish and Wildlife Service

1



dx.doi.org/10.17504/protocols.io.b3wgqpbw

jennifer_gumm

Protocols used to extra mRNA from frog retina, create cDNA libraries, and prepare sampled for sequence at the UT core facility under standard protocols.

DOI

dx.doi.org/10.17504/protocols.io.b3wgqpbw

Ryan K Schott, Leah Perez, Matthew A Kwiatkowski, Vance Imhoff, Jennifer Gumm 2022. Protocols from: Evolutionary analyses of visual opsin genes in frogs and toads: diversity, duplication, and positive selection. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b3wgqpbw>



protocol

Schott RK, L Perez, MA Kwiatkowski, V Imhoff, JM Gumm. 2022. Evolutionary analyses of visual opsin genes in frogs and toads: diversity, duplication, and positive selection. Ecology and Evolution (in press).

protocol ,









Jan 15, 2022








Jan 18, 2022

57000



Frog retinas.

mRNA Extraction

- 1 Transfer sample into a 1.5ml collection tube
- 2 Pipette off RNALATER
- 3 Add  **600 µL** Buffer RLT
- 4 Add  **6 µL** Beta-mercaptoethanol
- 5 Disrupt tissue with sterile pestle
- 6 Pipette into Qiashredder column- Spin  **00:02:00** @ 8,000rpm 2m
- 7 Remove Qiashredder column; Add cap - Spin  **00:03:00** @ max speed 3m
- 8 Add  **600 µL** 70% Ethanol to new collection tube
- 9 Transfer lysate to the collection tube- mix lysate and 70% Ethanol by pipetting
- 10 10. Transfer lysate to RNeasy column ( **700 µL** at a time)- Spin  **00:00:15** @ 9,800rpm;
Discard flow through/ Add rest of lysate; Spin  **00:00:15** @ 9,800rpm; Discard flow through 30s




- 11 11. Add  **700 µL** Buffer RWI- Spin  **00:00:15** @ 9,800rpm 15s
- 12 12. Transfer RNeasy column to new collection tube
- 13 13. Add  **500 µL** Buffer RPE- Spin  **00:00:15** @ 9,800rpm; Discard flow through 15s
- 14 14. Add  **500 µL** Buffer RPE – Spin  **00:01:00** @ 9,800rpm; Discard Flow through/ Spin ^{3m}
 **00:02:00** @ 13,000rpm
- 15 15. Transfer RNeasy column to new collection tube
- 16 16. Elute with  **30 µL** RNase-Free H2O- Spin  **00:01:00** @ 13,00rpm 1m
- 17 17. Elute with  **30 µL** RNase-Free H2O- Spin  **00:01:00** @ 13,00rpm 1m

cDNA Synthesis 1h 6m

- 18 1. Combine mRNA and RNase-free H2O to standardize all samples to aliquots containing 0.4µg mRNA total in 10µL.
- 19 2. Make 2 Master mixes
- 20 Master Mix 1: add  **1 µL** dntp mix and  **2 µL** dT primer per sample

- 21 Master Mix 2: add  **4 μ L** Buffer,  **4 μ L** DTT and  **0.5 μ L** RNAase inhibitor per sample
- 22 3. Pipette  **3 μ L** of Master Mix 1 into each sample.
- 23 4. Place sample on dry bath at  **65 °C** for  **00:05:00** . 5m
- 24 5. Put samples  **On ice** for  **00:01:00** . 1m
- 25 6. Pipette  **6.5 μ L** of Master Mix 2 into each sample.
- 26 7. Pipette  **1 μ L** Superscript into each sample.
- 27 8. Incubate samples  **Room temperature** for  **00:10:00** . 10m
- 28 9. Incubate samples at  **42 °C** for  **00:50:00** . 50m

PCR 17m 30s

- 29 1. Keep all reagents on ice at all times.
- 30 2. Make a master mix. Per sample add the following:
 **2.0 μ L** 10X Buffer
 **1.0 μ L** 50mM MgSO₄
 **0.5 μ L** dNTP mix (10mM each)

- ▢ **18.4 µL** ddH₂O
- ▢ **1 µL** forward primer (10µM)
- ▢ **1 µL** reverse primer (10µM)
- ▢ **0.5 µL** taq polymerase

31 3. Mix well by spinning.

32 4. Add ▢ **24 µL** of Master Mix to each PCR tube.

33 5. Add ▢ **1 µL** of sample for a total of ▢ **25 µL** per tube.

34 6. Program the thermocycler for the following program:

17m 30s

⌚ **95 °C** for ⌚ **00:10:00**

⌚ **94 °C** for ⌚ **00:02:00**

REPEAT FOLLOWING 3 steps 35-50 times

⌚ **94 °C** for ⌚ **00:00:30**

⌚ **45 °C** - ⌚ **50 °C** for ⌚ **00:01:00** *temp depends on primer

⌚ **72 °C** for ⌚ **00:02:00**

END Repeat

⌚ **72 °C** for ⌚ **00:02:00**

⌚ **4 °C** hold