



SEP 26, 2023

Polymerase Chain Reaction (PCR)

NUS iGEM¹¹National University of Singapore

NUS iGEM 2023



NUS iGEM

National University of Singapore

ABSTRACT

This is the common Polymerase Chain Reaction (PCR) protocol followed by the 2023 NUS-Singapore iGEM team.

MATERIALS

1. KOD One™ PCR Master Mix (Blue)
2. DNA Template
3. Primer(s)
4. DI Water

SAFETY WARNINGS



Proper lab PPE must be worn.

OPEN  ACCESS**DOI:**dx.doi.org/10.17504/protocols.io.261gedrndv47/v1

Protocol Citation: NUS iGEM 2023. Polymerase Chain Reaction (PCR).

protocols.io<https://dx.doi.org/10.17504/protocols.io.261gedrndv47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Sep 12, 2023**Last Modified:** Sep 26, 2023

Keywords: PCR, Polymerase
Chain Reaction, KOD One™
PCR master Mix, iGEM

- 1 Prepare a PCR tube.
- 2 Thaw the KOD One™ PCR Master Mix (Blue) from the -20°C fridge.

- 3 For a 50µL PCR sample, add the following into a PCR tube (follow the sequence):

Item	Volume
KOD One™ PCR Master Mix (Blue)	25µL
Each Primer (2 primers in total)	1.5µL
DNA Template	1µL
DI Water	21µL

- 4 Vortex the PCR tubes to mix well.

- 5 Place the PCR tube into the Thermal Cycler with the following conditions:

****Lid Temperature: 105°C; Volume: 50µL**

A	B	C
98°C	10s	
68°C	5s	Edit: -1°C per cycle
68°C	1 minute	
Go to step 1, repeat the cycle 13 times		

A	B	C
98°C	10s	
55°C	5s	
68°C	1 minute	
Go to step 5, repeat the cycle 25 times		
68°C	5 minutes	
12°C	Infinite Loop (∞)	

- 6** Upon finishing, cancel the current protocol and take out the PCR sample from the Thermal Cycler. Keep the sample in the 4°C fridge for storage. Otherwise, gel electrophoresis may be carried out to isolate the DNA fragment of interest.