

Jul 12, 2024

Immunostaining- Fluorescent

DOI

dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1

daniel.dautan daniel^{1,2}, Per Svenningsson^{1,2}

¹Department of Clinical Neuroscience, Karolinska Institutet, 171 76 Stockholm, Sweden;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Rese...

Kaplitt Protocols



Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1

Protocol Citation: daniel.dautan daniel, Per Svenningsson 2024. Immunostaining- Fluorescent. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk88y6l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 04, 2024

Last Modified: July 12, 2024

Protocol Integer ID: 101197

Keywords: ASAPCRN, immunofluorescence, mouse, brain,

Funders Acknowledgement:

Aligning Science Across







Parkinson's

Grant ID: 020608

Abstract

Staining of mouse brain sections with immunofluorescent visualization.



- 1 Wash freshly sectioned slices of tissue 2-3 times with 1X PBS to remove OCT.
- 2 Block with 5% Normal donkey/Goat serum in 0.3% TRITON X-100 in 1X PBS for  01:00:00 at  Room temperature . 1h
- 3 Wash sections 5 times in 1X PBS.
- 4 Following washes, transfer sections to a primary antibody solution containing a mix of primary antibody diluted in 0.3% TRITON X-100 in 1X PBS with 1% Normal Donkey/Goat serum  Overnight at  4 °C . 1h
- 5 The following day, wash sections 5 times in 1X PBS.
- 6 Transfer sections to secondary solution containing the adequate secondary antibody diluted in 0.3% TRITON X-100 in 1X PBS with 1% Normal Donkey/Goat serum for  04:00:00 at  Room temperature . 4h
- 7 Wash sections 3 times in 1X PBS.
- 8 Mount sections on microscope slides with medium (Vectashield with 4',6-diamidino-2-phenylindole (DAPI), H-1800-10).