

FEB 06, 2024

# 🌐 Cloning, Protein Expression, and Purification of 20S CPs and Assembly Intermediates

📁 In 1 collection

Frank Adolf<sup>1</sup>

<sup>1</sup>Department of Molecular Machines and Signaling, Max Planck Institute of Biochemistry, 82152 Martinsried, Germany

ASAP Collaborative Research Network

OPEN  ACCESS



Frank Adolf

Department of Molecular Machines and Signaling, Max Planck I...

## ABSTRACT

This protocol details methods for cloning, expression, and purification of 20S CPs and assembly intermediates for biochemical and structural analysis.

## GUIDELINES

Please familiarise yourself with the laboratory safety rules and guidelines and follow these while performing the experiment. Please wear appropriate PE while performing the experiment.

## MATERIALS

X-tremeGENE<sup>TM</sup> HP DNA Transfection Reagent - ROCHE  
cOmplete PROTEASE INHIBITOR COCKTAIL - ROCHE  
Strep-Tactin<sup>®</sup> Sepharose<sup>®</sup> resin - IBA

## SAFETY WARNINGS



Please refer to the Safety Data Sheets (SDS) for health and environmental hazards. Liquid nitrogen (LN2) and other cryogenics can cause severe damage to the skin and eyes. Always wear personal protective equipment when handling these cryogenics.

## DOI:

[dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1)

**Protocol Citation:** Frank Adolf 2024. Cloning, Protein Expression, and Purification of 20S CPs and Assembly Intermediates. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldmd59l5b/v1>

## MANUSCRIPT CITATION:

<https://www.biorxiv.org/content/10.1101/2024.01.27.577538v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 05, 2024

**Last Modified:** Feb 06, 2024

**PROTOCOL integer ID:** 94683

**Keywords:** ASAPCRN, proteasome, core particle, 20S proteasome, chaperone, molecular machine, multiprotein complex, POMP, PAC1, PAC2, PAC3, PAC4, propeptide, protease

**Funders Acknowledgement:**  
Aligning Science Across Parkinson's (ASAP)  
Grant ID: ASAP-000282

# Cloning of baculovirus transfer vectors

1w

1 Baculovirus transfer vectors were assembled utilizing a combination of the biGBac and MultiBac systems



## CITATION

Weissmann F, Petzold G, VanderLinden R, Huis In 't Veld PJ, Brown NG, Lampert F, Westermann S, Stark H, Schulman BA, Peters JM (2016). biGBac enables rapid gene assembly for the expression of large multisubunit protein complexes..

LINK  
<https://doi.org/10.1073/pnas.1604935113>

## CITATION

Vijayachandran LS, Viola C, Garzoni F, Trowitzsch S, Bieniossek C, Chaillet M, Schaffitzel C, Busso D, Romier C, Poterszman A, Richmond TJ, Berger I (2011). Robots, pipelines, polyproteins: enabling multiprotein expression in prokaryotic and eukaryotic cells..

LINK

<https://doi.org/10.1016/j.jsb.2011.03.007>

The MultiBac vector pACEBac1 was used as both library vector and acceptor vector for step 1 biGBac assembly reactions with the following primers:

ACEBac-1EC-BBA.rev

ATTTAAATCTTTAGACCATAGAGCGTTCTCGCGAATCGATACTAGTGTTTAAACTCGCTACCTTAGGACC

ACEBac-2EC-BBA.fwd

ATTTAAATAAACCTAATGATGCCTGATGTTTCCTAGGGTATACCCATCTAATTGGAACCAGATAAGTGAAATC

ACEBac-3EC-BBA.fwd

ATTTAAATAAACGGTTCACATAGCTTAGTTTCCTAGGGTATACCCATCTAATTGGAACCAGATAAGTGAAATC

ACEBac-4EC-BBA.fwd

ATTTAAATAAACACTGACATTGACTTGTTTCCTAGGGTATACCCATCTAATTGGAACCAGATAAGTGAAATC

ACEBac-5EC-BBA.fwd

ATTTAAATAAATCTATATCTCAATCGGGTTCCTAGGGTATACCCATCTAATTGGAACCAGATAAGTGAAATC

## 2 Clone all 20S CP subunits and assembly chaperones into pACEBac1 using Gibson assembly

1d



For affinity purification add C-terminal TEV-cleavable twin strep tags on  $\beta 2$  (PSMB7) and  $\beta 7$  (PSMB4)

### 2.1 Screen for positive clones and sequence verify by Sanger sequencing

1d



### 3 Preparation of bigBac assembly inserts:

6h



Amplify expression cassette from library vectors (step 2) by PCR with the following primers:

Cas1-ACEBac.fwd

AACGCTCTATGGTCTAAAGATTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas2-ACEBac.fwd

AAACTGGATACTATTGCACGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas3-ACEBac.fwd

AAACCTAATGATGCCTGATGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas4-ACEBac.fwd

AAACGGTTCACATAGCTTAGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas5-ACEBac.fwd

AAACACTGACATTGACTTGGTTTAAATCGACCTACTCCGGAATATTAATAGATCATGG

Cas1-ACEBac.rev

AAACGTGCAATAGTATCCAGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas2-ACEBac.rev

AAACATCAGGCATCATTAGGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas3-ACEBac.rev

AAACTAAGCTATGTGAACCGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas4-ACEBac.rev

AAACCAAGTCAATGTCAGTGTTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

Cas5-ACEBac.rev

AACCCCGATTGAGATATAGATTTATTTAAATGGTTATGATAGTTATTGCTCAGCGGTGG

#### 3.1 Purify all inserts by agarose gel extraction

2h



4



Preparation of biGBac step 1 assembly acceptor vectors:  
linearize the acceptor vector pACEBac1 by PCR with the following primers:

for cloning of 5 expression cassettes:  
ACEBac-1EC-BBA.rev and ACEBac-5EC-BBA.fwd

for cloning of 4 expression cassettes  
ACEBac-1EC-BBA.rev and ACEBac-4EC-BBA.fwd

for cloning of 3 expression cassettes  
ACEBac1-1EC-BBA.rev and ACEBac-3EC-BBA.fwd

**4.1** Purify all linearized vectors by agarose gel extraction



6h

5



Clone multi expression cassette Baculo transfer vectors listed below by biGBac step 1 assembly with expression cassette inserts from step 3 and acceptor vectors from step 4

pACEBac1-POMP-PSMG1-PSMG2-PSMG3-PSMG4  
pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4,  
pACEBac1-PSMA5-PSMA6-PSMA7,  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4,  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII,  
pACEBac1-PSMB5-PSMB6-PSMB7,  
pACEBac1-PSMB5-PSMB6-PSMB7-TEV-2xSTII.

**5.1** Screen for positive clones and sequence verify by Sanger sequencing



1d

1d

6



Assemble finale multi expression cassette baculo transfer vectors by multibac assembly

pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4-PSMA5-PSMA6-PSMA7  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII-PSMB5-PSMB6-PSMB7  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-PSMB5-PSMB6-PSMB7-TEV-2xSTII

2d

**6.1** Digest acceptor vectors listed below with I-CeuI, CIP and purify by agarose gel extraction 4h



pACEBac1-PSMA5-PSMA6-PSMA7  
pACEBac1-PSMB5-PSMB6-PSMB7  
pACEBac1-PSMB5-PSMB6-PSMB7-TEV-2xSTII

**6.2** Digest donor vector listed below with I-CeuI and BstXI, and purify inserts by agarose gel extraction 4h



pACEBac1-PSMA1-PSMA2-PSMA3-PSMA4  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII  
pACEBac1-PSMB1-PSMB2-PSMB3-PSMB4

**6.3** Ligate the following vector/insert pairs listed below and transform in DH5alpha *E. coli* 1d



pACEBac1-PSMA5-PSMA6-PSMA7 and EC1-4 PSMA1-PSMA2-PSMA3-PSMA4  
pACEBac1-PSMB5-PSMB6-PSMB7 and EC1-4 PSMB1-PSMB2-PSMB3-PSMB4-TEV-2xSTII  
pACEBac1-PSMB5-PSMB6-PSMB7 TEV-2xSTII and EC1-4 PSMB1-PSMB2-PSMB3-PSMB4

**6.4** Screen for positive clones and sequence verify by Sanger sequencing 1d



## Baculo virus amplification and insect cell expression 1w

**7** Culture Sf9 insect cells (Thermo Fisher Scientific) for virus amplification in serum-free Ex-cell 420 medium (Sigma-Aldrich)



Culture Trichoplusia ni (Thermo Fisher Scientific) in protein free ESF 921 insect cell culture media (Expression Systems LLC)

All steps were carried out according to standard protocols

## CITATION

Fitzgerald DJ, Berger P, Schaffitzel C, Yamada K, Richmond TJ, Berger I (2006). Protein complex expression by using multigene baculoviral vectors..

LINK

<https://doi.org/>

## CITATION

Bieniossek C, Richmond TJ, Berger I (2008). MultiBac: multigene baculovirus-based eukaryotic protein complex production..

LINK

<https://doi.org/10.1002/0471140864.ps0520s51>

**8** For bacmid preparation transform baculo transfer vector from the section above into EMBAC *E. coli* and plate on LB-agar plates with ampicillin, kanamycin, tetracyclin, gentamycin, IPTG, and XGal

1d



**8.1** Prepare Overnight cultures from white colonies in LB with ampicillin, kanamycin, and gentamycin

1d



**8.2** Prepare bacmids by alkaline lysis and subsequent isopropanol and 70 % (v/v) EtOH precipitation

2h



**8.3** Air dry bacmid pellets, resuspend DNA in milliQ water, and store at 4 °C until usage

1h



**9** P1 virus production

2d 12h



**9.1** For P1 virus production seed  $0.7 - 0.8 \times 10^6$  cells/well in **3 mL** medium in a six well plate and leave for **00:15:00** min at **27 °C** **15m**



**9.2** Prepare bacmids by diluting **1 µg** bacmid DNA in **20 µL** milliQ water, and add **200 µL** of medium **5m**



**9.3** For each bacmid mix **100 µL** medium with **15 µL** Extreme Gene HP DNA transfection reagent (Roche) **5m**



**9.4** Add **112 µL** of Medium Extreme Gene Mix to each bacmid DNA **5m**



**9.5** Add **156 µL** of the DNA-Extreme Gene Mix to each well in the six-well plate **5m**



**9.6** Wrap the plates with Parafilm and incubate for **60:00:00** h at **27 °C** **2d 12h**



**9.7** Transfer medium containing P1 virus from each well into a **15 mL** tube and store at **4 °C** until further usage **10m**





## 10 P2 and P3 virus amplification

2d



**10.1** Seed  $0.8 - 1.0 \times 10^6$  Sf9 cells/ml into an appropriate conical flask and infect cells with 0.5 - 1.0% P1 or P2, respectively 1h



**10.2** Incubate cell suspensions at  $27^\circ\text{C}$  at 80 rpm for 48:00:00 h 2d



**10.3** Harvest P2 and P3 by centrifugation, transfer medium containing virus into a 50 mL tube 15m  
250 mL bottle and store at  $4^\circ\text{C}$  until further usage



## 11 Protein Expression in High Five™ insect cells

2d 12h





**11.1** Seed  $2.0 \times 10^6$  High Five™ cells/ml into an appropriate conical flask and infect with 0.5 - 1.0% P3 virus 30m










**11.2** Incubate cell suspensions at  $27^\circ\text{C}$  at 80 rpm for 60:00:00 h 2d 12h



- 11.3** Harvest High Five™ cells by centrifugation, resuspend cell pellets in PBS, transfer to a **50 mL** tube, harvest cells by centrifugation, discard supernatant, snap freeze cell pellets in LN2, and store until further usage at **-80 °C** **30m**
-  

## Purification of 20S CPs and 20S CP assembly intermediates **8h**

- 12** Purification of mature 20S CPs together with their assembly intermediates **8h**

- 
- 12.1** Thaw cell pellets from section 2 - step 11, resuspend in buffer A 25 mM HEPES pH 7.5 (KOH), **30m**
-   **[M] 150 millimolar (mM) NaCl, [M] 1 millimolar (mM) DTT, and add 1 tablet cOmplete protease inhibitor mix (Merck) per 10 mL cell pellet resuspended in 40 mL buffer A in a 50 ml tube**
- 12.2** Lyse cells by sonication on ice **10m**
- 
- 12.3** Centrifuge lysate at 22k at **4 °C** for **01:00:00** h **1h**
- 
- 12.4** In parallel wash Strep-Tactin® Sepharose® resin (IBA) 3 times with buffer A **45m**
- 
- 12.5** Incubate supernatants from step 12.3 with resin for **01:00:00** h at **4 °C** **1h**
- 



12.6

Wash resin 3 times with buffer A by centrifugation at 2000 x g for 00:15:00 min at

15m



4 °C

12.7

Load resin on a gravity flow column and elute proteins with 2.5 micromolar ( $\mu\text{M}$ ) d-  
Desthiobiotin in buffer A

15m

12.8

Pool fractions of interest, concentrate in a spin protein concentrator with a MWCO of 100 kDa to  
max 5-8 mg/mL

1h

12.9

In parallel equilibrate a Superose 6 10/300 GL column (Cytiva) with buffer A

1h

12.10

Apply concentrated eluate from step 12.8 onto the SEC column at a flow rate of 1 ml/min, record  
absorption at 280 nm, and fractionate elution a 300  $\mu\text{L}$

1h

12.11

Analyze SEC fractions by SDS-PAGE, pool fractions of interest, snap freeze in LN2, and store until  
further usage at -80 °C for biochemical assays

2h



For structural analysis use fractions directly after SEC