



DEC 10, 2022

## HotSHOT genomic DNA extraction

This protocol is published without a DOI.

FishFloorUCL<sup>1</sup>

<sup>1</sup>University College London



Francois Kroll

COMMENTS 0

WORKS FOR ME

1

### ABSTRACT

How to extract genomic DNA from larvae or finclips using the HotSHOT method.

### PROTOCOL CITATION

FishFloorUCL 2022. HotSHOT genomic DNA extraction. **protocols.io**  
<https://protocols.io/view/hotshot-genomic-dna-extraction-ckbhusj6>

### LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Dec 09, 2022

### LAST MODIFIED

Dec 10, 2022

### PROTOCOL INTEGER ID

73801

## Materials

### 1 Prepare the BASE solution (50X):

- 14.03 g KOH crystals (1.25M final concentration)
- 4 mL of 0.5M EDTA (10 mM final concentration)
- ddH2O to 200 mL total volume

### 2 Prepare the NEUTRALISATION solution (50X)

- 63.04 g Tris-HCL (2M final concentration); also called Trizma HCl
- ddH2O to 200 mL total volume

## Procedure

- 3 Prepare fresh 1X BASE and NEUTRALISATION solution in nuclease-free H<sub>2</sub>O.
- 4 If larvae: add larvae one per well with Pasteur pipette. Take up all the liquid using a P200 pipette.  
If finclip: add finclip directly to well (this can be done while finclipping).
- 5 Add 50 µL of 1x BASE Solution into each well of the plate.
- 6 Seal the plate and place on PCR block:

🔥 95 °C

🕒 00:30:00

30m
- 7 Cool at

🔥 Room temperature
- 8 Add 50 µL of 1x NEUTRALISATION solution into each well of the plate.
- 9 Note, the extracted DNA concentration is often *too high* for downstream applications like PCR or KASP.

## Storage

- 10 Store at

🔥 4 °C

 if you will use the DNA in the next few weeks. Beware, the samples will slowly evaporate from sealed plate  
or  

🔥 -20 °C

 for long-term storage. This will also prevent the samples from evaporating.