




APR 15, 2024

Histology Protocol

 In 2 collections

Sasha Burwell¹

¹Duke University

ASAP Collaborative Research Network



Sasha Burwell
Duke University

OPEN  ACCESS



ABSTRACT

This protocol details the histological processing of a mouse brain.

DOI:

dx.doi.org/10.17504/protocols.io.q26g71bdqgwz/v1

Protocol Citation: Sasha Burwell 2024. Histology Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g71bdqgwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Apr 02, 2024

Last Modified: Apr 15, 2024

Histology Pro31hocol

1d 7h 5m 1s

1 Prepare 4% paraformaldehyde in [M] 0.1 Molarity (M) PB, and pH 7.4 . Let chill to 4 °C .



2 Deeply anesthetize mouse with isofluorane until completely under (no toe pinch response, low to no breathing).

3 If the mouse has an electrode bundle implanted, quickly attach their tongue with an alligator clip to the negative port of a 9V battery, and touch each electrode connection for 00:00:01 with a wire attached to the positive port of the battery. This will create a small lesion for histological viewing of electrode placement.

1s






4 Fix with a transcardial perfusion of 15 mL PBS followed by 50 mL of 4% PFA.


5 Carefully remove the brain from the skull and post-fix in 50 mL of 4% PFA at 4 °C Overnight


8h





6 The next day, wash 3x with PBS, with a 15 min agitation in a cold room for the first wash, and 15 mins or longer agitation in a cold room for the remaining washes.


- 7 The next day, wash with PBS, with a  00:15:00 agitation in a cold room for the first wash (1/3). 15m
- 8 The next day, wash with PBS,  00:15:00 or longer agitation in a cold room for the remaining washes (2/3). 15m
- 9 The next day, wash with PBS,  00:15:00 or longer agitation in a cold room for the remaining washes (3/3). 15m
- 10 Embed the brains in 4% agarose.
- 11 Using a Leica VT1200S vibratome, slice the brain along the coronal axis at  50 undetermined thick slices.
- 12 If performing tyrosine hydroxylase immunostaining:
 - 12.1 Place sections in a 24 well plate, 2-3 sections per well.
 - 12.2 Wash 1x in PBS.


12.3 Incubate for  02:00:00 in a blocking solution consisting of 5% goat serum, 3% bovine serum albumin, and 0.3% triton-x. 2h





12.4 Then transfer sections to a half-block solution containing 1:1000 rabbit α -TH, and leave 8h


 Overnight at  4 °C on a shaker.





12.5 Wash in  0.1 Molarity (M) PBS with 0.1% tween.




12.6 Incubate for  04:00:00 in a half-block solution containing 1:1000 goat α -rabbit 488 . 4h




12.7 Wash in  0.1 Molarity (M) PBS with 0.1% tween.





12.8 Wash in PBS.



13 Mount the sections onto SuperFrost glass slides. Let dry completely, then refresh with PBS for  00:20:00 . 20m

14 Rinse slides with MilliQ, then coverslip with  155 μ L of Vectashield Vibrance + DAPI mounting medium.



15 Let set at  Room temperature for at least 2 hours, and up to  Overnight . 8h

16 When set, image slides on a VS200 slide scanner using the following protocol:

16.1 Overview in brightfield.

16.2 Focus in 10% Tritc (make sure focus points are not over the VTA).

16.3 Image at:

1. 300ms 50% Cy5
2. 193ms 10% Tritc
3. 100ms 10% Fitc (only if GCaMP or TH present; otherwise, turn off this channel)
4. 20ms 10% DAPI