

FEB 29, 2024

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.3b4l4qpxrvo5/v1

Protocol Citation: Carlos Goller 2024. 16S Gene PCR Amplification and Sanger Sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3b4l4qpxrvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Feb 28, 2024

16S Gene PCR Amplification and Sanger Sequencing

Carlos Goller¹

¹North Carolina State University

BIT Metagenomics BIT-Protocols [1 more workspace](#) ↓



Carlos Goller

ABSTRACT

After genomic DNA (gDNA) has been isolated and quantified, it can be used for many downstream applications. While we will use our genomic DNA samples for whole genome sequencing using Nanopore sequencing technology, we can also use the gDNA to amplify one gene for Sanger sequencing. These different sequencing methods contrast one another:

Table. Comparison of Nanopore and Sanger DNA Sequencing.

Nanopore Sequencing	Sanger Sequencing
High Throughput (many genes at once)	Low Throughput (one gene at a time)
Sequences gDNA directly	Sequences PCR amplified DNA
Ion current base calling	Fluorophore base calling

PCR amplification of 16S

We will amplify the bacterial 16S gene, which encodes a subunit of ribosomes. Biologists use this [gene to compare evolutionary relationships](#) by sequence similarities and differences to organize biological taxonomies. Polymerase Chain Reaction (PCR) enables scientists to amplify many copies of a gene, and we will target 16S for amplification.

IMAGE ATTRIBUTION

Image created with BioRender.com by Carlos C. Goller

Last Modified: Feb 29, 2024

PROTOCOL MATERIALS

PROTOCOL integer ID: 95917

QIAquick PCR Purification Kit **Qiagen Catalog #28104** Step 3

Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0494S**

Keywords: metagenomics, Nanopore sequencing, WGS, 16S sequencing

In 2 steps


Funders Acknowledgement:

Biotechnology Program
Grant ID: NCSU

Introduction

1

Polymerase chain reaction - PCR



Polymerase chain reaction - PCR. Created with BioRender.

PCR Cycles Diagram indicating the three main steps: denaturation at 95-96°C, annealing at 68°C, and elongation at 72°C. Created with BioRender.com

You can learn about this reaction by watching the video "[Polymerase Chain Reaction.](#)"

After completing this lab, you will gain the following lab skills:

- Lab safety and proper personal protective equipment (PPE)
- Setup of two PCR reactions.
- Proper use of a thermocycler for PCR.

PCR Setup

2 Note Before You Begin:

Review the protocols and figures below to learn how the 16S PCR will help us amplify regions of this gene.

A	B	C
Reagent	Volume	Final Concentration
DNA polymerase master mix (contains proof-reading DNA polymerase, dNTPs, buffer, Mg ²⁺) NEB Q5 Hot Start High-Fidelity 2X Master Mix.	25 µL	1X
10 µM Forward Primer	2.50 µL	0.5 µM
10 µM Reverse Primer	2.50 µL	0.5 µM
Nuclease-free water	17.5 µL	-
gDNA (<100 ng/µL)	2.5 µL	<1,000 ng
FINAL REACTION	50 µL	-

PCR Table. Reagents and volumes needed for 16S PCR.

For each gDNA sample, set up **two reactions**:

- one LONG reaction with primers **27fwd + 1492rev**
- one SHORT reaction with **515fwd + 1492rev**

We will use the

Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0494S** and primers ordered from IDTA and prepared to be LabReady.

Q5 Hot Start High-Fidelity 2X Master Mix - 100 rxns **New England Biolabs Catalog #M0494S**


Primer Name	Primer sequence 5' → 3'
27fwd	AGA GTT TGA TCM TGG CTC AG
515fwd	GTG CCA GCM GCC GCG GTA A
1492rev	CGG TTA CCT TGT TAC GAC TT

Primer Table. Primer names and sequences for the amplification of 16S.

Source: [Identification of unknown bacterial isolates using Sanger sequencing of the 16S rRNA gene | CHMI services](#)

16S Amplification and Sanger Sequencing

Bacteria rRNA gene organization



Bacterial 16S Gene Schematic for Amplification and Sanger Sequencing. The structure of the 16S rRNA gene is presented along with the primers and their annealing sites. Created with BioRender.com

Temperature (°C)	Time	Repetitions
98°C	30 seconds	1X
98°C	10 seconds	35 X
63°C	30 seconds	
72°C	45 seconds	
72°C	2 minutes	1X
4°C	Hold	Hold



Thermocycler protocol table.

Thermocycler protocol table. Adapted from NEB Protocol for Q5 2X Hot Start Master Mix.



PCR Clean-up

- To remove the reagents from the PCR and isolate pure DNA for Sanger sequencing, we need to "clean up" our amplified DNA samples. In BIT 295, we use reagents from the

 QIAquick PCR Purification Kit **Qiagen Catalog #28104** Kit.

- Add all  50 µL of the PCR reaction to five volumes ( 250 µL) PB buffer

- Mix and transfer to pink spin column

-  Sample Centrifuge at  10.000 rpm, 00:00:30

30s


7 Remove flow through

8 Add 750 μ L PE to wash the column. Incubate at RT for 1 minute.

9 Centrifuge at  10.000 rpm, 00:00:30


30s



10 Remove the spin column from the collection tube and add to a new wash tube

11 Centrifuge at  10.000 rpm, 00:00:30 to remove all ethanol wash

30s

12 Remove the spin column from the collection tube and add it to a new collection tube

13 Add  50 μ L EB to elute. Incubate at RT for 1 minute.

14  2 μ L Centrifuge at  10.000 rpm, 00:00:30

30s

- 15 Remove the spin column and throw away
- 16 Quantify 2 µl using the Implen that has been blanked with the elution buffer (EB).
- 17 DNA can be stored at 4 °C or on ice while being used and should be stored long-term at -20 °C

Sanger Sequencing Sample Preparations

- 18 For Sanger sequencing, start by sending out 20 µL of ~100 ng/µL DNA of your LONG amplicon for three sequencing reactions with all three primers 27fwd, 515 fwd, and 1492rev. Prepare 10 µL a 1:10 dilution of your 100 µM primers.
- 19 If you do not get clear chromatograms to interpret the sequence, next send out your SHORT amplicon for Sanger sequencing with only primers 515fwd and 1992rev. Prepare 10µL a 1:10 dilution of your 100 µM primers.

[16S rRNA Sequencing: A PCR-based Technique to Identify Bacterial Species](#)

Adapted from: [Identification of unknown bacterial isolates using Sanger sequencing of the 16S rRNA gene | CHMI services.](#)