



MAR 30, 2023

WC Medium

Richard W Lambrecht¹

¹University of Konstanz

Richard W Lambrecht: This protocol is a digital version of the one used for a long time in the laboratory. It was not developed by the author.



Richard W Lambrecht

DISCLAIMER

This is a protocol used in the Becks Lab, 2023.

ABSTRACT

Protocoll description to prepare WC medium.

OPEN ACCESS

Protocol Citation: Richard W Lambrecht 2023. WC Medium. **protocols.io** <https://protocols.io/view/wc-medium-cr5dv826>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 30, 2023

Last Modified: Mar 30, 2023

PROTOCOL integer ID:
79749

Keywords: WC medium, Algal culture medium, Algal cultures

Stock Solution preparation

- 1 The stock solutions don't need to be autoclaved after preparation. Make sure to use volumetric flasks to ensure the right concentrations of the salts and transfer the solution to autoclavable flasks and store it in the fridge. Remember to start the preparation with less than 1L of MilliQ water, and after adding the salts complete the volume to 1L.

- 1.1** Stock solution **A: 36,80g of (CaCl₂+2H₂O)** in 1L of MilliQ water
- 1.2** Stock solution **B: 37g of (MgSO₄+7H₂O)** in 1L of MilliQ water
- 1.3** Stock solution **C: 12,60g of (NaHCO₃)** in 1L of MilliQ water
- 1.4** Stock solution **D: 11,40g of (K₂HPO₄ +H₂O)** in 1L of MilliQ water
- 1.5** Stock solution **E: 85 g of (NaNO₃)** in 1L of MilliQ water
- 1.6** Stock solution **F: 21,20 g of (Na₂SiO₃ + 5H₂O)** in 1L of MilliQ water
- 1.7** Stock solution **G (Supplement elements)**: add the following salts in 1L of MilliQ water:
- **4,36 g of (Na₂EDTA)**
 - **3,15 g of (FeCl₃ + 6H₂O)**
 - **0,010 g of (CuSO₄ + 5H₂O)**
 - **0,022 g of (ZnSO₄ + 7H₂O)**
 - **0,010 g of (CoCl₂ + 6H₂O)**
 - **0,180 g of (MnCl₂ + 4H₂O)**
 - **0,006 g of (Na₂MoO₄ + 2H₂O)**
 - **1,0 g of (H₃BO₃)**
- 1.8** Stock solution **H (Vitamin Mix)**: add the following in 1L of MilliQ water:
- **Thiamine hydrochloride (Vitamin B1)**: 0,1 g of (C₁₂H₁₇ClN₄O₄S)
 - **Biotin**: 0,0005 g of (C₁₀H₁₆N₂O₃S)
 - **Cyanocobalamin (Vitamin B12)**: 0,0005 g of (C₆₃H₈₈CoN₁₄O₁₄P)

Aliquot the solution in 15 mL falcon tubes, and freeze them until used at -20 °C.

Medium preparation (for 1L of WC)

- 2 Weigh **0,115 g of TES Buffer** ($C_6H_{15}NO_6S$) and add it to a volumetric flask with almost 1L of MilliQ water
- 3 Add **1 mL of each of the stock solutions A to G** to the medium preparation bottle.
- 4 **Complete** the volume until 1 liter with **MilliQ water**
- 5 Transfer the prepared medium to an autoclavable flask, labeled with autoclaving tape.
On the label put the information about the liquid contained (WC medium), the date, the person responsible for the preparation, and a box to check when adding the **solution H (Vitamin Mix)**.
- 6 **Autoclave** the medium at 121°C and let it cool down to room temperature.
- 7 In a clean bench: syringe filter 1 mL of the defrosted **solution H** (Vitamin Mix), add it to the medium (**1 ml per Liter**), and check the box on the tape.