



MAR 29, 2024

## DAB Staining of Fixed Mouse Brain Tissue Sections

Katerina Rademacher<sup>1</sup>, Ken Nakamura<sup>1</sup>

<sup>1</sup>Gladstone Institute of Neurological Disease



Haru Yamamoto  
UCSF

### ABSTRACT

This protocol describes steps for immunohistochemical 1 3,3'-Diaminobenzidine (DAB) staining of free floating fixed mouse brain tissue sections.

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1)

**Protocol Citation:** Katerina Rademacher, Ken Nakamura 2024. DAB Staining of Fixed Mouse Brain Tissue Sections.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldm127l5b/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** Jan 07, 2024

Last Modified: Mar 29, 2024

PROTOCOL integer ID: 93035

Keywords: ASAPCRN

Funders Acknowledgement:

ASAP

Grant ID: 020529

## MATERIALS

- 24-well netwell insert (Fisher NC9979302)
- dPBS (UCSF Stem Cell Core)
- Tris buffered saline (TBS)
- Methanol
- 30% H<sub>2</sub>O<sub>2</sub> (Millipore Sigma H1009)
- Endogenous peroxidase quenching buffer:
  - 10% MeOH, 3% H<sub>2</sub>O<sub>2</sub>, in 1X TBS
- TBS++++ blocking buffer:
  - 10mL 10X TBS
  - 10mL fetal bovine serum (UCSF Stem Cell Core)
  - 3g BSA (Millipore Sigma A3803)
  - 1g glycine
  - 30mg sodium azide (Millipore Sigma S2002)
  - Bring to 96mL with ddH<sub>2</sub>O
  - Store @-20°C in 12mL aliquots
  - Add 500 µL of 10% Triton x-100 (Millipore Sigma 93443) right before use
- Primary antibody (e.g. rabbit anti-TH, Abcam AB152)
- ABC complex - VECTASTAIN® Elite ABC-HRP Kit, Peroxidase (Rabbit IgG)
  - Manufacturer: Vector Laboratories
  - Cat#: PK-6101
  - Includes biotinylated goat anti-rabbit IgG secondary antibody (Vector Laboratories BA-1000-1.5)
- 50X 1 3,3'-diaminobenzidine (DAB) dissolved in 0.1M Tris
  - Millipore Sigma D12384
- 0.1 M Tris, pH 8.0
- Ethanol
- Xylene (ThermoFisher Scientific X51)
- Permout mounting medium (ThermoFisher Scientific SP15)
- Superfrost Plus slides (Fisher 12 550 15)

## Day 1

### 1 Blocking and primary antibody incubation.

#### 1.1 Pull all sections into 1x dPBS in 24-well netwells.

- 1.2 Wash 3 times with 1x dPBS, 10 mins each.
- 1.3 Wash 3 times with 1x TBS, 10 mins each.
- 1.4 Quench endogenous peroxidase activity with quenching buffer, 5 mins at RT. Just covering the sections is sufficient, and drain excess buffer by dabbing on kimwipes before and after.
- 1.5 Wash 3 times with 1x TBS, 10 mins each.
- 1.6 Block in TBS++++ for 75 mins at RT on rocker, 2 mL/well.
- 1.7 Incubate in primary antibody diluted in TBS++++ on rocker at RT O/N, 500 µL/well in a new 24-well plate without netwell insert.

## Day 2

- 2 Secondary antibody incubation and DAB development

- 2.1 Wash 3 times with 1x TBS, 5 mins each in netwell insert.
- 2.2 Incubate in secondary antibody at a 1:300 dilution in TBS++++ for 2 hrs, 500 µL/well in a 24-well plate without netwell insert.
- 2.3 15 mins prior to end of incubation, prepare ABC complex:
  - Both buffer A and B are 1:300 dilution in TBS++++
  - Incubate in RT beaker of water for 30 mins, in drawer
- 2.4 Wash 3 times with 1x TBS, 5 mins each in netwell insert.
- 2.5 Incubate 2 hours in ABC complex in TBS++++ at RT on rocker without netwell insert.
- 2.6 Wash 3 times with 0.1 M Tris buffer, pH 8.0, 5 mins each in netwell insert.
- 2.7 Prepare 1x DAB solution:
  - Thaw aliquots of 50x DAB, 1 tube is around 200 µL
  - Add 50x DAB into 0.1M Tris buffer
  - Vortex for 5 sec to mix well, sit in dark for 10 mins
  - Right before development, add 30% H<sub>2</sub>O<sub>2</sub> (1:10000 dilution), vortex to mix well
  - Distribute 2 mL to each well on a 12-well plate

- 2.8 Prepare an additional 12-well plate with Tris buffer to start washes immediately after development.
- 2.9 Develop DAB reaction, minimizing time difference between wells. Drain excess buffer by dabbing on kimwipes before and after.
- 2.10 Wash 2 times with Tris, 5 mins each.
- 2.11 Wash 2 times with TBS, 5 mins each.
- 2.12 Wash 2 times with PBS, 10 mins each.
- 2.13 Sections can be stored in PBS at 4C.
- 2.14 Mount onto Superfrost Plus slides labeled with pencil and allow to dry completely.

## Day 3

- 3 Dehydration and clearing.

- 3.1** Wash 2 x 2min in ddH<sub>2</sub>O.
- 3.2** Wash 2 x 2min in 70% EtOH.
- 3.3** Wash 2 x 2min in 95% EtOH.
- 3.4** Wash 2 x 2min in 100% EtOH.
- 3.5** Wash 2 x 2min in Xylene.
- 3.6** Coverslip using PermOUNT mounting medium.
- 3.7** Let dry in fume hood overnight then store in drawer in slide box until imaging.