



Jun 23, 2022

Silver Stain Gel for IP-TDMS Development

Lauren Adams¹¹Northwestern University1 *Works for me* Share

This protocol is published without a DOI.

Kelleher Research GroupTech. support email: kelleher-ofc@northwestern.edu**Kelleher KRG Research Group**

Northwestern University, National Resource for Translational...

ABSTRACT

Silver stain gels can be used to examine total protein in a given sample. This method is less sensitive than western blotting, but typically more sensitive than Coomassie blue staining.

PROTOCOL CITATION

Lauren Adams 2022. Silver Stain Gel for IP-TDMS Development. **protocols.io**
<https://protocols.io/view/silver-stain-gel-for-ip-tdms-development-cbwtspen>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 23, 2022

LAST MODIFIED

Jun 23, 2022

PROTOCOL INTEGER ID

65203

MATERIALS TEXT

Sample loading buffer with DTT

Microcentrifuge tube

Hot plate

Centrifuge

Gel box with power source

1 x MES running buffer

Invitrogen 4-12% BisTris gel (10, 12, or 17-well)

Protein MW Ladder

Clean gel box

MilliQ H₂O

Silver Nitrate Fixing Buffer (50% MeOH 45% H₂O, 5% acetic acid)

Silver Nitrate Washing Buffer (50% MeOH, 50% H₂O)

Silver Nitrate Sensitizing Buffer (0.2 g sodium thiosulfate/1 L)

Silver Nitrate Stain (0.05g AgNO₃ /50mL H₂O)

Silver Nitrate Developing Buffer (40 g Na₂CO₃, 800 µL formaldehyde/2 L)

Silver Nitrate Terminating Buffer (95% H₂O, 5% acetic acid)

- 1 Take aliquots from IP fractions and add sample loading buffer to a final concentration of 1X. Boil samples in hot plate at 95°C for 5-10 min. Briefly, centrifuge samples to ensure sample entire sample is towards bottom of the tube.
- 2 Open gel and place into gel box. Add 1 x MES running buffer, ensuring that there are no bubbles within the wells of the gel. Carefully load samples, without loss of sample from the well. Make sure that the gel wells have a large enough capacity to contain the entire sample volume. Load at least one lane with 2-3 µL of protein MW ladder. Run the gel at 120-150 V for approximately 40 minutes, or until the sample dye has reached the bottom of the gel.
- 3 Remove the gel from the cassette and gently wash with deionized water. Cut off the bottom of the gel (where there is a ridge). Place gel into clean gel box. Wash gel in MilliQ (MQ) H₂O 3 times for 5 minutes each.
- 4 Cover with Fixing Buffer for ≥ 20 minutes.
- 5 Wash in wash buffer for 10 minutes.
- 6 Wash in MQH₂O for 10 min.

- 7 Treat with sensitizing buffer for 1 minute.
- 8 Wash twice in MQH₂O for 1 minute each.
- 9 Stain AgNO₃ for 20 minutes; after this step, all waste must be collected in hazardous waste container.
- 10 Wash twice in MQH₂O for 20 seconds each.
- 11 Expose for 1 minute in developing buffer; watch the gel avoid overdeveloping.
- 12 Treat with terminating buffer for 1 minute; make sure the gel is no longer being developed.
- 13 Image the gel.