

May 19, 2020

SARS-CoV-2 Tailed Amplicon Illumina Sequencing

Daryl Gohl¹¹University of Minnesota

1 Works for me dx.doi.org/10.17504/protocols.io.bge5jtg6

Daryl Gohl

ABSTRACT

This protocol outlines how to process RNA for SARS-CoV-2 sequencing using tailed primers to generate tiled amplicons using the method described here: <https://www.biorxiv.org/content/10.1101/2020.05.11.088724v1>.

Best results are obtained for samples with N1 and N2 Ct values of <30 (based on the UMGC/MDL implementation of the CDC qRT-PCR diagnostic assay for SARS-CoV-2, see here: <https://www.biorxiv.org/content/10.1101/2020.04.02.022186v1.full>). For samples with N1 and N2 values between 30 and ~35, coverage and other sequencing metrics may be more variable and increased adapter dimer formation is expected.

MATERIALS TEXT

- 1) Fully skirted 96-well plate. (BioRad)
- 2) Semi-skirted 96-well plate (Thermo Scientific)
- 3) Nuclease-free water. (Fisher Scientific)
- 4) Microseal F foil seals. (BioRad)
- 5) Microseal B PCR seals. (BioRad)
- 6) SuperScript IV VIL0 master mix (Thermo)
- 7) Q5 Hot Start High Fidelity DNA polymerase. (NEB)
- 8) 10 mM dNTPs (NEB)
- 9) nCov-2019 pool 1.1, 1.2, 2.1, 2.2 primers. (IDT) – see Appendix
- 10) Indexing primers. (IDT) – see Appendix
- 11) Rainin Liquidator 96 pipette with p20/p200 tips. (Rainin)
- 12) Rainin single/multichannel pipettes with p20/p200/p1000 tips. (Rainin)
- 13) White Matrix troughs. (Thermo Scientific)
- 14) SequalPrep Normalization Plate Kit, 96-well. (Thermo Scientific)
- 15) AMPure XP beads. (Beckman Coulter)
- 16) Combinatorial Dual Indexing Primers:
For 384 sample barcoding scheme, see “Indexingprimers.xlsx”, from:
<https://protocolexchange.researchsquare.com/article/nprot-4831/v1>
- 17) Unique Dual Indexing Primers:
Available from Illumina (Nextera Unique Dual Indexing Primers, catalog number: 20027213, 20027214, 20027215, 20027216.

BEFORE STARTING

Tailed primers should be pooled to generate 4 primer pools (1.1, 1.2, 2.1, 2.2) according to the pooling scheme described in **Supplemental Data File 2** here: <https://www.biorxiv.org/content/10.1101/2020.05.11.088724v1.supplementary-material>.

Set up













10m

- 1 **Clean workspace and pipets by spraying with RNaseZAP or comparable product (such as RNase Away) and wiping down with KimWipes prior to beginning work.**











RNA samples should be stored at **-80 °C** and thawed on ice.

cDNA synthesis

- 2 Thaw RNA samples on ice then transfer **5 µl** of sample into a 96-well Thermo PCR plate.

- 3 Set up the following reverse transcription reaction master mix (multiply below volumes by number of reactions plus desired overage):
 -  **11 µl nuclease free water**
 -  **4 µl SuperScript IV VIL0 master mix**
- 4 Transfer  **15 µl** of reverse transcription master mix to each sample containing well.
- 5 Seal plate with a “B” seal, mix well by vortexing using a plate vortexer at 1900 rpm for  **00:00:10** s, and spin down briefly in a plate centrifuge ( **00:00:05** s at  **2500 rpm**).
- 6 Incubate in a thermocycler using the following conditions:
 -  **25 °C** for  **00:10:00**
 -  **50 °C** for  **00:10:00**
 -  **85 °C** for  **00:05:00**

Enrichment PCR

- 7 Transfer  **2.5 µl** of cDNA to each of 4 96-well Thermo PCR plates labeled: Project_Name_PCR1_1.1, Project_Name_PCR1_1.2, Project_Name_PCR1_2.1, and Project_Name_PCR1_2.2.
- 8 Set up the following four PCR master mixes, one for each of the four multiplexed primer pools (multiply below volumes by number of reactions plus desired overage):
 -  **14.75 µl nuclease-free water**
 -  **5 µl 5x Q5 reaction buffer**
 -  **0.5 µl 10mM dNTPs**
 -  **0.25 µl Q5 Polymerase**
 -  **2 µl primer pool (10 µM)** (Either pool 1.1, 1.2, 2.1, or 2.2)
- 9 Transfer  **22.5 µl** of master mix to each well of the appropriate PCR plate.
- 10 Seal plate with a “B” seal, mix well by vortexing using a plate vortexer at 1900 rpm for  **00:00:10** , and spin down briefly in a plate centrifuge ( **00:00:05** at  **2500 rpm**).
- 11 Amplify samples using the following PCR conditions:











⏸ **98 °C** for ⌚ **00:00:30**

35 cycles of:

⏸ **98 °C** for ⌚ **00:00:15**

⏸ **65 °C** for ⌚ **00:05:00**

Indexing PCR

- 12 For each sample, combine  **10 µl** of each of the four pools in a single Bio-Rad fully-skirted 96 well plate.
- 13 Seal plate with a “F” seal, mix well by vortexing using a plate vortexer at 1900 rpm for ⌚ **00:00:10** , and spin down in a plate centrifuge (⌚ **00:00:30** at ⚙ **2500 rpm**).
- 14 In a 96-well Thermo plate, add  **2 µl** of each sample to  **198 µl** of nuclease free water (1:100 dilution).
- 15 Seal plate with a “F” seal, mix well by vortexing using a plate vortexer at 2500 rpm for ⌚ **00:00:10** , and spin down in a plate centrifuge (⌚ **00:00:30** at ⚙ **2500 rpm**).
- 16 Transfer  **5 µl** of 1:100 diluted PCR 1 sample to a 96-well Thermo PCR plate.
- 17 Transfer  **2 µl** of 5 µM indexing primer mix to the 96-well Thermo PCR plate containing the samples.
- 18 Set up the following PCR master mix (multiply below volumes by number of reactions plus desired overage):
 -  **0.7 µl nuclease-free water**
 -  **2 µl 5x Q5 reaction buffer**
 -  **0.2 µl 10 mM dNTPs**
 -  **0.1 µl Q5 Polymerase**
- 19 Transfer  **3 µl** of master mix to each well of the appropriate PCR plate.
- 20 Seal plate with a “B” seal, mix well by vortexing using a plate vortexer at 1900 rpm for ⌚ **00:00:10** , and spin down briefly in a plate centrifuge (⌚ **00:00:05** at ⚙ **2500 rpm**).

21 Amplify samples using the following PCR conditions:

⌄ 98 °C for ⌚ 00:00:30

35 cycles of:

⌄ 98 °C for ⌚ 00:00:20

⌄ 55 °C for ⌚ 00:00:15

⌄ 72 °C for ⌚ 00:01:00

⌄ 72 °C for ⌚ 00:05:00

Normalization

22 Normalize samples using a SequalPrep plate according to manufacturer's instructions.

☐ [sequalprep_platekit_man.pdf](#)

23 Elute in ☐ 20 µl of SequalPrep Elution Buffer.

Pooling

24 Pool ☐ 10 µl of each sample in a trough, mix well and transfer material to a ☐ 1.5 ml non-stick tube.

25 Purify using AMPureXP beads at a 0.7x ratio. Elute library in ☐ 20 µl of EB.

Library QC

26 Perform final QC on pool by determining concentration (PicoGreen or Qubit assay). Prepare 2 nM pool dilution, based on the sample concentration as determined by PicoGreen and fragment size (expected size is ~555 bp).

Sequencing

27 Dilute pooled sample to 8 pM in HT1, following MiSeq loading instructions, spike in 5% 8 pM PhiX, and load in MiSeq 2x250 or 2x300 reagent cartridge.